
LER 2788
Ages 3+ | Grades PreK+

Great for centers and exploration, these jumbo test tubes will have little scientists making big discoveries!

- Test tubes are calibrated in 25 ml increments
- The top line indicates 100 ml mark
- Colored lids are write-on/wipe-off for labeling

Care and Use

- Always test write-and-wipe markers on the inside corner of the lid to make sure that marks can be erased
- Use a soft cloth or paper towel to wipe away dry-erase markings
- Hand-wash test tubes and lids in soapy warm water
- Allow tubes to dry in the storage rack without lids
- Do not use test tubes with boiling water!
- The ABS plastic test tubes are designed for use with liquid and air temperatures ranging from 25°F to 200°F (-6.5°C to 93°C)
- If freezing water inside the test tube, allow room for frozen water to expand before adding the lid
- Caution: Dropping frozen or extremely cold test tubes may crack the tubes

- What to do:**
1. Have students fill each of the 6 test tubes three-quarters (75 ml) full of water.
 2. Ask students to predict and sketch what each test tube of water would look like if they were to add one material to each: sand, liquid soap, vegetable oil, sugar, food coloring, and small rocks or pebbles.
 3. Have students add enough sand to one of the tubes to raise the water level to the top line (100 ml). Similarly, add enough pebbles to one of the tubes to raise the water level to the top line. Add a spoonful of sugar to one of the tubes, and a spoonful of vegetable oil to another. Add just one drop of the liquid soap and food coloring to their respective tubes. Make sure each tube has water plus only one other added material.
 4. Label the tube lids using a dry-erase marker to identify each water mixture.
 5. Tighten labeled lids on each matching tube.
 6. Have students observe, sketch, and discuss how the water in the tubes looked after they added each material. Have students compare their discoveries to their original predictions.

Additional supplies you will need:

- Water
- Vegetable oil
- Sugar
- Plastic spoons
- Liquid soap
- Food coloring
- Paper and pencil

Starter Experiment—Mix and Match

4. Fill two more test tubes to the first line (25 ml) with water.
5. Add three drops of yellow food coloring to one of the new tubes.
6. Add three drops of blue food coloring to the other tube.
7. Place lids tightly on both tubes and shake for 30 seconds. Put both tubes in the rack.
8. Remove the lids from all 3 colored-water tubes.
9. Pour some of the red water into both of the yellow and blue water tubes so that the water rises just above the halfway lines (50 ml) in each.
10. Place the lids back on the tubes and shake for 30 seconds each. What new colors were made by combining colored water?
11. Try mixing other combinations of colored water. Guide students to make a chart to see if they can think of simple formulas for creating new colors.

Look for these other great products from Learning Resources®:

LER 2774 Jumbo Magnifiers
LER 2781 Big View Bug Jars
LER 2784 Primary Science Set

www.LearningResources.com

Visit our website to write a product review or to find a store near you.

© Learning Resources, Inc., Vernon Hills, IL (U.S.A.)
Learning Resources Ltd., King's Lynn, Norfolk (U.K.)
Please retain our address for future reference.
Made in China. LRM2788-GUD
Fabriqué en Chine.
Made in China.
Informations à conserver.
Bitte bewahren Sie unsere
Adresse für spätere
Nachfragen auf.
Conservar estos datos.
Hecho en China.

Color Mixing

Additional supplies you will need:

- Water
- Food Coloring

What to do:

1. Fill one test tube to the top line (100 ml) with water.
 2. Add five drops of red food coloring into this tube.
 3. Place the lid on tightly and shake for 30 seconds.
- Put the tube in the rack for later.

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