

INTENDED USE

Management of edema and fracture healing, acute ankle sprains, soft tissue injuries of the lower leg, stress fractures of the lower leg, stable fractures of the foot and ankle, stable delayed union of non-union fractures of the distal tibia and fibula, Achilles tendon repairs, Bunionectomies, metatarsal fractures, forefoot, midfoot and hindfoot injuries.

This product was designed to complement the variety of medical treatments common to the above afflictions. The frequency and duration of use should be determined by your prescribing healthcare professional.

CLEANING INSTRUCTIONS

- Remove the foam liner from the walker frame prior to washing.
- Hand wash the foam liner in cold water with mild soap.
- Do not dry clean or tumble dry.
- Do not use bleach.
- Air dry only and do not expose this brace to any heat source.



UNITED-ORTHO.COM | 800-227-8748

FITTING INSTRUCTIONS

Air Stabilizer Standard



Please read carefully.
Correct application
is vital to the proper
functioning of the brace.



Step 1:

Open all the straps.
Open liner, leaving
liner in the boot.

*Additional comfort pads
provided to fill gaps in liner
and provide additional
padding where needed.



Step 2:

Step into boot liner.
Slide heel all the
way back.

*Sewn in comfort pad should
rest on the top of the foot.



Step 3:

Close boot liner
starting at toes and
working up.

Remove pull tab
strips from uprights
to secure liner.



Step 4:

Fasten the straps
starting at toes and
working up.

*Toes may extend past end of liner.
If toes extend past rigid shell, a
larger size is recommended.

*If numbness, tingling or other
discomfort is experienced, loosen
the strap until comfortable. If any
of these symptoms persist, contact
your medical professional.



Step 5:

To inflate: Turn
"inflate/deflate" knob
clockwise until it stops
turning. Press and
release air bulb several
times until the air
liner provides
desired support.

*Liner may expand during
initial use. Additional pumps
may be required to return to
desired support.



Step 6:

To deflate: Turn
"inflate/deflate"
knob counter-
clockwise until it
stops turning.

*If air liner is fully deflated,
it may take upwards of 50
pumps to feel the desired
support.