BeneCheck™ Uni

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
BSU010,BSU01D

**BeneCheck Meter Kit Meter (Front Side & Back Side) / BeneCheck Messgerät
Messgerät (Vorderseite & Rückseite) / Kit Medidor BeneCheck Medidor (Parte delantera y
parte trasera) / Kit de l'Appareil BeneCheck Appareil (Face Avant et Face Arrière)**

LCD Screen Display: Information and test result display

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Chapter 1	Introduction	1
Chapter 2	Setting the Meter	
2.1	Installing/ Replacing the Batteries.....	2
2.2	Set the Date and Time	2
2.3	Code the Meter.....	3
Chapter 3	How to Perform the Test	
3.1	Perform the Test.....	4
Chapter 4	Meter Memory Function	7
Chapter 5	Control Solution Test	8
Chapter 6	Transmission Function for BSU01D	9
Chapter 7	Care and Maintenance	
7.1	Storing Your Meter and Strip	10
7.2	Cleaning and Caring for Your Meter	10
Chapter 8	Error Message and Trouble Shooting.....	11
Chapter 9	Specification	13

Introduction

Please read carefully before using meter kit, and consult healthcare professional before making any important medical decision. Please contact your local customer service for further assistance with the product.

Normal Environmental Conditions

This meter designed under the following conditions:

- Indoor use
- Overvoltage category II
- Pollution degree 2

Intended Use

- *In vitro* diagnostic use only.
- Measuring uric acid in fresh capillary whole blood from fingertip.
- The meter can be used by laypersons or healthcare professionals.
- Not for neonatal use.

Principles of the Examination Method

- Electrochemical biosensor technology.

The meter is plasma-calibrated by reference instruments, which are traceable to the following standard reference materials and methods.


Test	Standard	Method
Uric Acid	NIST SRM 913	Uricase / UV

Chapter 2 Setting the Meter

2.1 Installing/ Replacing the Batteries

This meter uses AAA battery*2. Please remove the plastic tab under the battery before using.

Note: Dispose the batteries according to your local environmental regulations.

2.2 Set the Date and Time

Press "S" Button (3 secs) → One "Beep" Sound → Setting Mode → Turn Off Automatically after Setting




- Setting order: Year/ Month/ Date/ Hour/ Minute/ Unit
- Press "M" button to advance one unit, "S" button to enter next setting.

Note: Correct setting is important while managing your health records.

2.3 Code the Meter

- Code your meter when you first use it or open a new retail box of strips.
- Make sure the meter is off before you insert the code strip.
- Make sure the code on screen, code strip, and strip vial label or single foil pack are the same.


Chapter 3 How to Perform the Test

Materials you need to perform the test:


BeneCheck Meter/ BeneCheck Test Strip/ Lancing Device/ Lancets/ Tissue or Cotton Ball with 75% Ethanol or Disinfection Wipes

3.1 Perform the Test


- Wash and clean your hands with disinfection wipes, and make sure your hands are dry before testing.


Insert lancet firmly.



Remove and save the protective cap.


Recap lancing device and adjust penetration depth.


Take a strip and close the vial immediately.


Insert the strip.


Make sure the code number is correct.


Pull the sliding barrel till it clicks.

Press the release button to sample.


Wipe off first drop of blood and start sampling.

Touch the blood sample with strip.

Uric Acid


Fill up ok


Fill up FAILED


Note: Insufficient fill up of blood sample could lead to inaccurate or failed test result. Do not refill the test strip.


Result will show after countdown. Then discard the used strip to biohazard container.


The meter will turn off automatically.

Pierce used lancet into protective cap.

Discard used lancet to biohazard container.

Recap lancet device and storage.

Note: Please refer to Lancing Device Instructions for detail procedure.

Note:

- Please finish the test within 3 minutes or the meter will turn off automatically.
- The meter will not turn on if you insert the wrong end or wrong side of strip.
- Marked open date on new open strip vial, do not use expired strip.

- Dropping, bumping or other violent impact will damage the meter or cause malfunction.
- Do not use the meter in an environment with possible magnetic, electromagnetic, and radioactive interferences.
- Do not inhale or swallow.
- Education level at least 15 years old and 8 years intensive reading experience (school), no maximum.
- The electromagnetic environment should be evaluated prior to operation of the device.

Warning:

- Please follow local regulations to discard used test strips and lancets.
- Used test strips, lancets and any other material that has been in contact with blood should be treated as potential biohazards.
- If user has infectious disease, the used test materials could be sources of infection.
- Lancets cannot be reused.
- Always use certified lancets to ensure safety.
- Keep the system away from children and pets.
- Use of this instrument in a dry environment, especially if synthetic materials are present (synthetic clothing, carpets etc.) may cause damaging electrostatic discharges that may cause erroneous results.
- Do not use this instrument in proximity to sources of strong electromagnetic radiation, as these may interfere with the proper operation.
- This equipment is designed for use in HOME HEALTHCARE ENVIRONMENT and PROFESSIONAL HEALTHCARE FACILITY ENVIRONMENT. If it is suspected that performance is affected by electromagnetic interference, correct operation may be restored by increasing the distance between the equipment and the source of the interference.

Chapter 4

Meter Memory Function


- Uric Acid - 500 results

The latest test result will replace the oldest when the records exceed maximum memory capacity. The memories start record from M1 to M500.

Directions for Checking Memories:

No Test Strip in the Meter → Press "M" Button (3 secs) → Full Display on Screen → A Short "Beep" → Enter the Memory Mode → Press "S" Button to Confirm → Press "S" Button for Next Test Record/ Press "M" Button for Previous Record → Press "M" Button (3 secs) to Turn Off

Memory Records are shown as follows:

Chapter 5 Control Solution Test


Control solution is used to check the performance of the kit.

The meter kit should be checked:

- When the meter and strip do not work properly.
- When the test result is unusual or inconsistent.

Control solution range is shown as follow:

Please check your strip vial label or single foil pack for exact range.

Perform a Control Test:


Insert an Unused Strip → Make Sure the Code Number is Correct → Press "M" Button (3 secs) → Enter Control Solution Mode (Shown CL on Screen) → Shake the Control Solution Well → Discard First Three Drops → Put One Drop onto a Clean Surface → Touch the Control Solution with Strip Sample Inlet → Test Result Display after Count Down → Compare the result with the control range on test strip vial or single foil pack

Note:

- Do not reuse the test strip.
- Marked open date on new open control solution.
- Do not use the meter if the control test is out of range.
- If the control test keeps result wrongly, please contact an authorized distributor.
- Control solutions are not included. Please contact an authorized distributor for purchasing.

Chapter 6 Transmission Function for BSU01D

The meter provides Bluetooth transmission function. It can transmit test results to connected device wirelessly via Bluetooth.

The screen show "Bluetooth" means Bluetooth transmission function is working.

- The meter with Bluetooth 4.0 can apply to IOS 6, Android 4 and above.
- This meter complies with the emission and immunity requirements of IEC 61326-2-6.
- The meter and the transmission function may be interfered while other device is operating nearby. Ex: mobile phone, wireless internet, etc.
- If the transmission has been interfered. Please keep the meter away from the source of interference or turn off the interfering device.
- Please make sure the meter and the receiving device are placed within a reasonable distance (less than five meters) during transmission via Bluetooth.
- Please do not share the meter with other people if you are using Bluetooth transmission function. The test results from other people will also transmit to the receiving device, and it will influence your test record.
- The Bluetooth transmission function may not work on certain types of mobile phone due to the compatibility of Android systems.
- Frequency range: 2402 – 2480MHz. RF power: -2.21 dBm @ room temperature Typical.

Warning: If the meter is paired with a mobile device and Bluetooth is on, your test results will be sent to the paired device. Please check that the meter results are the same as those displayed on the paired device.

Chapter 7 Care and Maintenance

7.1 Storing Your Meter and Strip

Meter:

- Avoid bump or violent behavior.
- Do not use in extremely dry environment. It may cause static discharges.
- Do not use under electromagnetic radiation, ex: electrical equipment.
- Do not disassemble the meter for any reason.
- Keep meter clean by wiping the exterior appearance with tissues or lint-free cloth.
- No modification of this equipment is allowed.
- Do not expose the meter under lint, dust, sunlight, heat or humid environment.

Strip:

- Do not store in high humidity environment, or expose directly to sunlight.
- Do not freeze or refrigerate the meter and strips.
- Keep your hand dry and clean while handing the strips and performing the test.
- Do not bend, cut or fold the strips.

7.2 Cleaning and Caring for Your Meter

Please use soft cloth slightly damp with one of the following solution to clean meter surface after finish the test.


- 75% alcohol
- Super Sani-Cloth disposable wipes
- Mild dishwashing liquid with water
- 10% household bleach solution and 90% water

Note:

- Do not allow any other wet cloth or liquid.
- Do not allow any liquid run in or around the test port and battery cover.
- Make sure the meter is completely dry before use.
- Protection impairment if used in a manner not specified by the manufacturer.

Chapter 8

Error Message and Trouble Shooting

Message	Cause	Solution
	• Problem with code strip. • Problem with test strip. • Insert strip improperly.	Repeat the coding procedure (Chapter 2.3) and insert the strip again. If the problem persists, please contact local distributor for service.
	• Low battery.	Replace with new battery.
	• Problem with code strip or meter. • Insert code strip improperly.	Repeat the coding procedure (Chapter 2.3). If the problem persists, please contact local distributor for service.
	• Incorrect meter operating temperature.	Repeat the test after meter return to operating temperature. If the problem persists, please contact local distributor for service.
	• Used strip. • Damped strip.	Follow Chapter 3.1 and repeat the test with a new strip. If the problem persists, please contact local distributor for service.

A digital display showing the text "E-9" in a black, segmented font.

- Test incomplete due to removing the strip during measuring.

Follow Chapter 3.1 and repeat the test with a new strip. Do not remove the strip before the test is completed.

A digital display showing the text "E-1" in a black, segmented font.

- Improper code strip.

Repeat the coding procedure (Chapter 2.3). If the problem persists, please contact local distributor for service.

A digital display showing the text "HI" in a black, segmented font.

- Test result is higher than the range listed on Chapter 9.

Follow Chapter 3.1 and repeat the test with a new strip. If the problem persists, please contact local distributor for service.

A digital display showing the text "Lo" in a black, segmented font.

- Test result is lower than the range listed on Chapter 9.

Follow Chapter 3.1 and repeat the test with a new strip. If the problem persists, please contact local distributor for service.

Chapter 9 Specification

Test Sample	Fresh Capillary Whole Blood from fingertip
Measuring Time	15 seconds
Measuring Range	3-20 mg/dL (0.18-1.19 mmol/L)
Sample Volume	1 μ L
Storage & Transportation Condition	Meter: 0-60 °C (32-140 °F) Strip: 4-30 °C (39-86 °F)
Operation Temperature	10-40 °C (50-104 °F)
Storage & Transportation Relative Humidity	10-90%
Open Vial Relative Humidity	10-90 %
Memory	500 Test Results
Battery Type	1.5V AAA battery * 2
Battery Life	Approximately 1,000 tests
Dimensions	62*88*24 mm
Weight	About 64g (with battery)
Altitude	10,000 feet (3048 m) (700~1013 hpa)
Expected Service Life	5 years (Approximately 10,000 times)
Transmission Function	Bluetooth 4.0 for BSU01D

Note:

- Please refer to the strip insert for accuracy, precision, limitation, and other important information.

Contents of the Kit (please check the meter outer box for exact detail)

- BeneCheck Uni Uric Acid Monitoring Meter (BSU010,BSU01D) (with AAA battery*2)
- User's Manual
- Quick Guide
- Carrying bag
- Lancing Device

(Please refer to package for manufacture information.)

Optional (not included in the standard kit package, please contact your local distributor for ordering)

- BeneCheck Uric Acid Test Strip (BK-U1) (with code strip & Insert)
- BeneCheck Uric Acid Control Solution
- Lancets (Please refer to package for manufacture information.)

Note: After purchasing, if the contents are damaged, please contact authorized distributor immediately.

Hereby, General Life Biotechnology Co., Ltd. declares that the radio equipment type BSU01D is in compliance with Directive 2014/53/EU.

The full text of the EU declaration of conformity is available at the following internet address:

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PHARMADOC GmbH
MDSS GmbH