SMARTHEART™

Blood Pressure Monitor with Adult Wrist Cuff

Model # 01-523 series

INSTRUCTION MANUAL

ENGLISH AND SPANISH

Please read this instruction
manual completely before
operating this unit.

STOP!
PLEASE ENSURE YOU HAVE ALL OF
THE FOLLOWING COMPONENTS
BEFORE USING YOUR DIGITAL
BLOOD PRESSURE MONITOR


Digital Monitor

Storage Case

Charging Cable

Instruction Manual
Quick Start Guide
Blood Pressure Reading Tracking Log
FAQ Information Sheet

IF YOU ARE MISSING ANY PARTS,
INCLUDING INSERTS OR INSTRUCTION
MANUALS, DO NOT RETURN TO PLACE OF
PURCHASE. CONTACT CUSTOMER CARE AT
866-326-1313.

INDEX

Care & Safety Information	4-9
Introduction	10
Monitor Features	11
Getting Started	12-13
Charging the Battery	14
Setting Date/Time	15-17
Fitting & Applying Cuff	18-19
Taking Your Blood Pressure Reading	20-23
Interpreting Your Results	24-25
Memory Functions	26-27
Care & Maintenance	28-29
Device & Label Symbols	30
Display Symbols	31
Troubleshooting	32
Error Codes	33
FCC Statement	34
Electromagnetic Compatibility	35-37
Specifications	38-39
Warranty	40-41
Instrucciones en Español	43-88

**Toll-Free Customer
Care Help Line:
1-866-326-1313**

**Monday – Friday
8:30 a.m. –
4:30 p.m. CST**

SmartHeart™
Manufactured for
Veridian Healthcare
1175 Lakeside Drive
Gurnee, IL 60031

Made in China

#93-1276 03/22
©2022 Veridian Healthcare, LLC

CARE & SAFETY INFORMATION

NOTE: Read all instructions carefully before use. The following basic precautions are needed when using an electrical product.

CAUTION: Failure to read and observe all precautions could result in personal injury or equipment damage.

Improper care or use of your blood pressure monitor may result in injury, damage to the unit or ineffective treatment. Following these instructions will ensure the blood pressure monitor's efficacy and long life.

GENERAL CAUTIONS AND WARNINGS

- Measurement results should be discussed with your physician or healthcare professional; never self-diagnose or attempt treatment as this may be dangerous.
- Intended for adult use only; this device is not approved for infant or child.
- Keep out of reach of children.
- Do not use the unit if it has any damaged parts, if it has been submersed in water or dropped.
- If any abnormality occurs, discontinue use until the unit has been examined and repaired.
- Only use included attachments and parts; do not use attachments from other brands or models as these may not be properly calibrated for use with this device and may result in measurement error.

CARE & SAFETY INFORMATION

- Do not share with an infected person.
- Do not disassemble the unit or attempt to repair it; substitution of a component different from that supplied may result in measurement error and will void manufacturer warranty.
- Always follow local regulations for proper disposal of the monitor, cuff and batteries.
- This device is intended for non-invasive measuring and monitoring of arterial blood pressure.
- It is not intended for use on extremities other than the arm or for functions other than obtaining a blood pressure measurement.
- To verify the calibration of this device, please contact the manufacturer.
- This device may be used only for the purpose described in this booklet. The manufacturer cannot be held liable for damage caused by incorrect application and use.
- The service life of the cuff may vary by the frequency of washing, skin condition of the user, and adherence to proper storage conditions. The typical service life is 10000 measurements.
- It is recommended that the performance should be checked every 2 years and after maintenance and repair, by retesting at least the requirements in limits of the error of the cuff pressure indication and air leakage (testing at least at 50mmHg and 200mmHg).

CARE & SAFETY INFORMATION

- The device doesn't need to be calibrated within two years of reliable service.
- Keep the unit out of reach of infants, young children or pets to avoid inhalation or swallowing of small parts. It is dangerous or even fatal.

OPERATING CAUTIONS AND WARNINGS

- Prolonged over-inflation may cause congestion, swelling or bruising in some people.
- Wait a minimum of 2 minutes between measurements (20-minutes is recommended); excessive frequent measurements may restrict blood circulation.
- Do not use this device while connected to any monitoring equipment on the same limb.
- Consult with your healthcare professional for guidance for use if you have any of the following: intravascular shunt; mastectomy on the arm corresponding to the measurement side of the body.
- Consult with your healthcare professional before using the device in pregnancy including pre-eclampsia, or if diagnosed with arrhythmia or arteriosclerosis.
- If you experience any discomfort or an abnormality occurs, stop using the device immediately; to stop

CARE & SAFETY INFORMATION

the device during operation, press the **U1 or U2 (POWER)** buttons and the cuff will automatically deflate.

- Operation outside of stated operating temperature may result in measurement error or device malfunction; operation environment temperature is: 41°F – 104°F (5°C – 40°C); Humidity: 15% to 90% RH (non-condensing).
- Never use this unit while operating a vehicle or in the bathtub or shower.
- To avoid measurement errors, please avoid the condition of strong electromagnetic field radiated interference signal or electrical fast transient/burst signal.
- The blood pressure monitor, its adaptor, and the cuff are suitable for use within the patient environment. If you are allergic to polyester, nylon or plastic, please don't use this device.
- During use, the patient will be in contact with the cuff. The materials of the cuff have been tested and found to comply with requirements of ISO 10993-5:2009 and ISO 10993-10:2010. This cuff will not cause any potential sensation or irritation reaction.
- Wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walkie-talkies can affect this equipment and should be kept at least a distance d away from the equipment. The

CARE & SAFETY INFORMATION

distanced is calculated by the MANUFACTURER from the 80MHz to 5.8 GHz column of Table 4 and Table 9 of IEC 60601-1-2:2014, as appropriate.

- The equipment is not AP/APG equipment and not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.
- On the rare occasion of a fault causing the cuff to remain fully inflated during measurement, open the cuff immediately. Prolonged high pressure (cuff pressure >300mmHg or constant pressure >15mmHg for more than 3 minutes) applied to the arm may lead to an ecchymosis.
- If the cuff pressure reaches 40 kPa (300 mmHg), the unit will automatically deflate. Should the cuff not deflate when pressure reaches 40 kPa (300 mmHg), detach the cuff from the arm and press any button to stop inflation.
- The operator shall not touch output of batteries / adapter and the patient simultaneously.
- The plug/adapter plug pins insulates the device from the main supply. Do not position the device in a position where it is difficult to disconnect from the supply mains to safely terminate operation of ME equipment.

CARE & SAFETY INFORMATION

STORAGE CAUTIONS AND WARNINGS

- Storage outside of stated storage temperature may result in measurement error or device malfunction; storage environment temperature is: -4°F – 140°F (-20°C – 60°C); Humidity: ≤93%RH (non-condensing) at water vapor pressure to 50hPa.
- Keep the unit out of reach of children.
- After charging, unplug the adapter; do not store while plugged in or with adapter attached.
- When not in use, store the device with the adapter in a dry room and protect it against extreme moisture, heat, lint, dust and direct sunlight. Never place any heavy objects on the storage case.

CLEANING CAUTIONS AND WARNINGS

- Never immerse the unit in water to clean as it may damage the unit.
- Follow the 'Care & Maintenance' portion of this manual for instruction on how to clean and care for your monitor.
- Don't use any abrasive or volatile cleaners.

INTRODUCTION

It is recommended that you first seek the advice and recommendation of your physician or healthcare professional when using home diagnostic devices, including blood pressure monitors.

The Blood Pressure Monitor is digital monitors intended for use in measuring blood pressure and heartbeat rate with wrist circumference ranging from 5.3" to 8.4" inches (about 13.5 – 21.5 cm). It is intended for adult indoor use only.

Automatic digital blood pressure monitors use the oscillometric method to electronically measure your blood pressure. The monitor detects your blood's movement through the artery in your arm and converts the movements into a digital reading. The oscillometric method does not require a stethoscope, making the monitor ideal for home use.

Blood pressure readings determined with this device are equivalent to measurements obtained by a trained healthcare professional using the cuff/stethoscope auscultation method, within the limits prescribed by the American National Standard for Electronic or Automated Sphygmomanometers.

Thank you for purchasing a SmartHeart Automatic Digital Blood Pressure Monitor. With proper care and use, your monitor will provide you with many years of reliable readings.

**PLEASE READ THIS INSTRUCTION MANUAL
COMPLETELY BEFORE OPERATING THIS UNIT.**

MONITOR FEATURES

DETAILS OF THE DIGITAL DISPLAY SCREEN AND INTERPRETATION OF THE SYMBOLS ARE INCLUDED THROUGHOUT THIS MANUAL.

GETTING STARTED

Before taking a measurement

**NOTE: THIS RECHARGEABLE DEVICE IS
'ASLEEP' FOR SHIPPING AND TRANSPORT.**

To activate your monitor, press and hold the **SET** button for several seconds. The display will illuminate and your monitor is ready for set-up and use.

1. Your readings should only be interpreted by your physician or healthcare professional with access to your individual medical history. Regular use of a home blood pressure monitor will allow you to track and record your readings for discussion with your physician.
2. Conduct your measurement in a quiet place while seated in a relaxed position. Rest for 15 minutes before taking your reading. If an error occurs or you desire to take a second reading, allow 15 minutes between readings for your blood vessels to return to normal.
3. Avoid smoking, eating, taking medication, alcohol, physical activity or any other stressful activity for 30 minutes prior to taking a reading.
4. Always remove any jewelry or constrictive clothing that may interfere with the cuff placement.

GETTING STARTED

5. Keep yourself and the monitor still during measuring; do not talk during the reading.
6. It is recommended that you take your readings at the same time each day to better monitor any indications in your results.
7. Record your daily measurement on the included chart or some other written document to share with your physician.


Contraindications

1. The device should not be used by any person who may be suspected of, or is pregnant.
2. The device is not suitable for use on patients with implanted, electrical devices, such as cardiac pacemakers, defibrillators.

CHARGING THE BATTERY

This device includes an internal rechargeable battery within the electrical stimulator unit and USB cable. A USB wall charger/adapter is required for charging; not included.

It is necessary to charge the battery when 'Lo' + low battery symbol appears on the display, when the display dims or when the display does not illuminate.


Before first use, the device should be fully charged. Charging time may take up to 90 minutes from depleted battery.

1. Connect the wide end of the USB cable to a wall adapter.
2. Connect the smaller mini USB to the wrist monitor.
3. During charging the digital display will show a battery flashing to indicate the current charge level. When the device is fully charged, the battery symbol will be stop flashing and remain solid on the display.
4. Unplug the charger and stimulator when fully charged; never operate the unit while charging.

SETTING DATE/TIME

The date/time feature is on the bottom right of the display screen. This unit stores the date, time and year allowing for a date/time stamp on all readings in the memory recall function. The **SET**, **MEM** and **START STOP** buttons are located on the bottom center of the monitor.

SETTING DATE/TIME

HOW TO SET DATE AND TIME

1. When the unit is off (a blank display screen), press and hold the **SET** button for 3-5 seconds to enter clock setting mode.

2. **HOURS/MINUTES**—first the hour will flash on the display. Press the **MEM** button to advance the display to the desired hour, press the **SET** button to confirm your selection and advance to the minute setting. Press the **MEM** button to advance to the correct minute, press **SET** to confirm.

NOTE: the time will show PM when the set time is PM, but the same area will be blank for AM.

3. **MONTH**—next the Month will flash on the display. Press the **MEM** button to advance the display to the desired month, press the **SET** button to confirm your selection and advance to the date settings.
4. **DATE**—repeat the same sequence of step 2 to confirm the MONTH.
5. **YEAR**—repeat the same sequence of step 2 to confirm the YEAR.
6. After confirming YEAR, the LCD will display "dONE" and then shut off.

SETTING DATE/TIME

PLEASE SET DATE AND TIME BEFORE INITIAL USE TO ENSURE YOUR READINGS ARE LABELED AND RECORDED PROPERLY.

		
---	---	---

Hour → Minutes → Month →

		
---	---	---

Date → Year → Complete →

FITTING & APPLYING YOUR CUFF

PROPER CUFF FIT AND PLACEMENT IS CRITICAL IN OBTAINING ACCURATE BLOOD PRESSURE MEASUREMENTS. PLEASE READ THESE INSTRUCTIONS CAREFULLY AND CONSIDER CONSULTING WITH YOUR PHYSICIAN FOR FURTHER DIRECTION.

Before applying your blood pressure cuff, be sure your wrist is within the appropriate cuff range:


Wrist Circumference Range Suitable for 5.3" – 8.4"
(13.5 cm – 21.5 cm).

Reminder, it is important to avoid smoking, eating, taking medication, alcohol consumption or physical activity 30 minutes prior to taking a reading.

If for any reason you are unable to or cannot not use your left wrist, please modify the instructions for cuff application to your right wrist. Your physician can help you identify which wrist is best for you to take measurements from.

Applying Your Cuff

1. Remove any constrictive clothing or jewelry that may interfere with cuff placement.

FITTING & APPLYING YOUR CUFF

2. Be seated at a table or desk with your feet flat on the floor.
3. Hold your left arm in front of you with your palm facing upward.
4. Apply the cuff to your left wrist so that the digital display face is positioned on the inside area of your wrist facing you, Fig. 1.
5. Adjust the cuff approximately 0.5" from the edge of the head of the ulna bone, Fig. 2.
6. The cuff should fit comfortably, yet snugly around your wrist.


Fig. 1


Fig. 2

If you are not comfortable with applying your cuff, please seek the assistance of another member of your household or work with your physician to practice the cuff application. Incorrectly applied cuffs may result in inaccurate readings.

TAKING YOUR BP READING

Please read the preceding portions of this manual prior to taking your first reading.

Reminder: It is important to avoid smoking, eating, taking medication, alcohol consumption or physical activity 30 minutes prior to taking a reading.

Select the desired user, 1 or 2, before measurement.

The user icon displayed represents the memory bank where the measurement result will be stored.

1. When the monitor is OFF, press and hold the **MEM** button. The display will illuminate and the current user, 1 or 2, will flash on the display.
2. Press **MEM** again to toggle between user 1 or 2.
3. Press **SET** to confirm the selection.

Proceed to the measurement steps below.

TAKING YOUR BP READING

Taking the measurement.

1. Rest your elbow on a solid surface with your palm facing upward. Elevate your arm so that the cuff is at the same level as your heart, Fig 1. Be certain that cuff is elevated to the proper level; you may choose to support your hand with a folded towel or book. Relax your left hand. Be seated with your feet flat on the floor and legs uncrossed.

Fig. 1

2. Press and release the **START STOP** button, the unit will run a self-test, Fig. 2.
3. As the inflation begins, the display will briefly show 0.
4. The cuff will automatically inflate until the reading is detected. You will see the numbers increase on the screen and feel the cuff inflate slowly.


Fig. 2

TAKING YOUR BP READING

YOU CAN STOP THE INFLATION OR DEFLATION PROCESS AT ANY TIME BY PRESSING THE START/STOP BUTTON.

- When the measurement is complete, the cuff will automatically deflate, and your blood pressure measurement and pulse reading results will show on the display screen, Fig 3.
- The Hypertension Indicator will indicate your reading range on the left display in comparison to the color bars on the face of the monitor, Fig 4. See page 24 for information regarding interpreting the Hypertension Indicator results.


Fig. 3

Fig. 4

TAKING YOUR BP READING

7. If an irregular heartbeat was detected during the reading, the Irregular Heartbeat Detector symbol will appear on the display, Fig 5. 

See page 25 for more information on the Irregular Heartbeat Detector feature.

8. The reading will automatically be stored in memory, up to 60 readings.
9. Press the **START STOP** button to turn the unit off to conserve energy and battery life. The unit will also automatically shut-off after 1 minute of non-use.




Fig. 5

INTERPRETING YOUR RESULTS

Irregular Heartbeat Detector

An irregular heartbeat is defined as a heartbeat rhythm that is more than 25% slower or 25% faster than the average measurement detected while the monitor is conducting the measurement. When this device detects an irregular rhythm two or more times during the reading, the IHB symbol will appear on the display screen.

This monitor will operate if the icon appears, but the results may be compromised, especially if this symbol appears often. Please consult with your physician or trained healthcare professional for further information regarding an irregular heartbeat and if this symbol appears frequently.

INTERPRETING YOUR RESULTS

Hypertension Indicator

This unit features our unique Hypertension Indicator. The color bars on the left side of the monitor display correspond with an icon on the digital display, indicating where the measurement results fall within the American Heart Association standards.

This chart should only be considered as a guideline; always consult with your physician to interpret your individual results. Never adjust your medication(s) dosage, or severely alter your diet or exercise routine without consulting with your physician.

Hypertension Indicator

Category	Systolic (upper #)	Diastolic (lower #)
Hypertensive Crisis	Higher than 180	&/or Higher than 120
High Blood Pressure, Stage 2	140 or higher	or 90 or higher
High Blood Pressure, Stage 1	130-139	or 80-89
Elevated	120-129	& Less than 80
Normal	Less than 120	& Less than 80

MEMORY FUNCTIONS

This monitor stores and recalls up to 60 readings per user. These measurements can be shared with your physician or trained healthcare professional.

Select the desired user, 1 or 2, for review.

The user icon displayed represents the memory bank where the measurement results are stored.

1. When the monitor is OFF, press and hold the **MEM** button. The display will illuminate and the current user, 1 or 2, will flash on the display.

2. Press **MEM** again to toggle between user 1 or 2.
3. Press **SET** to confirm the selection.

Proceed to the measurement recall steps below.

Recalling Measurements in Memory

1. Press the **MEM** button; the display will show the average of the last three measurements in that memory bank, Fig 1.
2. Press the **MEM** button again to begin scrolling through the saved measurements. Measurements will appear on the display from most current to oldest; the date / time shown will indicate when the measurement was taken.

Fig. 1

MEMORY FUNCTIONS

3. All results for a given measurement will display, including measurement results, pulse rate, Hypertension Indicator and Irregular Heartbeat alert (if applicable), Fig 2.
4. When the number of readings exceeds 60, the oldest data will be replaced with the new record.
5. Press **START STOP** to turn the monitor off at any time during review of the stored measurements.

Fig. 2

Clearing Measurements From Memory

To delete all results in a given bank:

1. When the monitor is off, press and release the **MEM** button to access the memory bank feature.
2. Press and hold both the **MEM** and **SET** buttons for 3 seconds.
3. The display will show 'dEL' and 'dONE', Fig 3.
4. The memory bank is now empty, Fig 4, and the monitor will turn off.

Fig. 3


Fig. 4

CARE & MAINTENANCE

Proper care and maintenance of your blood pressure monitor will help ensure long life for the product; improper use of the monitor may void the included warranty.

Before taking a measurement

- At least 30 min required for ME equipment to warm from the minimum storage temperature between uses until it is ready for intended use. At least 30 min required for ME equipment to cool from the maximum storage temperature between uses until it is ready for intended use.
- Before use, make sure the device functions safely and is in proper working condition.

Monitor Care & Maintenance

- Do not use any liquids on the monitor; use a soft, dry cloth to clean the monitor as needed.
- Storage outside of stated storage temperature may result in measurement error or device malfunction; storage environment temperature is: -4°F – 140°F (-20°C – 60°C); Humidity: ≤93%RH (non-condensing).
- Avoid high temperatures and direct sunlight.
- Do not attempt to disassemble your monitor; disassembling your monitor may void the manufacturers warranty.
- Do not subject the monitor to strong shocks; take care not to drop the monitor.
- No component can be maintained by user in the monitor. The circuit diagrams, component part

CARE & MAINTENANCE

lists, descriptions, calibrations instructions, or other information which will assist the user's appropriately qualified technical personnel to repair those parts of equipment which are designated repairable can be supplied.

Cuff Care & Maintenance

- The cuff may be spot cleaned with a mild-detergent as needed.
- Never attempt to iron or dry the cuff; if the cuff becomes wet, allow the cuff to air dry completely before use or storage.
- To disinfect the cuff, wipe the inside (skin contact side) of the cuff using a cloth lightly moistened with 70% - 90% ethyl rubbing alcohol and allow to air dry; it is suggested to disinfect the cuff frequently when used in clinical sittings or if the cuff is contaminated in any way.
- It is recommended the performance should be checked every 2 years or after a repair. Please contact customer care at 866-326-1313.
- It is recommended to clean the cuff every 200 uses.
- The monitor can maintain the safety and performance characteristics for a minimum of 10,000 measurements of three years.
- The cuff integrity is maintained after 1,000 open-close cycles of the closure.

DEVICE & LABEL SYMBOLS

These symbols may appear on your device, instructions or packaging and may vary by make and model.

Symbol	Meaning
	Read This Manual—All included manuals should be read prior to first use.
	Warning—Symbol indicates a warning, prohibition or mandatory action that mitigates a risk that is not necessarily obvious to the device operator.
	Type BF Applied Parts—Indicates that a part of this unit comes in contact with the patient in order to carry out it's intended function; in the case of this device the cuff is the Type BF Applied Part.
	Environment Protection—dispose of this product properly; consult with your local recycling ordinances for proper recycling and disposal.
	Keep Dry—This device should be kept dry; never submerge the unit or cuff. Consult with the Care and Maintenance section of this manual for information on cleaning your monitor.
	Serial Number.
	Direct Current.

DISPLAY SYMBOLS

Symbol	Meaning
SYS	Systolic Blood Pressure—The 'top' blood pressure result
DIA	Diastolic Blood Pressure—The 'bottom' blood pressure result
PUL	Pulse Rate—Pulse rate per minute during the measurement
mmHg	Millimeter of Mercury—Unit of blood pressure measurement for this device
	Low Battery—Battery life indicator
	Pulse Reading—Appears when pulse reading is displayed.
	Irregular Heartbeat Symbol—An irregular heartbeat was detected during the measurement
	User 1/User 2—Indicates which user is active
LAST 3 AVG.	Average—The display shows the average of the last three measurements in memory
MEMORY REVIEW	Memory—The display shows the measurements stored in memory

TROUBLESHOOTING

Problem	Possible Cause	Solution
Display shows an abnormal result.	Cuff position was not correct or was not tight enough.	Apply the cuff properly and attempt a new measurement.
	Body posture was not correct during measurement.	Review the 'Taking a Measurement' portion of this manual and attempt a new measurement.
	Speaking, excessive movement, angry, excited or nervous anxiety during measurement.	Wait a period of time and attempt a new measurement after calm has returned; do not speak or move during measurement.
Display is blank when power is on.	Battery charge is low.	Recharge the device following recharging instructions and attempt a new measurement.
Low battery symbol + Lo appear on the display.	Battery charge is low.	Recharge the device following recharging instructions and attempt a new measurement.

ERROR CODES

Code	Meaning	Corrective Action
E1	The cuff is not secure.	Re-apply the cuff and attempt new measurement.
E2	The cuff is too tight.	Re-apply the cuff and attempt new measurement.
E3	The cuff pressure is in excess.	Relax for several minutes and attempt new measurement.
E10 or E11	The monitor detected motion while measuring.	Sit comfortable and remain still; attempt a new measurement.
E20	The measurement process does not detect the pulse signal.	Re-apply the cuff and attempt new measurement.
E21	The unit measured incorrectly.	Relax for a moment and then measure again.
EExx	A calibration error occurred.	Retake the measurement; if the problem persists, contact the Customer Care.
OUT	Out of measurement range.	Relax for a moment and then measure again.

FCC STATEMENT

NOTE:

POTENTIAL FOR RADIO/TELEVISION INTERFERENCE (for U.S.A. only)

This product has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules.

These limits are designed to provide reasonable protection against harmful interference in a residential installation. The product generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If the product does cause harmful interference to radio or television reception, which can be determined by turning the product on and off, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna
- Increase the separation between the product and the receiver.
- Connect the product into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

ELECTROMAGNETIC COMPATIBILITY

Table 1

The ME EQUIPMENT or ME SYSTEM is suitable for home healthcare environments

Warning: Don't near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances is high.

Warning: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Warning: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation."

Warning: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the equipment, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Technical description:

1. All necessary instructions for maintaining BASIC SAFETY and ESSENTIAL PERFORMANCE with regard to electromagnetic disturbances for the expected service life.
2. Guidance and manufacturer's declaration -electromagnetic emissions and Immunity.

Guidance and manufacture's declaration - electromagnetic emissions	
Emissions Test	Compliance
RF emissions CISPR11	Group 1
RF emissions CISPR11	Class B
Harmonic distortion IEC 61000-3-2	Class A
Voltage fluctuations / flicker emissions IEC 61000-3-3	Comply

ELECTROMAGNETIC COMPATIBILITY

Table 2

Guidance and manufacture's declaration - electromagnetic immunity

Immunity Test	IEC 60601-1-2 test level	Compliance level
Electrostatic Discharge IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV signal input/output 100 kHz repetition frequency	± 2 kV for power supply lines ± 1 kV signal input/output 100 kHz repetition frequency
Surge IEC 61000-4-5	± 0.5 kV, ± 1 kV differential mode ± 0.5 kV, ± 1 kV, ± 2 kV common mode	± 0.5 kV, ± 1 kV differential mode ± 0.5 kV, ± 1 kV, ± 2 kV common mode
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % UT; 0,5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0 % UT; 1 cycle and 70 % UT; 25/30 cycles; Single phase: at 0°. 0 % UT; 250/300 cycle	0 % UT; 0,5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0 % UT; 1 cycle and 70 % UT; 25/30 cycles; Single phase: at 0°. 0 % UT; 250/300 cycle
Power frequency (50Hz/60Hz) magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz
Conducted RF IEC 61000-4-6	3 V 0,15 MHz – 80 MHz 6 V in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	3 V 0,15 MHz – 80 MHz 6 V in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz
Radiated RF IEC 61000-4-3	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz

NOTE UT is the a.c. mains voltage prior to application of the test level.

ELECTROMAGNETIC COMPATIBILITY

Table 3

Guidance and manufacture's declaration - electromagnetic immunity

Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	Test Frequency (MHz)	Band a (MHz)	Service	Modulation	Maximum Power (w)	Distance (m)	IEC 60601-1-2 Test Level (V/m)	Compliance level (V/m)
	385	380-390	TETRA 400	Pulse modulation 18Hz	1.8	0.3	27	27
	450	380-390	GMRS 460, FRS 460	FM $\pm$ 5kHz deviation 1kHz sine	2	0.3	28	28
	710 745 780	704-787	LTE Band 13, 17	Pulse modulation 217Hz	0.2	0.3	9	9
	810 870 930	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18Hz	2	0.3	28	28
	1720 1845 1970	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217Hz	2	0.3	28	28
	2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217Hz	2	0.3	28	28
	5240 5240 5785	5100-5800	WLAN 802.11 a/n	Pulse modulation 217Hz	0.2	0.3	9	9

PRODUCT SPECIFICATIONS

Name	Blood Pressure Monitor with Adult Wrist Cuff
Model Number	01-523 series
Display System	Digital display/LCD
Measuring Method	Oscillometric Method, Automatic Air Inflation and Measurement
Power Source	3.7V 420mAH Built-in rechargeable lithium-ion battery, 5V  1A USB cable
Measuring Range	Rated cuff pressure: 0mmHg~299mmHg(0kPa ~ 39.9kPa) Measurement pressure: SYS: 60mmHg~230mmHg (8.0kPa~30.7kPa) DIA: 40mmHg~130mmHg (5.3kPa~17.3kPa) Pulse value: (40-199)beat/minute
Accuracy	Within operating environment 41°F – 104°F (5°C – 40°C): Pressure: ±3 mmHg; Pulse: ±5% of readings
Pressurization	Automatic pressurization by pump
Deflation	Automatic pressure release valve
Memory	Built-in memory with 2 banks, 60 measurements per bank
Automatic Power-Off	Approximately 1 minute after last button operation
Battery Life	Approximately 200 measurements
Operation Environment	A temperature range of :+5°C to +40°C A relative humidity range of 15% to 90%, non-condensing, but not requiring a water vapour partial pressure greater than 50hPa An atmospheric pressure range of: 700 hPa to 1060 hPa

PRODUCT SPECIFICATIONS

Storage Environment	Temperature: -20°C to +60°C A relative humidity range of ≤ 93%, non-condensing, at a water vapour pressure up to 50hPa
Monitor Dimensions	3.14" x 2.85" x 0.5"
Monitor Weight	3.8 oz. (without batteries)
Wrist Circumference	Suitable for 5.3" – 8.4" (13.5 – 21.5 cm)
Accessories	USB cable, instruction manual, storage case, Quick Start Guide, Blood Pressure Log
Degree of protection	Type BF applied part
Protection against ingress of water	IP22, It means the device could be protected against solid foreign objects of 12.5mm and greater and against vertically falling water drops when ENCLOSURE tilted up to 15"
Software Version	A01

Specifications are subject to change without notice

5-YEAR LIMITED WARRANTY

Congratulations on your purchase of a Digital Blood Pressure Monitor. Your Digital Blood Pressure Monitor is covered by the following limited warranty commencing upon the date of purchase, and subject to the following terms and conditions:

The Warranter warrants that its Digital Blood Pressure Monitor will be free from defects in materials and workmanship under normal consumer usage for a period of five years for the original purchaser.

Periodic maintenance, repair and replacement of parts due to normal wear and tear are excluded from coverage. Defects or damage that result from: (a) improper operation, storage, misuse or abuse, accident or neglect, such as physical damage (cracks, scratches, etc.) to the surface of the product resulting from misuse; (b) contact with liquid, water, rain, extreme humidity or heavy perspiration, sand, dirt or the like, extreme heat, or food; (c) use of the Digital Blood Pressure Monitor for commercial purposes or subjecting the Digital Blood Pressure Monitor to abnormal usage or conditions; or (d) other acts which are not the fault of the Warranter, are excluded from coverage. This warranty does not cover batteries or other power sources that may be provided with, or used with the Digital Blood Pressure Monitor.

If the Digital Blood Pressure Monitor fails to conform to this limited warranty, return the Digital Blood Pressure Monitor postage prepaid to: **Attn: Repair Department, 1175 Lakeside Drive, Gurnee, IL, 60031**. When returning a product, please also include: (i) a copy of your receipt, bill of sale or other comparable proof of purchase; (ii) a written description of the problem; and (iii) your name, address and telephone number. Carefully package the product to avoid any damage that may occur while in transit; shipping insurance with returned receipt is recommended. At our option, the Warranter will repair or replace the unit found to be defective in materials or workmanship under normal consumer usage. The purchaser will be notified of any additional repairs required prior to completing the repair, and will

5-YEAR LIMITED WARRANTY

be responsible for parts charges, if any, and repair charges not covered by this limited warranty.

EXCEPT AS PROVIDED FOR IN THIS LIMITED WARRANTY, ALL EXPRESS AND IMPLIED WARRANTIES AND CONDITIONS ARE DISCLAIMED, INCLUDING WITHOUT LIMITATION THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. THE REPAIR OR REPLACEMENT AS PROVIDED UNDER THIS LIMITED WARRANTY IS THE EXCLUSIVE REMEDY OF THE CONSUMER, AND IS PROVIDED IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. IN NO EVENT SHALL THE WARRANTER BE LIABLE, WHETHER IN CONTRACT OR TORT (INCLUDING NEGLIGENCE) FOR DAMAGES IN EXCESS OF THE PURCHASE PRICE OF THE PRODUCT, OR FOR ANY INDIRECT, INCIDENTAL, SPECIAL OR CONSEQUENTIAL DAMAGES OF ANY KIND, OR FOR DAMAGES TO, OR LOSS OF, OTHER PROPERTY OR EQUIPMENT OR PERSONAL INJURIES TO THE FULL EXTENT THESE DAMAGES MAY BE DISCLAIMED BY LAW.

Some states and jurisdictions do not allow the limitation or exclusion of incidental or consequential damages, or limitation on the length of an implied warranty, so the above limitations or exclusions may not apply to you. This warranty gives you specific legal rights, and you may also have other rights that vary from state to state or from one jurisdiction to another.

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