

Digital Basal Thermometer

DMT-3018

User Manual

- EN -

Intended Use:

The Digital Basal Thermometer (DMT-3018) is used by women of childbearing age for the measurement of basal body temperature (BBT).

What is basal body temperature, and why do you need a digital basal thermometer?

Basal body temperature is your temperature when you are fully at rest. It is taken immediately upon awaking before any activity. Unlike traditional thermometers, basal body thermometers can measure with increased accuracy to easily detect your fertile window. The DMT-3018 Digital Basal Thermometer is sensitive enough to measure subtle changes in body temperature with the high accuracy of $\pm 0.05^{\circ}\text{C}$ ($\pm 0.10^{\circ}\text{F}$), and only takes about 60 seconds for the reading with proper use.

Note: It is important to record your temperature every day to identify your BBT trend for ovulation prediction.

How does your BBT help you pinpoint ovulation? What do you need to know?

Your BBT temps will likely surge right after ovulation, with a 0.20°C - 0.60°C (0.40°F - 1.00°F) or more shift.

What's Included:

- 1 Thermometer
- 1 Manual
- 1 Protective case

Storage and transportation condition:

-20°C ~ 55°C (-4°F ~ 131°F), 15%~95%RH, 70kPa~106kPa
Ingress Protection Rating: IP 27
Classification: Type BF (A)
Contraindication: No contraindication

Explanation of Symbols

---	Direct Current		Manufacturing date
	Type BF Applied Part	REF	Catalog Number
	Refer to instruction manual		Manufacturer
	Storage and Transportation Humidity limitation: 15% ~ 95% RH		Storage and Transportation Temperature Limit: -20°C ~ 55°C (-4°F ~ 131°F)
LOT	Batch Code		Recyclable
	Authorized representative in the European Community		
	CE 0197 The product is in compliance with the requirements of MDD 93/42/EEC, "0197" is the identification number of notified body.		
	WEEE (Waste Electrical and Electronic Equipment) Disposal of this product and used batteries should be carried out in accordance with the national regulations for the disposal of electronic product.		
	Ingress Protection Rating The first num.2: Protected against solid foreign objects of 12.5mm ø and greater. The second num.7: Protected against the effects of temporary immersion in water.		

Electromagnetic Compatibility Information

The device satisfies the EMC requirements of the international standard IEC 60601-1-2. The requirements are satisfied under the conditions described in the table below. The device is an electrical medical product and is subject to special precautionary measures with regard to EMC which must be published in the instructions for use. Portable and mobile HF communications equipment can affect the device. Use of the unit in conjunction with non-approved accessories can affect the device negatively and alter the electromagnetic compatibility. The device should not be used directly adjacent to or between other electrical equipment.

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
Operating Instructions

1. Disinfect the probe before each use with 75% ethyl alcohol on a damp cloth.
2. Press the ON/OFF button to switch the device on.

Step 1)
Initial Display: all of the symbols appear.

Step 2)
First Use Vs Second Use or Later-Display Differs Here.

Step 3)
Calibration temperature 36.50°C ~ 37.50°C (97.70°F ~ 99.50°F) with flashing $^{\circ}\text{C}$ or $^{\circ}\text{F}$ is displayed in both instances.

Table 1

Electromagnetic Compatibility information - electromagnetic emission	
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.	
Emissions test	Compliance
CISPR 11 Radiated emission	Group 1 Class B
Harmonic emissions IEC 61000-3-2	Not applicable
Voltage fluctuations/flicker emissions IEC 61000-3-3	Not applicable

Table 2

IMMUNITY test	IEC 60601 test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
Electrostatic transient/burst IEC 61000-4-4	± 2 kV for power supply lines 100 kHz repetition frequency ± 1 kV for input/output lines	N/A
Surge IEC 61000-4-5	± 0.5 kV, ± 0.1 kV differential mode line-line	N/A
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0%UT (100 % dp in UT) for 0.5 cycle at 0°; 40°; 90°; 135°; 180°; 225°; 270° and 315° 0% UT (100 % dp in UT) for 1 cycle at 0° 70%UT (30 % dp in UT) for 2500 cycles at 0° 0% UT (100 % dp in UT) for 25000 cycles at 0°	N/A
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m; 50Hz	30 A/m; 50Hz
NOTE: UT is the a. c. mains voltage prior to application of the test level.		

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Step 4)
When the $^{\circ}\text{C}$ or $^{\circ}\text{F}$ symbols flashes with **Lo** sign, the temperature is ready for use. If ambient or baseline temperature is over 32.00°C (90.00°F), that real temperature will be displayed instead without a flashing $^{\circ}\text{C}$ or $^{\circ}\text{F}$. In both cases, you can now measure.

Step 5)
Temperature measurement will commence and the $^{\circ}\text{C}$ or $^{\circ}\text{F}$ symbol will flash.

Step 6)
The $^{\circ}\text{C}$ or $^{\circ}\text{F}$ symbol will stop flashing and stay on when the reading is ready.

3. Place the thermometer tip under the tongue and close your mouth.

The thermometer tip should remain in constant contact with the tissue under the tongue, and the mouth should remain closed during measuring. The average measuring time is 60 seconds.

NOTE: Avoid eating or drinking just before measuring the temperature, as this will result in an inaccurate temperature reading.

Table 3

IMMUNITY test	IEC 60601 test level	Compliance level
Conducted RF IEC 61000-4-6	3Vrms 150kHz to 80 MHz 6 Vrms 150 kHz to 80 MHz outside ISM bands	N/A
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz	10V/M
Proximity magnetic fields IEC 61000-4-39	30kHz,8 A/m, CW 134.2kHz, 65A/m Pulse modulation 2,1 kHz 13,56 kHz, 7,5A/m, Pulse modulation 50kHz	30kHz,8A/m, CW 134.2kHz, 65A/m Pulse modulation 2,1 kHz 13,56 kHz, 7,5A/m, Pulse modulation 50kHz

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a The ISM(industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz,13,553 MHz to 13,567 MHz, 26,957 MHz to 27,283 MHz, and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz,7 MHz to 7,3 MHz,10,1 MHz to 10,15 MHz,14 MHz to 14,2 MHz,18,07 MHz to 18,17 MHz, 21,0MHz to 21,4MHz,24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.
b The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2,7 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 has been incorporated into the formulae used in calculating the recommended separation distance for transmitters in these frequency ranges.
c Field strengths from fixed transmitters, such as base stations for radio(cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.
d Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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4. Consecutive beeps for approximately 10 seconds indicate that the measurement is complete.

Also, the $^{\circ}\text{C}$ or $^{\circ}\text{F}$ symbol will stop flashing when the measurement is completed. The displayed temperature can increase slightly if the measurement continues after the beep. So it is recommended to keep the thermometric sensor in mouth approximately 2 minutes regardless of beep sound. The highest measured temperature from that reading will appear on the LCD display. If the measurement is over 37.80°C (100.04°F), the beeping will sound with a more rapid consistent buzzer sound to indicate fever. The thermometer will automatically shut off 10 minutes after the measurement is complete or you can manually shut the thermometer off with the ON/OFF button.

NOTE: Make sure to wait until the minimum measurement time is reached before reading your temperature. This is when you hear the last beep of the consecutive beeps. This is to ensure accuracy.

How do you switch between Celsius and Fahrenheit?

Turn the thermometer off. Press the ON/OFF button for 2 seconds to switch the measurement unit between $^{\circ}\text{F}$ and $^{\circ}\text{C}$.

Memory Recall

The digital basal thermometer can only store the last measured value. The memory temperature reading may not save properly if following situations have occurred:

- Switched the measurement unit between $^{\circ}\text{C}$ or $^{\circ}\text{F}$ recently.
- The battery needs replacement.

Operating Mode

Direct mode (not predictive)

Table 4

Recommended separation distances between RF wireless communications equipment				
The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between the wireless communications equipment and the device as recommended below, according to the maximum output power of the communications equipment.				
Frequency MHz	Maximum Power W	Distance	IEC 60601 Test Level	Compliance Level
385	1,8	0,3	27	27
450	2	0,3	28	28
710				
745	0,2	0,3	9	9
780				
810				
870		0,3	28	28
930				
1720				
1845	2	0,3	28	28
1970				
2450	2	0,3	28	28
5240				
5900	0,2	0,3	9	9
5785				
Note 1: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.				

WARNINGS!

- This device should not be used in the vicinity or on the top of other electronic equipment such as cell phone, transceiver or radio control products. If you have to do so, the device should be observed to verify normal operation.
- The use of accessories and power cord other than those specified, with the exception of cables sold by the manufacturer of the equipment or system as replacement parts for internal components, may result in increased emissions or decreased immunity of the equipment or system.
- PORTABLE RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30cm (12 inches) to any part of the ME EQUIPMENT or ME SYSTEM, including cables specified by the MANUFACTURER. Otherwise, degradation of the performance of this equipment could result.

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