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(EN) USER INSTRUCTIONS: Shoe Heel Inserts

Dear Customer, thank you for purchasing shoe heel inserts!

Please read the following instructions carefully to use the product correctly and enjoy its benefits.

It is recommended to keep this manual for future reference.

IMPORTANT: KEEP FOR FUTURE REFERENCE!

Read the entire manual before starting to use the product and follow the recommendations provided.

Do not use the product if any part is missing or damaged.

Store the manual for reference when needed.

1. Preparation:

Before applying the heel inserts, thoroughly clean the inside of the shoes, especially around the heel area, to ensure proper adhesion and hygiene.

2. Placing the heel inserts:

Remove the protective film from the adhesive layer of the heel insert (if applicable).

Place the insert in the back part of the shoe, ensuring it adheres well to the heel area. Make sure the insert is laid evenly and does not create bumps.

3. Adjustment:

Put on the shoes and check if the heel inserts provide the expected comfort. If discomfort or slippage occurs, adjust their position.

4. Care and replacement:

Regularly inspect the condition of the inserts. Replace them if they show wear, lose elasticity, or no longer provide cushioning.

Gel or silicone inserts may require cleaning. Use a damp cloth and mild detergent, then dry thoroughly before reuse.

5. Notes:

Heel inserts are designed for use in closed-heel shoes.

If you experience foot issues or pain, consult a doctor or orthopedist before using the inserts.

Using heel inserts can significantly improve walking comfort, absorb shocks, and correct minor posture problems. Always ensure proper placement and regularly check their condition.

WARNINGS!

- **Do not use with damaged skin:** Avoid using inserts if you have open wounds, burns, or irritation in the heel area.
- **Do not use in inappropriate footwear:** The inserts are intended for use in closed-heel shoes. Avoid using them in sandals, flip-flops, or other open footwear types as they may slip and cause discomfort.
- **Do not use if allergic to materials:** If you are allergic to silicone, gel, or other materials used in the inserts, do not use them to avoid allergic reactions.
- **Do not use in poorly fitting footwear:** Using inserts in shoes that are too tight or too loose may cause increased friction, discomfort, or foot injuries.
- **Do not modify the product:** Do not cut, split, or alter the inserts. Modifying their shape or structure may affect their function and safety.
- **Do not use if pain occurs:** If you experience pain or discomfort while using the inserts, discontinue use immediately and consult a doctor or specialist.
- **Do not use if visibly damaged:** Do not use inserts if they are cracked, torn, or have damaged adhesive surfaces, as they may shift in the shoe and cause injuries.
- **Do not store in unsuitable conditions:** Avoid storing inserts in high temperatures (e.g., sunlight, cars) or in high humidity, as this may degrade their performance.

- **Do not share with others:** Inserts are personal items. Do not share them to avoid transferring bacteria or fungi.

Thank you for choosing our product!

We wish you comfort and satisfaction during use. 😊

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