User Manual

Aster Replacement Electric Toothbrush Heads

Installation of Electric Toothbrush Head

Make sure the electric toothbrush is turned off.

Unplug the old toothbrush head and insert the new toothbrush head into the connection port to ensure a firm connection.

Replacement of Electric Toothbrush Head

It is recommended to replace the toothbrush head every 3 months, or replace it in time according to personal usage.

Unplug the old toothbrush head and insert the new toothbrush head.

Precautions for using electric toothbrush head

Before using the toothbrush head, read the product manual carefully.

During use, please avoid brushing your teeth hard, just let the electric toothbrush head slide gently on the surface of the teeth.

Please clean the toothbrush head in time after use.

Avoid placing the electric toothbrush head in a hot or humid environment to avoid affecting its service life.

Safety Warnings:

Avoid accidental ingestion: When replacing the toothbrush head, keep small parts out of reach of children to prevent children from accidentally ingesting or misusing it, causing suffocation or other dangers.

Children's use: Children should be supervised by adults when replacing electric toothbrush heads to ensure safe operation.

Matching model: Make sure the replacement toothbrush head you purchase matches the model of your electric toothbrush. Do not use a toothbrush head that does not match the model, so as not to affect the use effect or safety.

Regular inspection: Regularly check the wear of the toothbrush head. If it is found to be severely worn or deformed, please replace it in time to avoid damage to the oral cavity.

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