

Steaming Steps

1. Connect the air inlet of the barrel and the air outlet of the fumigant with the rubber tube.
2. Put the medicine can into the barrel and connect it with the air inlet L of the barrel with the rubber tube, and then put the massage stool in the proper position in the barrel.
3. Add proper amount of water to the fumigant, and then turn on the main power switch of the fumigant and the power switch on the control panel. The fumigant will generate steam in about 5 minutes. The specific steam generation time will vary according to the ambient temperature and gear changes.
4. Sit on the stool and put your feet on the massage stool in the bucket naturally and comfortably. And cover the barrel mouth with towel.
5. Use the key on the remote control or fumigation instrument control panel to control and adjust the time and temperature in real time.

Panel Functions

1. Power button (1 or A): press this button to control whether the machine is in standby or working state.
2. Timing decrease button (2 or) : press this button to control the decrease of timing time of the machine.
3. Timing increase button (3 or B) : press this button to control the increase of the machine's timing time.
4. Gear decrease button (4 or D): Press this button to control the decrease of the machine gear.
5. Gear increase button (5 or E) : press this button to control the increase of the gear of the machine.
6. Safety valve of fumigant (cannot be taken out).

Attentions

1. Do not add water for medicine change, inject water and clean when power on the machine. After 5 minutes, open the cover after the machine cools down, and restart the fumigant after operation. Do not open the cover or pull out the steam pipe when working, so as to avoid scalding.
2. If the power cord is to be replaced, in order to avoid danger, it must be replaced by professionals from the manufacturer's maintenance department or similar departments.
3. The rubber tube must be kept unobstructed and shall not be bent or blocked. When the fumigant is dry, the machine will automatically cut off the power. At this time, the machine can not be started until the inner tank of the fumigant is cooled. Generally, it is better to cut off the power manually when the water is not dry, which is good for machine maintenance.
4. When putting the traditional Chinese medicine, it must be wrapped with a sandbag to prevent the medicine from blocking the steam outlet, otherwise it may cause machine failure or safety accident.
5. After using the fumigant, wipe off the water stains, keep it clean, and there should be no scaling. It is strictly forbidden to wash the fumigant shell directly to avoid short circuit caused by moisture. When idle, it should be placed in a cool place to avoid sun exposure and moisture corrosion. Sharp objects are forbidden to collide to prevent steam leakage caused by damaging cylinder block.
6. Do not use the fumigant as other heaters, and it is not allowed to repair or modify the main engine by yourself.
7. The elderly and children should use the machine under the care of their families. It is not suitable for those who are over hungry, who just finish meals 30 minutes ago, who suffer

from serious body failure, acute severe water loss, high fever, massive blood loss, severe hypertension and heart disease, purulent ulcer and wound, female menstrual period, pregnancy, drunkenness and unclear consciousness.

Usage

1. Leakage protection device must be installed in the room, qualified power socket should be used, the load capacity should be above 10A, the voltage should be 110V, and the maximum working current of single fumigant should be 5A.
2. The fumigant and the barrel shall be placed on the same plane, the power switch and the display board of the fumigant shall avoid water splashing, and it is strictly prohibited to use in the toilet and the place with high humidity.
3. Do not add too much water, and the water level shall not exceed the maximum water level, nor be lower than 50% of the water volume in the tank; when heating, adjust to the 9th gear, and then adjust to the temperature that is suitable for the individual after the steam comes out, or directly add hot water for heating.
4. Pay attention to the safety valve of the fumigant. If the safety valve is found to be out of water or out of steam, turn off the power immediately and check whether the steam outlet pipe is blocked. If there is a block, the cover cannot be opened immediately. Only after the machine is cooled can the cover be opened for inspection.

After Use

1. After the steam health care system of traditional Chinese medicine is finished and the water temperature in the standby room drops, the accumulated water in the machine shall be discharged in time, and the water shall be poured into the fumigant to clean the inner tank.
2. After one month of use, pour "edible citric acid" into the machine, heat it for a while, and remove the scale in the machine.
3. Pay attention to soaking (1 / 4 water level of barrel) before and after use to prevent water leakage and dry crack of wood.
4. It is better to use it at least twice a week to keep moisture in the barrel.
5. If it is not used for a long time, it shall be placed in a wooden barrel with a wet sponge and sealed with a plastic bag.

Package List

1. 1 steam fumigant
2. 1 ultrathin remote controller
3. 1 set of steaming medicine cans
4. 1 silicone steam pipe
5. 1 user manual