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1. Die Stromversorgung

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Instructions for use

1. the power supply

The device has an integrated rechargeable battery with a capacity of 300mAh. Use a 5V 1A USB adapter (e.g. from a cell phone) and a USB-C cable (included) to charge. The charging time is about 2 hours. With a fully charged battery, the device can work for about 1.5 hours with white light and 2 hours with blue light. It is important to check the battery charge before using the glasses.

Briefly press the power button to see whether the power indicator lights up green (battery is sufficiently charged) or red (battery is low and the glasses can be used for a maximum of 15 minutes).

If the power indicator flashes red and an alarm sounds, this means that the glasses are almost empty and need to be recharged.

2. adjustment of the wearing comfort:

The glasses can be adjusted for a comfortable fit. This includes adjusting the temples inwards or outwards to adjust the angle.

For smaller head circles, there is a non-slip strap for a comfortable fit.

It is important to set the correct height for wearing the glasses so that they do not interfere with daily activities such as reading and office work. For this purpose, there are three holes for the nose pad, between which the height of the glasses can be freely adjusted.

If you need to use the Blue-3 Pro glasses together with normal dioptric glasses, you can wear both at the same time, but make sure that the dioptric glasses do not block the light from the light therapy glasses.

Do not use the light therapy glasses with your eyes closed.

For safety reasons, it is not recommended to use the light therapy glasses while charging.

3. operation of the glasses

The device is switched on by pressing and holding the button for 3 seconds, after which the light switches on.

The blue/white light can be easily switched with a button within 5 seconds of switching on. The light module flashes three times when the button is pressed to confirm the selected light color.

Once you have set the light to the desired color, the device prepares for the therapy. You can set a timer and put on your glasses to start the light therapy.

The light is gradually brightened up to maximum brightness within 60 seconds so that your eyes can get used to it.

The light therapy stops automatically after the set time has elapsed.

4. set the duration of use:

The duration of use of the therapeutic goggles can be set to 15, 30 and 45 minutes, depending on the user's current situation.

Press the power button to select the duration of use. The number of displays indicates the duration of the therapy (one display = 15 minutes, two = 30 minutes, three = 45 minutes).

After the set time has elapsed, the device warns you with an acoustic signal and switches off automatically. You can resume therapy by pressing and holding the on/off button.

5. how to use

Use it to stabilize the circadian rhythm: Use the blue light every morning at around the same time for at least 15 minutes. If you need to adjust your circadian rhythm to get up later and go to bed later, use the device about 3-4 hours before going to bed.

To manage the negative effects of jet lag, you should use glasses so that your body can adapt more quickly to the time change. For example, if you are traveling west, it is advisable to use glasses in the afternoon so that you are tired later and go to sleep in the evening. It is then advisable to use the glasses in the morning after waking up in order to anchor the new sleep rhythm more quickly.

To reduce seasonal (winter) depression/waves, use the device with blue light directly after waking up for 15 - 30 minutes and with white light during the day/afternoon for 30-45 minutes.

CAUTION, do not use the device 1.5 hours before bedtime as this would reduce melatonin production and worsen sleep. We recommend wearing orange/red glasses within 1.5 hours before bedtime to block blue and green light.

6. precautions for use:

Please note that the eyes should be open when using the glasses so that the light can act on the photoreceptors in order to achieve the expected function of the light therapy.

The blue light emitted by the glasses is of low intensity and complies with photobiological safety standards.

7. warning

Keep the appliance out of the reach of children.

Use the device in a dry environment, the device is not waterproof.

If the device makes suspicious noises or behaves unusually, do not use it and return it.

Never repair or dismantle the appliance yourself.

People with health complications such as diabetes or various eye diseases should first consult a doctor about the use of light therapy.

Do not use this product if you suffer from eye diseases such as cataracts, glaucoma, eye inflammation, optic neuritis, etc.

In rare cases, wearing glasses for light therapy can cause headaches, eye strain or nausea. If you notice these problems, discontinue light therapy.

INSTRUCTIONS AND INFORMATION ON HANDLING USED PACKAGING
Dispose of the used packaging material at a place designated by the municipality for waste disposal.

DISPOSAL OF USED ELECTRICAL AND ELECTRONIC EQUIPMENT
This symbol on the products or in the accompanying documents means that used electrical and electronic equipment must not be disposed of with normal household waste. Take these products to the designated collection points for proper disposal, reprocessing and recycling. In some countries of the European Union or other European countries, you can also return your products to your local retailer if you purchase an equivalent new product. By disposing of this product properly, you will help to conserve valuable natural resources and prevent possible negative effects on the environment and human health that may result from improper waste disposal.

This product fulfills all essential requirements of the applicable EU directives.

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