

Infinite Relax

Smart Cupping Therapy Set

User Manual



01 How to better use InfiniteRelax cupping set?

Please refer to the instructional video for a better understanding of how to use the cupping device :

https://www.amazon.com/live/video/0a23fd9827174ac295c747644d2da9a3?ref_=dp_vse_ibvc1

02 Why do you feel it has no suction or weak suction?

The cupping device has 12 levels of suction, ranging from weak to strong. The sufficient suction should meet your needs.

If you feel that there is no suction or weak suction, please refer to the instructional video or the images below.

1. If you don't align and tighten the attachments properly when installing them, it can result in the cupping device lacking suction.



2. When leaving the cupping cups on the skin, please slowly lift the cupping device in a vertical motion. Incorrect technique can lead to unsuccessful cupping.



03 InfiniteRelax Cupping Set Introduction

Brief Control Panel



- Warm Tips:**
1. Please use this equipment correctly according to the use manual
 2. It is normal to have red marks after cupping, it will subside naturally within 3-7 days.
 3. First time users should choose low suction
 4. Do not use on one area for more than 15 minutes(first time users less than 3 minutes)
 5. Do not let essential oils or water get inside the machine

6 cupping cups with massaging, safe, easy to use, and dynamic suction. It has a 20-minute intelligent timing stop, and the unique one-touch pressure relief mechanism instantly releases the cup pressure and eliminates the pain during cup removal.

Cupping Device effectively reduces tension, pain, cellulite, and muscle soreness in targeted areas, providing quick relief and relaxation.

Brand	InfiniteRelax
Product Size	5.70inch*3.9inch*1.8inch
Charging Voltage	5V 2A
Suction Strength	-60~0kpa
Rated Power	2.6W

04 Preparation Before Use

1. Please read the safety precautions carefully when using this product.
2. Check whether the instrument is complete and normal.
3. It is recommended to clean and moisturize the skin before using.
4. Tighten and install the cupping cup.

05 Operational Guidelines



Power/Mode Key

1. Press and hold for 1.5 seconds to start the machine, and The LCD screen on the control panel will display the 20 minute countdown.
(It is important to only start with a few minutes as a beginner to test skin sensitivity)
2. Press and hold for 1.5 seconds to shutdown, and the LCD screen will go out.
3. In the power on state, short press to change the mode. The machine has two modes, deep negative pressure mode and soothing mode Breathing pattern.



Suction Strength Key

1. The power on defaults to the third gear strength.
2. Short press to adjust the suction strength (1-4 gear, low suction. 5-8 gear, mid suction. 9-12 gear, strong suction.)
3. Each short press will increase 1 gear upward, and the LCD will display the corresponding strength gear.
4. It is recommended to use it from gear 1 for the first time or when the skin is not adapted.



One Click Pause Release Button

During cupping, if you want to stop or remove the product briefly, press this key to pause cupping. The cupping tank adsorption function and pressure relief can easily lift the tank and enter the standby state.



Red Light Irradiation Key

Short press this button to turn on or off red light irradiation.



Operation Steps of One Main Tank

1. take out the massage instrument, inspect the equipment and wipe it clean. Select acupoints and clean the skin according to massage needs.
2. Press the power button for 1.5s to turn on the device, and select the appropriate gear through suction strength button.
3. If you need red light therapy, you can turn on the red light irradiation button.
4. After cupping, use the pause/release button to release the pressure, remove the massage instrument, and long press the power button to shutdown.



Operation Steps of Five Supporting Can

1. take out the massage instrument, inspect the equipment and wipe it clean. Select acupoints and clean the skin according to massage needs.
2. Screw the conversion tank into the machine, and then insert the supporting can into the conversion tank, so that the supporting can and the conversion tank are assembled together on the machine.
3. Press the power button for 1.5s to turn on the device, and select the appropriate gear through suction strength button.
4. Suck on the skin for 10s, the skin feels tight, pull the machine directly upwards, leaving the supporting can on the skin surface of the body.
5. Short press the red and blue light button on the supporting can, to turn on the red or blue light irradiation.
6. After cupping, pull the pressure relief value on the supporting can to remove it, and long press the power button to turn it off.

06 Charging Instructions

1. Connect the type-C charging cable to the adapter socket, and connect the other end to the product charging interface to start charging Electricity. Adapter size 5V=2A.
2. Please keep the type-C charging port dry and clean during charging. The charging time will vary slightly due to environmental conditions, Subject To the actual time. Do not charge for more than 8 hours.
3. During charging, the indicator flashes and goes out after being fully charged.
4. The product cannot be used when it is charged.

07 Safety Precautions

DANGER(contents related to serious personal injury caused by wrong operation)
Do not use with the following medical electronic equipment at the same time, which will cause errors or damage to medical instruments, Seriously life-threatening!

- 1.Implantable medical electronic devices such as cardiac pacemaker;
- 2.Life sustaining medical electronic devices such as artificial heart and lung;
- 3.Electrocardiograph and other medical electronic instruments;
- 4.Patients with other medical implants.

Please do not use under the following circumstances to avoid accidents or discomfort.

- 1.Patients with pregnancy, menstruation, damaged epidermis, varicose veins and taking drugs.
- 2.Patients with infectious diseases, skin sensitivities,postoperative patients, high fever or critical illness are prohibited.
- 3.Patients with malignant tumor or local damage and ulceration, traumatic fracture, acute suppurative inflammation, skin Patients with loss of elasticity.
- 4.Patients with abnormal heart, brain nerve and blood pressure.
- 5.Severe neuroticism, convulsion, spasm and manic anxiety.
- 6.Bleeding prone diseases, such as thrombocytopenia, leukemia, hemophilia, Henoch Schonlein urpura and other diseases.
- 7.Near the heart, facial features, private parts or large blood vessels on the body surface.
- 8.Patients with skin perception disorder or abnormal skin.
- 9.A patient who is being treated by a doctor or feels physically abnormal.
- 10.Excessive fatigue, drunkenness and staying up too late (2 hours before and after meals are also prohibited).
- 11.Patients whose skin is allergic to metal and red light.
- 12.The product has a heating surface, which is not sensitive to heat, children,unable to express their wishes freely or unable to self Management personnel.

08 Cleaning and Maintenance

Please observe the cleaning and maintenance methods to prolong the service life of the instrument:

1. Please wipe the fuselage with dry non-woven cloth. Do not use acid-base detergent, alcohol, diluent, gasoline. Clean with solvent such as nail

remover.

2. After each use, please turn off the power supply and clean it (with wet towel). 3. Please store this product in a safe and dry place.

4. Do not leave the product idle when there is water residue.

5. Do not use a hard brush or stainless steel ball to wipe, so as not to damage the surface.

6. Do not use water and other liquids to directly inhale the instrument for flushing.

7. When the product is idle for a long time, please charge the product regularly to ensure its normal operation.

Do not use or store this product in the following places:

Places vulnerable to vibration and impact; Places with direct sunlight; Hot and humid places; Places with water; Strong electromagnetic field; Near fireworks; Places accessible to children, etc.

09 Fault Analysis and Troubleshooting

Fault phenomenon	Fault judgment and solution
Unable to power on	1.Press and hold the power button for more than 1.5 seconds to start up. 2.Please charge the machine and try again.
Unable to charge	1.Replace the adapter with a row socket and try again. 2.Please pull out the charging head and plug it in again and confirm that it is inserted in place. 3.Please check the charging head and charging port for dirt or sundries. If yes, please clean up. 4.Replace the adapter or charging cable.
Weak suction	1.Remove the transparent tank and re tighten it for installation. 2.Check whether the air inlet and exhaust holes are blocked. 3.Check whether the filter cotton needs to be replaced.

If the problem cannot be solved,do not disassemble the machine without permission.

Reminder:

- 1.When using it for the first time, please pay attention to step-by-step adjustment to find a suitable mode and intensity.
- 2.In the deep negative pressure mode, if the suction reaches a certain degree, please press the pressure relief /pause key in time to stop working or relieve pressure to avoid skin discomfort.

10 FREQUENTLY ASKED QUESTIONS

1. Is cupping good for tight back muscles?

yes, cupping is perfect for tight back, neck shoulder, and many other large muscle groups.

2. Does cupping get rid of muscle knots?

Absolutely! Cupping works wonders on releasing muscle knots and tension.

3. How long does it take to feel the benefits of cupping therapy?

Most people begin to feel improvements after using it just once, but some take two or three sessions to experience real results. cupping therapy effects tend to become more realized over time.

4. How do I use the suction cupping for cellulite removal?

It is recommended to apply a light coat of essential oils on the skin. Featuring dynamic suctioning, adjustable pressurized cupping, and red blue light therapy, it is easy and effective to reduce cellulite.

5. What does cupping do to your body?

cupping therapy increases blood circulation to the area where the cups are applied. This relieves muscle tension and improves overall blood flow to the area to stimulate cell repair.

6. Are they easy to use at home, by yourself?

our devices are intuitive right out the box and very easy to apply to enjoy a healthy therapy session from home.

7. Does it leave a mark?

Everyone skin is different slight pain with high levels of pressure is common and so are the classic red circles that come with cupping though this is temporary, these marks will gradually disappear within a few days and they are nothing to worry about!

8. Can the same device be used on multiple people? Are they easily sanitized?

Sure, our devices are easy to sanitize by wiping down or hand washing the Cup with your preferred disinfectant.

9. Is there any scientific evidence on cupping reducing cellulite?

yes! There are many scientifically-proven studies which have proved significant decrease in female cellulite and improved skin elasticity from a direct result of cupping therapy.