




Prepare

Use the razor provided along with the device to remove all visible hair on the skin before using the device.



Connect the power

Plug the device while using it to do hair removal.



Goggles

Wear the goggles before the treatment starts.



Get a comfy set up

Using mirror can help you see better when you zap your Brazilian wax area, armpits and legs.

*Please refer to the user manual for safety precautions before use.

Ulike Support & Warranty:
amazonsupport@ulike.com

CONSISTENCY IS KEY TO EFFECTIVE HAIR REMOVAL RESULTS

ARMPIT

3/4/5
Spot Method (2-3 flashes/area of skin)

* Don't use Ulike device

- ✗ near your eyes or ears
- ✗ on damaged skin e.g. sunburns, cuts, infections, cold sores etc
- ✗ on implants like pacemaker, neurostimulator, insulin pump or silicone implants
- ✗ on the testicles, genitals and nipples
- ✗ on warts, moles or large freckles
- ✗ on tattoos or permanent make-up
- ✗ if you have epilepsy with flashlight sensitivity or light sensitivity
- ✗ if you are on medicines that make you sensitive to light e.g. Retin A, isotretinoin Accutane or other topical retinoids

LEG

3/4/5
Spot Method (2-3 flashes/area of skin)

FACE

3/4/5
Spot Method (2-3 flashes/area of skin)

ARM

3/4/5
Spot Method (2-3 flashes/area of skin)

BRAZILIAN WAX AREA*

3/4/5
Spot Method (2-3 flashes/area of skin)

*If you want to be more comfortable in the brazilian area, it is recommended to turn on the device for 1-2 mins before applying it on skin.

*After the initial treatments, you may notice the emergence of pimple on your skin. This is a common histamine response and is entirely normal. With continued use, the pimple will gradually subside.

Treatment Schedule

Weeks 1~4: treat with the device 3 times a week. For better results during weeks 1~4, we recommend treatments **EVERY 2 DAYS**.

From week 5 on: you can adjust the treatment schedule to once every two weeks or once a month, depending on your hair growth.

Intensity: 5 > 4 > 3 > 2 > 1

The mode are just suggestions. Use the highest mode that feels comfortable for you. Higher modes can give you faster and better results. Start with the lowest mode and work your way up. You can use different modes for different areas. Your comfort is important.