
3-IN-1 LIGHT THERAPY LAMP

HAPPY LAMP & SUN LMAP & SAD LAMP

Sun Lamp

5 Color
Temperature

Night Light

5 Brightness
Adjustment

Detachable

Timing
Function

Remote
Control

WHAT'S INSIDE?


Thank you for purchasing Light Therapy Lamp!
Bringing the sunlight into your home or office relieves the winter blues, boosts mood, increases energy, and makes you feel revitalized. Please read the user manual carefully before use, and check whether the product is damaged or not when you receive the product.

1×Lamp Cap


1×I-Beam Clip

1×Lamp Hose

1×Threaded Connector

1×Adapter

1×User Manual

1×Remote Control

PRODUCT SPECIFICATIONS

Lamp Type	Light Therapy Lamp
Product Model	GLD-YP5C-G2V3
Product Color	White
Texture of Material	ABS+Aluminumtube+Silicone Gel
Control Method	Touch Control+ Remote Control
Remote Range	5M/16.4FT
Timing Function	60S/10MIN
Number of Lamp Beads	108PCS-3000K, 108PCS-6000K
Color Rendering Index	>90Ra
Light Source Type	Direct Illumination(SMD LED)
Light Source Power	18W
Light Intensity	10000LUX
Brightness Levels	20%/40%/60%/80%/100%
Color Temperature	3000K/3500K/4200K/5000K/6000K
Table Thickness for Installation	≤5.8CM

HOW TO WORKS?


Long-term lack of sunlight or healthy light can negatively affect our health, gradually leading to problems such as low mood, seasonal affective disorder (SAD), winter blues and sleep disturbance caused by circadian rhythm disturbances.

Therefore, we require a device that can consistently provide bright and ample light to meet our needs.

HAPPY LAMP adjusts our brain chemistry and biological clock by mimicking the sun's rays, increasing the frequency of photostimulatory hormones and neurotransmitters. While boosting energy levels, it greatly enhances positive emotions and well-being, thereby alleviating symptoms caused by lack of sunlight.

INSTALLATION INSTRUCTION


Therapy Light Mode

1. Thread the hose (B) into the I-Beam Clip(A).
2. Clip light head (D) onto hose (B) until you hear the "clicking" sound.
3. Insert the adapter (AA) into the hole on the back of the control panel and the electrical outlet.

- Please hold the lamp post to prevent it from swinging while adjusting the angle.
- Do not over-adjust the lamp neck to avoid damage to the connector parts.


PRODUCT INFORMATION

ICON GUIDE

Icon Legend					
Uses	Night Light Mode Button	5 Steps Brightness Button	5 Color Temperature Button	Timing Mord Button	Turn ON/OFF Button
Function	Enter Night Light Mode	Adjust 5 Brightness Levels	Adjust 5 Color Temperature	 ·The entry into the 60S timing mode. ·The light will automatically turn off after the time expires.  ·The entry into the 10min timing mode. ·The light will automatically turn off after the time expires.	 Button shows a white light, it indicates that the lamp is not in use now and the adapter is functioning properly.  Button shows a blue light, it indicates that the lamp is in use now and the adapter is functioning properly.
Special Instructions		When afted insert the Adapter (AA) into the electrical outlet, button didn't have white or bule light Please check whether the adapter is connected successfully or replace the outlet If you still can't use it after checking, please contact us in time support@sunlamlux.com			
About Adapter	If you find that the adapter is damaged, you can buy our original adapter or buy it yourself according to the following information				
	Model	Input	Output	Constant Voltage Output	
	LDGLD-1300900U	100-240V-50/60HZ 0.35A Max	13.5V 1.77A	SELV 23.8W	

REMOTE CONTROL INFORMATION

- The remote control range is within 10M/32FT, and you need to **AIM** the light therapy lamp.
- Batteries are not included in the remote, you need to have **AA batteries**.

RECOMMENDATIONS FOR USE

1. It is recommended that you adjust the light to the darkest when you use it for the first time, and gradually get used to it.
2. Pay attention to personal feelings, if any discomfort or adverse reactions occur, please reduce the use time or consult a doctor for advice.
3. Do not look directly at the light for a long time, do not tilt the light.
4. Try to place the light on your side or peripheral viewing angle, and adjust the angle that is comfortable for you. You can start your normal activities while using it.
5. The most effective time to use it is when you wake up in the morning, it is like sunlight entering your eyes at this time, and you can use it as needed at other times.

If you are pregnant, please seek medical advice before using happy lamps.

EFFECTS OF DIFFERENT LIGHT COLORS

Color Temperature	Therapeutic Effect
Warm Light	Warm light typically has a lower color temperature, giving a warm and cozy feeling. It has positive effects on relaxation, stress relief, and anxiety reduction. Therefore, warm light may be helpful for conditions such as mood disorders, anxiety, and insomnia.
Cool Light	Cool light simulates natural light and contains a range of wavelengths. White light can improve alertness, enhance attention, and cognitive function. It is often used in the treatment of seasonal affective disorder (SAD) to help regulate the circadian rhythm and improve mood.
White Warm Light Blend (Nature Light)	The Nature light blend provides a warm and gentle illumination, offering both soothing effects and promoting alertness and focus. It can enhance alertness, boost mental energy, and improve mood.
ATTENTION	Each individual may have a different response to light therapy. The therapeutic effects vary from person to person and are influenced by individual physiological and psychological differences. Therefore, specific treatment outcomes may vary based on personal conditions, methods of use, and other factors during the treatment period. It is recommended to consult a professional healthcare provider or light therapy specialist before using light therapy to obtain personalized advice and guidance.

FREQUENTLY QUESTIONS

Q: WHAT ARE THE EFFECTS OF USING HAPPY LAMP IN THE MORNING AND EVENING?

A: DAY: By stimulating the body to produce serotonin to relieve anxiety, improve mood and well-being.

NIGHT: Regulating melatonin to promote sleep.

Q: DOES THE EFFECT OF THE THREE COLOR TEMPERATURES HELP IN THERAPY?

A: All three color temperatures of light have a certain therapeutic effect. Among them, the effect of white color temperature is the best.

Q: WHAT IS THE LIGHT INTENSITY FOR EACH LEVEL?

A: Under the 25-30CM distance test, the light intensity of the three lights are respectively

Level 1: 2000 LUX

Level 2: 4000 LUX

Level 3: 6000 LUX

Level 4: 8000 LUX

Level 5: 10000 LUX

Q: HOW FAR DOES IT TAKE TO ACHIEVE THE HEALING EFFECT OF 10000 LUX?

A: Recommended distance is about 25-30CM.

Q: CAN IT BE USED AS A READING LIGHT?

A: Yes. This therapy lamp has 5 brightness settings and movable lamp head. You can adjust the brightness and tilt head down to use as an everyday lamp.

Q: IS LIGHT THERAPY SAFE FOR EYES?

A: Phototherapy lamps without UV and blue light are safe for eyes.

SAFETY TIPS & INFORMATION

Before starting to operate the product make sure that all parts are ready, do not attempt to install or run the product if any parts are missing or damaged. Contact customer service for replacement parts.

1. Do not attempt to replace the plug.
2. Do not use extension cords.
3. Do not look directly at the LED lamp to avoid eye damage.
4. Do not cover the product with cloth, paper, or combustible materials.
5. Do not put your hand on the LED light, because it will become hot when you use it.
6. If not in use for a long time, unplug the adapter.
7. Place on a stable surface.
8. Please use indoors only.
9. If the light is off or flashing, check that the adapter and the DC connection socket are secure.

CARE AND MAINTENANCE

Package contents ensure that the lamp is unplugged before re-cleaning.

1. Use a rag to clean the device.
2. Do not use detergent or abrasive for cleaning.
3. Do not immerse the device in water.

When the light doesnot work , Please check those

1. Make sure the switch and power supply are on.
2. Please make sure the adapter is connected successfully or replace the outlet.
3. Check that the lamp cap and lamp hose connect successfully.

For proper installation and safety, you must use our supplied parts to reduce the risk of fire, electric shock, or injury.

HOW SHOULD WE USE IT?

1. Try to place the light on your side or peripheral viewing angle, and adjust the angle that is comfortable for you. You can start your normal activities while using it.
2. This table lamp can be used at any time, the most effective time is 6 am to 6 pm, and you can use it according to your needs.
3. Do not stare directly into the light throughout use as this may cause retinal overstimulation, allowing for brief staring into the light.

AFTER-SALES SERVICE

Thank you again for purchasing our products.

We hope the light therapy lamp can bring you health, sunshine, optimism, and happiness. If you have any questions, please send an email to **support@Sunlamlux.com** and we will deal with it in time.

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