

COSLUS

ACTIVATE YOUR **24-MONTH** WARRANTY FOR FREE



* Please contact customer service if you're unable to provide the order number

www.coslus.com/warranty



Win Free Accessories



Covered by 24-Month Warranty



Friendly Responsive Customer Service

Expire time: 7 days after the delivery date


 +1 855 594 3950 (9am-5pm Mon. - Fri. (EST))

 support@coslus.com


FLIP →

QUICK GUIDE


Scan QR code to watch the >> operation video.




Charge Instruction




Install nozzle



Fill water




Adjust water pressure




Adjust mode


Note: Press mode button to select mode.



Put nozzle tip into mouth



Start irrigating



90° to gum line.



Tilt and shake the main unit to drain the water from the speed control knob.