

Nekteck



Shiatsu Foot Massager



Congratulations

Purchasing this massage device reflects your desire for comfort and wellness of your health

In order to benefit from the advantages of this device for a long time, we would ask you to read and observe the safety instructions carefully.

We using hope you'll love and enjoy your personal [Shiatsu Foot Massager](#)

Table of contents

Safety instructions	2
Scope of delivery	3
Specifications	3
Starting up	3
Control panel	4
LED Display	5
Specials	5
Foot reflex zones	6
Troubleshooting	8
Storage Maintenance	8
Warranty conditions	9

Safety instructions




Please read the following instructions carefully before putting your massage device into use to ensure trouble-free operation and optimum efficiency. Please retain these operating instructions for further use!

- The massage device conforms to recognized technical principles and the latest safety regulations.
- Any possible repairs may only be carried out by authorized expert staff. Im proper use and unauthorised repair are prohibited for safety reasons and may lead to loss of warranty.
- Never touch the power plug with wet hands.
- Please avoid contact of the device with water, exposure to high temperature and direct sunlight.
- To avoid the danger of short-circuit, disconnect the power cord when not used for long period of time.
- Do not use any damaged cables, plugs or loose sockets. If the plugs or cords are damaged, they must be replaced by the manufacturer, a service representative or by qualified personnel.
- In case of malfunction, disconnect immediately from the power outlet . Misuse or incor-rect use excludes any liability for damage.
- In order to avoid excessive stimulation or issues with the muscles and nerves, the recom-mended massage time of 30 minutes should not be exceeded.
- Never use pointed or sharp objects.
- please refrain from using the massager or any massage during preg-nancy or if one or more of the following complaints are present in the massage area: recent injuries, thrombotic diseases, all kinds of inflammation and Swelling, and cancer. It is recommended to consult the doctor before taking a massage for the treatment of illness or infirmity.
- Do not use this device while driving.
- If you are dependent on electrical aids e.g. pacemakers, please be sure to consult your doctor for medical advice before taking a massage.
- This appliance can be used, by children aged 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning the use and hazards of the appliance.

- Children are not allowed to play with the device.
- Children without supervision should, not do Cleaning and maintenance.
- Persons who are sensitive to heat should be careful when using the unit.
- Use the device only with the supplied equipment.

Scope of delivery

- Shiatsu Foot Massager
- Main adapter




Specifications


Dimensions:	36.5 x 34.0 x 17.0 cm
Weight:	2.3 kg
Voltage:	Input: AC 100-240V ~ 50/60Hz Output: DC 12V, 2000mA
Nominal power:	24 watt
Runtime:	15 minutes

Putting the device into use

Before starting up, check whether the main voltage matches the voltage specified on the main plug once you have confirmed that the connect the power cable to the massage device and then plug unit matches the required power voltage the power plug into the socket.

Control panel



Button	Function
	Press this button to turn on the device. The unit starts in kneading MASSAGE MODE with a low intensity level
MODE	Once the unit is turned on MASSAGE MODE+ AIR PRESSURE MODE will be in it's default setting (Kneading+ Rolling+ Air Pressure are activated) Press the Mode section (once) and it will change from default to MASSAGE MODE Only (Kneading+ Rolling Only) Press the Mode section (2X) it will activate AIR MODE Only (Air Pressure Activated)
HEAT	You can activate the foot sole heat therapy by pressing the HEAT button on the device.
AIR	By pressing the AIR button on the device one of the 2 intensity levels can be selected.

The automatic massage lasts 15 minutes. The massage device automatically switches off after this period of time in order to prevent overheating. Unplug the power cable from the socket after every use.

LED Display

AIR & MASSAGE MODE	If this function is activated, the LED illuminate "green" on the device.
MASSAGE MODE	If this function is activated, the LED illuminate "yellow" on the device.
AIR-MODE	If this function is activated, the LED illuminate "red" on the device.

FEATURES



COMPRESSION

Rising and receding airbags have a pumping motion which produces natural muscle reflexes. This is particularly beneficial in the calf and foot area.



HEATING

A deep penetrating and heat-radiating effect which provides relaxation and relieves is particularly beneficial for muscle tension due to the deep reaching effect.



REFLEXOLOGY

Reflexology regenerates foot muscles and has a harmonizing effect on almost all organs inside your body.



SOLE ROLLER

The foot massage stimulates the nerve endings and certain organs of the body, which also relaxes the muscles.



HEEL MESSAGE

For a particularly beneficial foot massage, the heel area of the Shiatsu Foot Massager with air compression is equipped with special acupressure heel airbags and heel rollers.

Foot reflex zones

Reflexology has the advantage over an all-body massage due to the fact that there are approx. 72,000 nerve endings of different organs in the feet. Precisely performed massages can be very beneficial to our health and well-being.

Precision massages of the foot reflex zones can stimulate metabolism and blood circulation among other things, it also stimulates digestion and support the self-healing process of the body.

The more we care about the vitality of our feet, the better it is for our body: The feet are known a microsystem in which each organ is said to be connected and reflected and can be treated individually.





- > Neglecting the feet can affect the whole body.
- > Muscle tension can arise due to poor posture of the body. This in turn can have negative effects on the spinal column and the nerve tracts and result in other physical complaints.
- > Sleeping problems and cold feet are usually results circulatory disorders.
- > Swollen legs, muscle cramps and blood congestion can be caused by excessive sitting and physical inactivity.
- > Provide your feet with sufficient energy so that the veins embedded in the muscles can transport the blood throughout the entire body again.

Troubleshooting

FAULT

The device is connected to the power outlet, however neither the device nor the control panel works.

REMEDIAL ACTION

- 1) Check the connection to the adapter.
- 2) The device can overheat during prolonged use. For this reason, you should switch it off and let it cool down for 30-50 minutes.

Storage | Maintenance

With proper usage maintenance you can increase the life span of your massage device and at the same time ensure safe operation at all times.

1. Recommended storage: Dry, free of dust.
2. it away from any moisture, liquid or combustible materials and open flame from the device.
3. Clean the device with a soft, dry cloth.
4. Unplug the device from the power outlet when not used

Warranty conditions

We provide all our customer 18 months warranty from the date of its original purchase. if the device fails due to manufacture defects, please contact us by email: support@nekteck.com to launch a warranty claim, we will instruct you on how to return the product or repair /replacement. Please do not return the product without seller/authorization.

Nekteck

Email: support@nekteck.com

WEB: www.nekteck.com

32520

Nekteck

The background features a teal color palette with white, flowing, organic shapes that create a sense of movement and depth. The shapes are layered, with some appearing to flow from the top right towards the bottom left, while others curve in the opposite direction. The overall effect is clean, modern, and minimalist.