BODIFY®

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USER MANUAL • MANUEL D'UTILISATION • MANUALE D'USO • MANUAL DEL
USUARIO • GEBRUIKERSHANDLEIDING • ANVÄNDARMANUAL • INSTRUKCJA
OBSŁUGI • KULLANIM KILAVUZU

MAX

Deutsch

English

Français

Italiano

Español

Nederlands

Svenska

Polski

Türkçe

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Table of Contents

Welcome to the Team	19
Let's Get Started: Your First Training Session	19
Contents of this Set	19
Legend of Symbols	20
Warnings and Safety Instructions	20
Charging the controller	24
Moistening of Electrode Pads and Skin	24
The correct attachment of the trainer	25
Handling of the controller	26
Meaning of the LED & Screen Symbols	26
Bodify® Training Programs	27
Screen Mirroring	27
Process of the training	28
How often should you exercise?	28
What intensity level should you choose?	28
Training example	29
Customer Support	29
Frequently asked questions & answers	29
Troubleshooting Instructions	30
Cleaning and Maintenance	30
Disposal	30
Standards and Norms (Declaration of Conformity)	31
Technical Specifications	32
Manufacturer Information	32

Welcome to the Team

Here you will find all the important information about your new EMS trainer. To help you get started as quickly as possible, you will find a quick overview in the next section.

If you have any questions, feel free to contact our support at any time. We will always help you quickly and find a solution that will make you happy - we promise!

✉ hallo@bodify.me

☎ +49 176 8525 4881 (weekdays, 11-16 Uhr)

Let's Get Started: Your First Training Session

1. Familiarize yourself with the user manual .
2. The controller is also the battery of the trainer and should be charged before first use. Only charge the controller with the cable provided and use a commercially available USB adapter.
3. Before you start with the EMS training, the electrode pads (pads) and your skin must be moistened with water before each session.
4. Your skin should not be sweaty or greasy before use. Attach the trainer so that the controller is above your belly button. Then attach the strap to your back. The strap should hold well but not be too tight.
5. Your trainer offers four different training programs with 15 intensity levels. You can choose between a 15-minute or 25-minute program. Initially, start with a lower level (level 1-3).
6. Your trainer recognises whether there is skin contact and always starts in the first level (after switching on).

Contents of this Set

1 EMS Controller, 1 belt, 1 extension belt, 1 USB-C cable

Legend of Symbols

We use the following symbols in this user manual, on the product, and on the packaging.

	Warning: Important warning notice for potential injuries or health risks.
	Caution: Safety notice for potential damage to the device.
	Manufacturer
	CE Marking: This product complies with European and national regulations.
	UK Conformity Mark
	Observe user manual
	BF-Type applied part
	Batch number
	Disposal notice for packaging and product components. Dispose according to regulations.
	Environmentally friendly disposal of packaging
	Disposal according to the WEEE Directive (Waste Electrical and Electronic Equipment)

Warnings and Safety Instructions

Your health is important to us. Therefore, use your trainer responsibly and only for its intended purpose. In general, only train if you are an adult, healthy, and feel well.

General Warnings

Use Bodify products only for their intended purpose:

- For human use only and for the purpose described in this manual.
- For external use only.
- Use only with original accessories.

EMS should not be used in the following areas:

- EMS should not be applied to the head, eyes, mouth, or carotid artery.
- In the genital area.
- Avoid using in humid environments, such as bathrooms, showers, or heavy rain.

⚠ Absolute Contraindications

If any of the following conditions apply to you, do not use EMS to prevent health damage:

- You have a pacemaker or another electrical device implanted in your body.
- You have metallic implants.
- You use an insulin pump.
- You have a high fever (above 38°C).
- You suffer from known heart rhythm disturbances or other heart problems.
- You have a neurological disorder causing seizures or other sensory disturbances.
- You are pregnant or have recently given birth.
- You have cancer or are undergoing tumor treatment.
- You have recently had surgery.
- You are connected to a high-frequency surgical device.
- You use other electrical devices that affect muscles or nerves.
- You have an acute or chronic gastrointestinal disease.
- You suffer from acute muscle cramps or muscle pain.
- You have an acute infection or inflammation.
- You have Type 1 diabetes.
- You have arteriosclerosis or circulation disorders (e.g., high blood pressure).
- You have injuries, inflammations, or wounds in the area of the electrodes.
- You are under 18 years old.
- You have recently consumed alcohol or drugs.

⚠ Relative Contraindications

We recommend that you consult a doctor if any of the following apply. The doctor should explicitly approve the training.

- You suffer from pain or have received a specific diagnosis.
- You have had a herniated disc.
- You have implants that have been in your body for more than 6 months.
- You have problems with internal organs.
- You have cardiovascular issues.
- You are prone to thrombosis or embolisms.
- You suffer from chronic pain, regardless of where it occurs.
- You take appropriate medications.

- You wear an intrauterine device (e.g., for contraception).
- You are breastfeeding.
- You have a sensitivity disorder (e.g., due to metabolic disorders), which leads to reduced perception of pain.
- You have irritated or inflamed skin.
- You are prone to bleeding, or there is a risk of injury.
- You have movement limitations.
- You have fluid buildup, such as edema or swelling.
- You have open skin wounds, eczema, burns, or injuries.
- For any other conditions, procedures, or sensitivities that are not mentioned here.

Possible Side Effects:

If you train with EMS despite contraindications or without following the warnings and safety instructions, it may lead to health problems. Possible side effects include overtraining, partial strength loss, or delayed muscle soreness. Therefore, pay attention to rest phases for your muscles. If you follow our training instructions, these problems should not occur. If you notice any of the following symptoms, consult a doctor.

Important Training Considerations:

Start with a low intensity at first. Only increase the intensity when you get used to the device. Train each muscle group a maximum of 2-3 times per week. Give your muscles 1-2 days of rest after training.

- Each training session should not exceed 25 minutes, depending on intensity and load.
- Listen to your body and do not exceed your limits.
- Drink enough water to stay hydrated during training.
- If you feel unwell during training, stop immediately and consult a doctor.

⚠ General Safety and Precautionary Measures

- Keep devices away from heat sources like heaters or direct sunlight. Avoid use near microwave ovens.
- Do not expose the device to moisture. Keep it dry and protect it from dirt.
- Never immerse the controller in water or use it in humid environments.
- The products should not be used by people with contraindications or by animals.

- If the device does not work correctly or causes pain, stop the training immediately.
- Do not use the device while sleeping or when operating vehicles or machinery.
- Never use your trainer during activities requiring your full attention, agility, and fine motor skills.
- During stimulation, avoid wearing metallic objects such as belt buckles or necklaces. If you wear jewelry or piercings, remove them before use.
- Avoid using EMS if you have a tattoo with metallic particles under or near the electrodes.
- Do not use EMS products near flammable substances, such as gases or explosives.
- If you experience dizziness or weakness during training, stop and consult a doctor.
- EMS products are not suitable for children or individuals with physical or cognitive impairments.
- Keep EMS devices away from children and animals, especially the packaging material, to avoid accidents.
- Only turn on your controller once it is properly positioned and placed in the correct area.
- Our EMS products should only be used by individuals who have read and understood all safety and contraindications in this manual.

General Note:

This product is neither a medical device nor an active therapeutic product under Regulation (EU) 2017/745 (MDR). This product is not intended for the treatment or alleviation of a disease, injury, or disability, nor does it serve to maintain, alter, replace, or restore such conditions.

Product Damage:

Do not use the device if it is damaged, and contact Bodify customer support.

- To ensure that the device functions properly, it should not be dropped or disassembled.
- Regularly inspect the device for damage. If you notice any damage or if the device has been improperly used, return it to the manufacturer for inspection or repair.
- Immediately turn off the device if a defect or malfunction occurs.
- Do not attempt to open or repair the device yourself.
- The manufacturer assumes no liability for damages caused by improper or incorrect use.

Battery Instructions

- If liquid from the battery comes into contact with your skin or eyes, rinse the affected area immediately with water and seek medical attention.
- Keep the controller away from excessive heat and moisture.
- Explosion risk! Do not throw the controller into a fire.
- Do not open, disassemble, or crush the controller.

Charging the controller

The controller has an integrated and rechargeable battery. You can easily charge the controller with the included USB-C cable. (The controller does not need to be opened for charging!) To do this, plug the cable into the connector on the controller. Use a standard USB adapter. The controller charges correctly when the LED blinks red. As soon as the LED lights up blue, the controller is fully charged. A full charging cycle takes about 3 hours.

- ✔ USB-C charging port
- ✔ After starting, the controller shows the current battery level.
- ✘ Do not leave the controller plugged in for an unnecessarily long time

The controller must not get wet or be cleaned wet. If the controller has water damage, do not continue to use it! Do not expose the controller to excessive sunlight! Never open or repair the controller yourself!

Moistening of Electrode Pads and Skin


Before you start with the EMS training, the **electrode pads and your skin must be moistened with water before each session**:

1. First, remove the controller to avoid accidental activation. Thoroughly moisten a cloth with water. Tip: Kitchen paper works great. Then wipe over the pads and the corresponding area on your skin. A noticeable film of water should remain on the pads and your skin.

2. Make sure **all three pads and your skin are evenly moistened**. This is the

only way the impulses can be transmitted evenly, ensuring optimal EMS training.

⚠ Use your trainer only when the pads and your skin are moistened before training. Do not hold the belt under running water or immerse it in water. Ensure your skin is clean and free of creams or oils. Keep the pads clean and avoid contact with creams, cleaning agents, sharp objects, etc.

The correct attachment of the trainer


Position the trainer so that the controller is above your belly button. The pads should cover your abdominal muscles well. After placing the trainer on your abdomen, attach the strap to your back. The belt should fit so that the trainer is firmly held in place. The EMS belt should not be fastened too tightly or constrict. If the belt of the trainer is not sufficient, use the belt extension provided.

- Check if the electrode pads and your skin are moistened.
- The controller should be at the same height as your belly button.
- Do not switch on the trainer until it is correctly placed.


Handling of the controller

The controller can be easily removed from the trainer and reattached to it. Hold the trainer and carefully pull off the controller. The illustration shows the simple operation of the controller.

Meaning of the LED & Screen Symbols

Note: The different colors (blue, yellow, white & green) represent the various training programs. Once the controller is turned on, you can switch between the training programs.

-  Off: The LED is not lit. Your trainer is therefore switched off and not active.
-  Blue, yellow, white, or green lit: The trainer is active and running the respective training program. Press the "M" button to switch between training programs. As soon as the timer reaches 00:00, the trainer will automatically turn off.
-  Flashing and Δ symbol: Your trainer is turned on but there is no skin contact. Or the controller is not properly attached. The trainer will then automatically turn off after 30 seconds.
-  Flashing: The trainer is active and running the respective training program. The battery is almost empty and should be charged soon.
-  Off and  symbol: Trainer does not start. The battery is empty and needs to be charged.

Bodify® Training Programs

With the various training programs, you can make your workout diverse. Once your trainer is activated and attached to your skin, you can switch between training programs using the “M” button. The color of the LED indicates which training program is active. **Note: If the training program is changed during the workout, the intensity level will reset to 1.**

Blue: The strongest and most advanced abdominal training. Comparable to an intense workout with many varied repetitions.

- ✔ All-in-One Training
- ✔ 25 minutes

Yellow: Intense strength training. Designed for long and sustained muscle contractions. Comparable to sit-ups for the muscles.

- ✔ Strength Training
- ✔ 25 minutes

White: Short power training. For muscle contractions with a focus on the lateral abdominal muscles. Comparable to a short power workout.

- ✔ Power Training
- ✔ 15 minutes

Green: Massage program for muscle relaxation. Tailored to gently loosen and relax the muscles. Comparable to a pressure massage for the muscles.

- ✔ Massage
- ✔ 25 minutes

Screen Mirroring

The screen of the controller can be rotated by 180°. To do this, press the “+” and “-” buttons at the same time for a total of 5 seconds. The rotation is only possible when the trainer is attached and a training program is active. Note: After switching off, the screen will rotate back to its original position.

Process of the training

Once your trainer is activated and attached to your skin, it will start the first of four training programs. The programs are dynamic and run through a specific routine to stimulate your muscles optimally. **Tip: You can move freely during the training.**

Soft - Level: 1

Strong - Level: 15

Note: After turning on, your trainer automatically starts at the first intensity level and in the first of four training programs. Pressing the plus button will switch you to the next higher intensity level. **For example, to reach level 3, turn on the trainer and press the plus button twice. The screen will display the number "3".**

How often should you exercise?


Train one muscle group 2-3 times a week. This is enough to train your muscles sustainably. After you have trained a muscle group, you should give your muscles 1-2 days to recover. As a general rule, only train as much as you feel comfortable with! Start at a low level 1-2. Do not use the EMS belt every day or several times on the same day.

What intensity level should you choose?

Initially, start on a low level (level 1-2) to get accustomed to the training. EMS training is many times more effective than conventional training. Remember: The training should always be comfortable for you. Reduce the intensity if you find the training too intense. It is normal to train at different levels in different training programs.

Training example

To make it easier to get started with EMS training, you can follow the following training intervals. Of course, you can adapt the training individually to your daily routine.

Customer Support

If you have any questions or need help, please contact us. We will help you quickly and easily:

✉ hallo@bodify.me

☎ +49 176 8525 4881 (working days, 11am -16pm)

Frequently asked questions & answers

- 1. How often should I train?** It's best to train the same body area 2-3 times a week. In between, you should give your muscles 1-2 days of regeneration.
- 2. What intensity should I use and when should I increase it?** Increase the intensity slowly, starting at level 1 or 2. Remember, EMS training is many times more effective than conventional training.
- 3. What can I do during training?** During your EMS training, you can move freely and continue your everyday activities such as working, cooking, or doing housework. Alternatively, you can supplement your training with Bodify by adding additional exercises.
- 4. When should I use the massage function?** The massage function is used for activating your muscles and can be specifically used to massage individual muscle groups. Note: Use the trainer no more than 2-3 times per week, including times when you use the massage function.
- 5. I think there's something wrong with my trainer, what should I do?** Contact our customer support and do not use your trainer until then. We are happy to assist you and find a quick and easy solution.
- 6. Can I use my trainer for a different area?** No, the trainer has been specifically designed for the abdominal muscles and should not be used otherwise.
- 7. The LED and screen no longer light up, even after (or during) charging?** This could be due to: The battery is empty or there is a defect in the controller. Possible solution: Charge the controller and supervise the charging process for 20 minutes. Check if the controller turns on again. If not, contact our support.

- 8. I have an uncomfortable feeling under the pads?** Possible Reasons: The pads and your skin are not sufficiently moistened. The contact between the pads and your skin is not optimal. Too many training sessions. Give your muscles a few days of rest. If the uncomfortable feeling persists, contact our customer support.
- 9. The contractions are weak/not noticeable even at higher levels?** Possible Reasons: (1) The pads and your skin are not sufficiently moistened. (2) The contact between the pads and your skin is not optimal.
- 10. Is abdominal hair a problem?** The pads should have good skin contact despite body hair. It helps to moisten your skin with a bit more water.

Troubleshooting Instructions

- **The controller won't turn on?** Possible reason: Low battery. Solution: Charge the controller.
- **The controller is blinking and showing the following symbol .** Possible reason: The controller is attached to the belt but does not detect skin contact. The controller will blink and turn off after approximately 30 seconds. Solution: Make sure the controller is properly attached to the belt, and the belt is placed on your body as described.
- **The contractions are very weak?** (1) Possible reason: The intensity level is too low. Solution: The trainer always starts at level 1. Gradually increase the intensity. (2) Possible reason: Insufficient skin contact. Solution: Moisten your skin before training in the corresponding area.
- **The stimulation is uncomfortable?** (1) Possible reason: The intensity is set too high. Solution: Lower the intensity on the controller. (2) Possible reason: The skin contact between the pads and your skin is not optimal. Solution: Moisten your skin and the pads again with water.
- **Other problem?** Contact our customer service team.

Cleaning and Maintenance

Only clean the controller with a soft, dry microfiber cloth. Do not use chemical cleaners. The belt, charging cable, and controller must never be washed in a washing machine. Store your EMS trainer and its accessories in a cool, dry place.

Disposal

For the sake of the environment, the device should not be disposed of in household waste at the end of its service life. It can be disposed of at the appropriate collection points in your country. Please dispose of the device in accordance with local regulations. In the European Union (EU), the WEEE Directive (Waste Electrical and Electronic Equipment) applies. If you have any questions, contact your local disposal authorities. Used and fully discharged batteries should be disposed of in special collection containers, at hazardous waste points, or at electrical retailers. Proper battery disposal is legally required.

Standards and Norms (Declaration of Conformity)

This device is certified according to EU directives and marked with the CE conformity mark "CE". We hereby declare that this BODIFY product, Art. No.: 3528, complies with the essential requirements of European directives 2014/30/EU, 2014/35/EU, and 2011/65/EU. The full EU declaration of conformity can be requested from Fassmer GmbH, Weingartsweg 17D, 40668 Meerbusch, Germany.

Technical Specifications

Classification	Device with internal power supply
Purpose	Abdominal muscle stimulator
Art. Nr.	3528
Type	Bodify EMS Abs Belt Max
Storage and transport conditions	Temperature (0-55°C) / Humidity (15-75%)
Operating Conditions	Temperature (10-30°C) / Humidity (15-75%)
Weight (controller only)	27g
Intensity (adjustable)	1-15
Dimensions (controller only)	Ø 5.10 cm x 1.55 cm
Pulse form	Biphasic rectangular pulses
Pulse frequency	1-97Hz
Pulse duration	70-200µs
Battery	Lithium Polymer, 3.7V, 240mAh

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You can also request an electronic copy of this manual through our website or by contacting us. We will gladly send you the user manual as a PDF. Subject to errors and changes.

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Bodify® is a registered trademark of Fassmer GmbH
Fassmer GmbH • Weingartsweg 17d
Bodify®, Fassmer GmbH'nin tescilli bir markasıdır.