

Claritag[®]
ADVANCED

FOR THE REMOVAL
OF SKIN TAGS

READ ALL INSTRUCTIONS
BEFORE USING

CLARITAG ADVANCED •
SKIN TAG REMOVER

ENGLISH/SPANISH



TABLE OF CONTENTS

A: THE CLARITAG ADVANCED™ SKIN TAG REMOVAL KIT CONTENTS

B: WHAT IS A SKIN TAG & WHEN TO USE

C: WHEN NOT TO USE

D: DO YOU MEET THE TREATMENT CRITERIA

E: HOW IT WORKS & USER GUIDE

F: FOLLOW-UP-TREATMENT

G: WARNINGS

H: CAUTIONS

I: RISKS AND SIDE EFFECTS

J: TROUBLESHOOTING


K: PICTORIAL GUIDE OF SKIN GROWTHS

L: DEVICE CARE FOR AFTER USE

IMPORTANT!

Carefully review safety information on the following panels BEFORE USE. For Adults Age 21 years of age and older.

A: THE CLARITAG ADVANCED SKIN TAG REMOVAL KIT CONSISTS OF:



THE CLARITAG ADVANCED® SQUEEZE AND FREEZE™ SKIN TAG REMOVAL DEVICE IS INTENDED FOR THE REMOVAL OF SKIN TAGS. HERE'S WHAT'S IN THE BOX:

A 1 CLARITAG® DEVICE (WITH PRE-INSTALLED **A1** LIQUID COOLING GAS* AND **A2** PRE-ATTACHED TREATMENT HEADS)

*COOLING GAS ARE CONTAINED IN AEROSOL CANISTERS FILLED WITH 152A. THE MIXTURE IS FLAMMABLE BUT NOT HARMFUL TO THE ENVIRONMENT. CAN TREAT UP TO 10 SKIN TAGS.

B 1 DEVICE ACTIVATION BASE

C 1 SPARE SET OF TREATMENT HEADS


D INSTRUCTION MANUAL & GUIDE**

**THIS BROCHURE DESCRIBES ALL ABOUT THE PRODUCT, HOW IT WORKS AND HOW THE CLARITAG ADVANCED™ SQUEEZE AND FREEZE™ SKIN TAG REMOVER CAN BE USED TO REMOVE SKIN TAGS.

B: WHAT IS A SKIN TAG & WHEN TO USE:

IS IT A SKIN TAG?

This is a helpful guide ONLY and not intended to serve as a replacement for diagnosing your skin condition. If you are uncertain about what skin condition you have, you should be seen by a doctor.




Skin tags will typically cause no pain, unless they become irritated through friction, get caught on jewelry, nicked by a razor, or become twisted. A skin tag must feel soft and be easily movable.

You can recognize a skin tag as an overgrowth of normal skin cells, which form a "stalk". They are not harmful, and present as smooth and soft, movable flaps of skin that stick out from the surrounding area.

Although skin tags can occur anywhere on the body, they are typically located in warm, moist areas subject to frequent rubbing, such as skin creases on the neck, under the arm, under the breast, or in the groin area. Most are skin-colored and typically smaller 1-6 mm in diameter.

EXAMPLES OF COMMON SKIN TAG APPEARANCE:



For more examples of skin growths that are not skin tags, see Section L: Pictorial Guide of Skin Growths.

If you are not sure that the growth on your skin is a skin tag, see your doctor before you use the Claritag Advanced™ Squeeze and Freeze™ skin tag remover. If you treat a skin condition as a skin tag with the Claritag Advanced™ Squeeze and Freeze™ skin tag remover and it is not a skin tag, you delay proper treatment. This delay may lead to a more serious skin condition or disease.

C: WHEN NOT TO USE

DO NOT USE if you now have or ever had any form of skin cancer or skin lesions (basal cell, squamous cell or malignant melanoma), or been diagnosed as pre-cancerous anywhere on your body.

DO NOT USE on moles, warts, neurofibroma, seborrheic keratosis or other skin conditions, including irritated, infected or bleeding skin.

DO NOT USE if you have poor circulation or if you have been diagnosed with blood conditions affected by extreme cold.

DO NOT USE if you have a weakened immune system as you are more susceptible to infection. You may have a weakened immune system if you have Cancer, Diabetes, Human Immunodeficiency Virus (HIV), Hepatitis B or Hepatitis C, Epstein-Barr, Toxoplasmosis, Cryoglobulinemia, Malaria, Long COVID, Graft versus Host disease, major burns, undernutrition, or are geriatric.

DO NOT USE if you have been treated with chemotherapy, radiation, steroids such as prednisone, biologics or receive a transplant.

DO NOT USE if you are pregnant or breastfeeding.
DO NOT USE if you are less than 21 years of age.

D: DO YOU MEET THE TREATMENT CRITERIA

1 If you have a skin tag, you should use this device (see Section B pictures). Confirm that your skin tag is NOT a mole, wart or cancerous lesion (See pictures and description in Section K). If you have a mole, wart or cancerous lesion instead of a skin tag, DO NOT USE the device.

WARNING: If you are unsure or have a potentially cancerous lesion, consult your dermatologist.

2. DO NOT USE DEVICE IF:

- The skin tag is near the eyes or mouth, on the face (including inside or outside the nose or ears), on fingers or toes, near the genitals or on irritated skin.
- You have any of the following diagnoses? Pregnant or breast feeding, Diabetic, Skin cancer or precancerous skin lesions, Autoimmune Disease, Poor circulation or if you have been diagnosed with blood conditions affected by extreme cold, Hepatitis B or Hepatitis C, HIV, Epstein-Barr, Toxoplasmosis, Cryoglobulinemia, Malaria, Graft vs Host disease, Long COVID, major burns, old age.
- You been treated with any of the following: chemotherapy, radiation, steroids such as prednisone, biologics, or received a transplant.
- You are less than 21 years of age. If you meet all of the above criteria, you are ready to use the device to remove your skin tag. If you meet all of the above criteria, you are ready to use the device to remove your skin tag.

E: HOW IT WORKS



CLARITAG IS A COMPACT, HAND-OPERATED DEVICE FOR SKIN TAG REMOVAL. HERE ARE STEPS IN A NUTSHELL:
CLARITAG® USES CRYOTHERAPY (FREEZING) TO REMOVE THE SKIN TAG, USING A METHOD SIMILAR TO THAT WHICH IS COMMONLY USED BY DERMATOLOGISTS.



Detailed instructions are found in section F: Instructions for Use and Section F: Follow-up-Treatment.


F: FOLLOW-UP-TREATMENT

DO NOT USE THE CLARITAG DEVICE ON THE SAME SKIN TAG FOR MORE THAN TWO TREATMENTS (4 applications).

Post-treatment Process:

After treatment, the skin tag will change color from white to pink, red, bright red/dark red and black as seen below in the photo illustrations. The skin tag may exhibit mild to moderate swelling.

COLOR CHANGES FROM TREATMENT:



SKIN TAG WILL EVENTUALLY FALL OFF ON ITS OWN


- Keep the treated area clean. You may add antibacterial cream to the affected area twice daily.
- Do not pick or scratch the treated area. You may use a bandage.
- If a blister forms, do not puncture it.
- If your skin tag is inflamed or shows any type of infection, see your doctor.
- When the skin tag turns black, allow the skin tag to fall off on its own. Do not attempt to remove it.

CAUTION: Post-treatment it is important to prevent infection.


- Keep the treatment area clean. You may wish to wash with soap and water, swim or shower. However, keep the skin tag as dry as possible.
- Apply antibacterial cream to the affected area twice daily.
- Keep the tag dry except when cleaning or applying antibacterial cream. DO NOT apply anything other than antiseptic cream to the treatment area.
- Do not pick or scratch the treatment area. While sleeping, you can cover it with a bandage so as not to pick or scratch it inadvertently.
- To prevent any rubbing or friction on the treated area, a bandage may be applied.
- If a blister forms, do not puncture it. If it punctures accidentally, apply antibacterial cream to the affected area twice daily.
- When the skin tag turns black, allow the skin tag to fall off on its own. DO NOT attempt to remove it.
- If the skin tag is inflamed or shows any type of infection, see your doctor immediately.
- DO NOT apply lubricants such as petroleum jelly (Vaseline) on the skin tag.

STEP-BY-STEP USER GUIDE


1 PREPARE FOR TREATMENT




Check the "use by" date to confirm it's not expired, and if you see any damage to the device, do not use.



To prepare the skin tag for treatment, wipe clean and dry.



Place the Claritag® device (A) in base (B) on a flat surface.




Firmly press down on the device to activate the tips. Pressing down on the base sends the liquid cooling gas to the Claritag® "squeeze and freeze" tips.

You will hear a hissing sound as the treatment heads begin to freeze.


2 TREATMENT OF SKIN TAG

IMPORTANT: apply to skin tag within 5 seconds after activating device!




When successfully activated, both treatment heads will have partially changed color to frosty blue. The device is now ready for use - be sure to use within 5 seconds of activation!


DO NOT TOUCH FOAM PADS with bare hands when treatment tips are BLUE.




NOTE: If both treatment heads do not partially change color, the device has not been activated successfully. Repeat Step 1.4.



Place the Claritag device tips over sides of skin tag, enclosing the skin tag between the treatment tips. Get as close to the base of the skin tag as possible.



Squeeze the tips tightly together and hold for **20 seconds**. The device must be used within **5 seconds of activation** in order to work.



Release the pressure and let go of the skin tag. The first application is successful when the skin tag changes color to a frosty white.

NOTE: When the skin tag is successfully treated it will change color to a frosty white. If the skin tag does not turn frosty white, the skin tag did not freeze and has not been treated. Repeat the application a second time.

3 REPEAT TREATMENT OF SKIN TAG

REPEAT STEPS 1 & 2 AFTER ALLOWING THE SKIN TAG TO REST FOR 20 SECONDS. THE SECOND FREEZING FULLY TREATS THE SKIN TAG SO IT DOES NOT GROW BACK.

G: WARNINGS



- FLAMMABLE AEROSOL. Pressurized container. May burst if heated. DO NOT EXPOSE to temperatures exceeding 40°C/120°F.
- HIGHLY FLAMMABLE. The activated foam pads or gas are highly flammable. Keep away from open flame or spark. DO NOT SMOKE or use near an open flame.
- Contents under pressure. DO NOT puncture or destroy gas canister, even if empty.

- Use as directed. For external use only.
- KEEP OUT OF REACH OF CHILDREN.
- Do not breathe in spray, use only in well-ventilated areas.
- Do not use near eyes or on the eyelids as contact with the eye may cause blindness. Do not swallow or inhale. Do not rub eyes or bring the gas or foam pads activated with your Squeeze and Freeze™ into contact with eyes, it may cause blindness. If gas gets into your eyes, flush with water for 15 minutes and get medical help immediately.

- DO NOT USE on the face, including inside or outside the nose or mouth, ears (inside or outside), on the genitals, as use may cause burns or permanent scarring.
- DO NOT USE on fingers or toes; it could require immediate medical care.
- DO NOT USE in any area that is infected or irritated.
- The use of Claritag with light or dark skin persons may cause lightening or darkening of the skin.
- DO NOT USE on skin tags with a length greater than 6 mm (0.2 inches) (3 nickels stacked)

H: CAUTIONS

- Treat one skin tag per area at a time. If there are multiple skin tags, only treat ONE in each area at that time; treat each remaining nearby skin tag leaving a two-week interval between treatments. This minimizes discomfort in one area.
- Treating a skin tag is a two-step process. For effective skin tag removal, treatment requires two applications, which consists of freezing the skin tag for 20 seconds, allowing the skin tag to rest for 20 seconds and then re-freezing the skin tag for an additional 20 seconds.
- Do not share Treatment Tips with another person. They can become contaminated and spread infection. You can use the alternate set of Treatment Tips for another person. Mark the storage box clearly to keep them labeled by user.
- To avoid infection at the treatment site, please follow the advice in Section F.
- For best results, treatment should be administered by another adult.

I: RISKS AND SIDE EFFECTS


- Possible discomfort.** A stinging, itching or aching feeling during and after freezing that decreases after the area thaws.
- Possible changes in skin tag coloring** from white, to red/pink, to black, prior to falling off.
- Freezing treatments **may cause a change in skin color**. Dark skin may become lighter in appearance or light skin may become darker in appearance.
- Possible swelling**. Possible mild to moderate swelling.
- Possible blistering, flaky and dry skin** appearance.

J: TROUBLESHOOTING

Problem	Possible Cause	Solution
Treatment heads did not turn blue	Did not press down firmly enough to send the cooling gas to the head.	Firmly press the device down into the base.
The skin tag did not turn white	Both heads were not properly activated.	Press down firmly so both heads will partially turn blue.
	Waited too long to freeze the skin tag after activating.	Once activated, you have 5 seconds to squeeze and freeze the skin tag.
	Did not squeeze device over skin tags long enough.	You must squeeze and freeze the skin tag for 20 seconds.
Treatment heads are blue before use	The device was incorrectly stored.	Do not re-use the device until it returns to room temperature AND the heads become white again. Please read storage directions (section L).

K: PICTORIAL GUIDE OF OTHER SKIN GROWTHS

Skin Tags must feel soft and flexible, sticking out from the skin like a stalk, and are the color of the surrounding skin. Pictorial examples of skin tags are shown in Section B "When to Use".



L: DEVICE CARE FOR AFTER USE

Cleaning:

- This Claritag device may be wiped clean with warm soap and water.
- After use, put the device back into the provided case for the next use.

Storage:

- The device should be stored away from heat and sunlight. Do not expose to temperatures above 120°F (49°C).
- When properly stored and used before the "use by" date on the device, the Claritag will last for up to 2 years and provide a total of 10 skin tag treatments.

Disposal:

- Ensure the device is emptied of all gas. The empty units should be disposed in accordance with local regulations.
- Operating and Storage Environmental Specifications:
 - Operating and Storage Ambient temperature -4°F to 120°F (-20°C to 49°C).
 - Operating and Storage Humidity 5-95% non-condensing.

Warranty:

This product comes with a limited manufacturer warranty against defects. If this product is determined to be defective, DGI Technologies Corp. will replace it at its discretion. All exchanges must be authorized by DGI Technologies Corp. within 10 days of purchase.

This warranty does not cover damage caused by nature, misuse, negligence, accident, modification of or to any part of the product, improper use or attempted repair by anyone other than DGI Technologies Corp.

Customer Support:

US/Canada 1-800-353-2350
UK 0-800-088-5444
Monday-Friday 9AM - 5PM EST

Manufactured for DGI Technologies Inc
400 Candlewood Commons
Howell, NJ, 07731 USA

PN: AZ6034 R5