

IPL Hair Removal and Skin Rejuvenation Device

Thank you for purchasing our products. For most effective use of the product, please read the user manual thoroughly before use. Please read [PRECAUTIONS] before using the product.

Safety Precautions

- For their skin protection purposes, people with following characteristics are not allowed to use the product:
 - Dehydration of skin sensation
 - Local and systematic inflammation, allergic dermatitis, eczema or other skin diseases (photosensitive skin and persons using photosensitive drug, or allergic to cosmetics, clothing, metals, etc.)
 - In Pregnancy, menstruation or nursing.
- Do the sunbath or sunbed or use sunless tanning beauty product in the past two weeks.
- The product is not allowed to be used on:
 - Tanning skin
 - Tattoos, dark brown and black spot (see skin color comparison card)
 - Hemangioma
- Which cosmetic surgery has taken place

II. Technical Parameters

IPL Hair Removal and Skin Rejuvenation Device	
Tube	Xenon quartz tube
Light source	Intense Pulsed light
Energy frequency	1.4 ~ 4.2 / cm ²
Irradiation area	100 ~ 110cm ²
Wavelength	390nm±10nm
Life span of lamp	≥ 30000Times
Assest specification	AC100V ~ 240V AC 50/60Hz 3A
Weight	200g
Size of the main unit	165x65x45
Consumption power	160 ~ 35W
Working normally	relative humidity 30% ~ 75 %

User Manual

Gracias por aver adquirido i nostri prodotti. Per un uso più efficiente del prodotto, leggere attentamente il manuale dell'utente e di utilizzare l'apparecchio [PRECAUZIONI] prima di utilizzare il prodotto.

Precauciones de seguridad

- Para fines de protección de la piel, las personas con las siguientes características no pueden utilizar el producto:
 - Alteraciones de la sensación de la piel.
 - Inflammaciones locales o sistémicas, dermatitis alérgica, eczema o otras afecciones de la piel. Piel sensible o personas que usen drogas fotosensibilizantes o alergias a cosméticos, medicamentos, metales, etc.
 - Durante la gestación, la menstruación o el lactante.
- Si ha realizado alguna cirugía estética.

II. Parámetros técnicos

Dispositivo para la depilación IPL y rejuvenamiento de la piel	
Tubo	Tubo de cuarzo xenón
Fuente de luz	Luz pulsada intensa
Frecuencia energética	1,4 ~ 4,2 / cm ²
Zona de irradiación	100cm ² ~ 110cm ²
Longitud de onda	390nm ± 10nm
Durata della lampada	≥ 30000 Cicli
Specifiche dell'adattatore	AC100V ~ 240V AC 50/60Hz 3A
Peso	200g
Dimensione dell'unità di base	165x65x45mm
Consumo di potenza	160 ~ 35W
Intervallo di lavoro	umidità relativa 30% ~ 75%

Manuale dell'utente

III. Product Introduction

- The product is packed with the device, adapter, user manual and protective glasses.
- Appearance and Structure

① Button for power and pulse levels
② Flash launch button
③ LCD Display
④ Flash window
⑤ AC120V power socket
⑥ Level Display: 1-5 Levels
⑦ Fan
⑧ Step for use

Product Features

- Skin rejuvenation**
The product applies to the removal of lip hair, armpit hair, limbs and body hair, forehead hair line and other hairs affecting appearance. Painless hair removal applies to black hair and dark hairs only, not to white, gray or golden hairs. Please refer to following skin color comparison card for details.

Products use

- Before depilating, scrape off the hair on the depilating area with a hair scraper.
- Connect the power cord to the power outlet of the laser hair removal.
- Long press the switch for 3 seconds, the display light will be on, and the laser hair removal can be used normally.
- When operating, please wear the professional protective glasses attached to this product to protect your eyes.
- Place the product vertically against the skin and press the flash button to start hair removal. No more than 2 continuous flash removals to be conducted in the same place.
- Single press on the flash button is manual hair removal mode, long press on the flash button for 3 seconds is automatic hair removal mode. Automatic hair removal mode will automatically stop flash hair removal every 7 seconds, and the screen displays the automatic mode.

Note

Treatment with Lumea is not effective if you have had a light bleed hair, as light hairs do not absorb enough light. Below you see the hair colors for which Lumea is suitable and effective.

4 Steps for use

- 1. 4.1 Short press "Power" button to select an appropriate energy level. "1-5" indicates the five levels. Level 1 is minimum, level 5 is maximum. "1-5" indicates the five levels. Level 1 is minimum and level 5 is maximum.
- 2. 4.2 Wear the protective glasses and let lighting lens area apply to the skin vertically. When the display background light flashes, short press the flash launch button to send the light to the skin and remove the hair.
- 3. 4.3 Long press and hold launch button for 3 seconds when the light flashes, please stop operation in the area.
- 4. 4.4 The device is set to Auto mode.

4.5 After hair removal is completed, cool the skin with ice water or frozen water.

4.6 Moisturize the skin, use lotion to moisturize the skin.

4.7 When high energy level is used, the time of operations should be appropriately reduced.

4.8 The skin appeared red swollen and feels a sense of burn heat, please stop operation in the area.

4.9 Do not sunbath within 5 weeks before and after the hair removal operation.

Operation Level and Circle Table

Part	1-2	3-4	5-6	7-8	9-10	11-12	13-14	15-16	17-18	19-20	21-22	23-24	25-26	27-28	29-30	31-32	33-34	35-36	37-38	39-40	41-42	43-44	45-46	47-48	49-50	51-52	53-54	55-56	57-58	59-60	61-62	63-64	65-66	67-68	69-70	71-72	73-74	75-76	77-78	79-80	81-82	83-84	85-86	87-88	89-90	91-92	93-94	95-96	97-98	99-100
Area	1-2	3-4	5-6	7-8	9-10	11-12	13-14	15-16	17-18	19-20	21-22	23-24	25-26	27-28	29-30	31-32	33-34	35-36	37-38	39-40	41-42	43-44	45-46	47-48	49-50	51-52	53-54	55-56	57-58	59-60	61-62	63-64	65-66	67-68	69-70	71-72	73-74	75-76	77-78	79-80	81-82	83-84	85-86	87-88	89-90	91-92	93-94	95-96	97-98	99-100
Level	1-2	3-4	5-6	7-8	9-10	11-12	13-14	15-16	17-18	19-20	21-22	23-24	25-26	27-28	29-30	31-32	33-34	35-36	37-38	39-40	41-42	43-44	45-46	47-48	49-50	51-52	53-54	55-56	57-58	59-60	61-62	63-64	65-66	67-68	69-70	71-72	73-74	75-76	77-78	79-80	81-82	83-84	85-86	87-88	89-90	91-92	93-94	95-96	97-98	99-100
Area	1-2	3-4	5-6	7-8	9-10	11-12	13-14	15-16	17-18	19-20	21-22	23-24	25-26	27-28	29-30	31-32	33-34	35-36	37-38	39-40	41-42	43-44	45-46	47-48	49-50	51-52	53-54	55-56	57-58	59-60	61-62	63-64	65-66	67-68	69-70	71-72	73-74	75-76	77-78	79-80	81-82	83-84	85-86	87-88	89-90	91-92	93-94	95-96	97-98	99-100
Level	1-2	3-4	5-6	7-8	9-10	11-12	13-14	15-16	17-18	19-20	21-22	23-24	25-26	27-28	29-30	31-32	33-34	35-36	37-38	39-40	41-42	43-44	45-46	47-48	49-50	51-52	53-54	55-56	57-58	59-60	61-62	63-64	65-66	67-68	69-70	71-72	73-74	75-76	77-78	79-80	81-82	83-84	85-86	87-88	89-90	91-92	93-94	95-96	97-98	99-100
Area	1-2	3-4	5-6	7-8	9-10	11-12	13-14	15-16	17-18	19-20	21-22	23-24	25-26	27-28	29-30	31-32	33-34	35-36	37-38	39-40	41-42	43-44	45-46	47-48	49-50	51-52	53-54	55-56	57-58	59-60	61-62	63-64	65-66	67-68	69-70	71-72	73-74	75-76	77-78	79-80	81-82	83-84	85-86	87-88	89-90	91-92	93-94	95-96	97-98	99-100
Level	1-2	3-4	5-6	7-8	9-10	11-12	13-14	15-16	17-18	19-20	21-22	23-24	25-26	27-28	29-30	31-32	33-34	35-36	37-38	39-40	41-42	43-44	45-46	47-48	49-50	51-52	53-54	55-56	57-58	59-60	61-62	63-64	65-66	67-68	69-70	71-72	73-74	75-76	77-78	79-80	81-82	83-84	85-86	87-88	89-90	91-92	93-94	95-96	97-98	99-100

Note

During hair removal operation, do not repeat or miss any area in a single treatment. When the operation is completed, the operation is completed, it is suggested to use again for a total of 2-3 times.

4.9 After hair removal is completed, cool the skin with ice water or frozen water.

4.10 Moisturize the skin, use lotion to moisturize the skin.

4.11 When high energy level is used, the time of operations should be appropriately reduced.

4.12 The skin appeared red swollen and feels a sense of burn heat, please stop operation in the area.

4.13 Do not sunbath within 5 weeks before and after the hair removal operation.

4.3 Clean the flash window and its glass with a cotton swab or plug-free cloth dipped in medical alcohol.

- 4.4 Plug the adapter to power supply and then plug its other end to the rear socket of the device. LCD display shows "OFF" and LCD display background light turns on.
- 4.5 Press down "Power button level button" and hold for 3 seconds till the product is turned on. The display shows current remaining amount of lighting, energy level and other information and operation mode.
- 4.6 Short press "Power" button to select an appropriate energy level. "1-5" indicates the five levels. Level 1 is minimum, level 5 is maximum. "1-5" indicates the five levels. Level 1 is minimum and level 5 is maximum.
- 4.7 Wear the protective glasses and let lighting lens area apply to the skin vertically. When the display background light flashes, short press the flash launch button to send the light to the skin and remove the hair.
- 4.8 Long press and hold launch button for 3 seconds when the light flashes, please stop operation in the area.
- 4.9 The device is set to Auto mode.

4.10 After hair removal is completed, cool the skin with ice water or frozen water.

4.11 Moisturize the skin, use lotion to moisturize the skin.

4.12 When high energy level is used, the time of operations should be appropriately reduced.

4.13 The skin appeared red swollen and feels a sense of burn heat, please stop operation in the area.

4.14 Do not sunbath within 5 weeks before and after the hair removal operation.

Operation Level and Circle Table

Part	1-2	3-4	5-6	7-8	9-10	11-12	13-14	15-16	17-18	19-20	21-22	23-24	25-26	27-28	29-30	31-32	33-34	35-36	37-38	39-40	41-42	43-44	45-46	47-48	49-50	51-52	53-54	55-56	57-58	59-60	61-62	63-64	65-66	67-68	69-70	71-72	73-74	75-76	77-78	79-80	81-82	83-84	85-86	87-88	89-90	91-92	93-94	95-96	97-98	99-100
Area	1-2	3-4	5-6	7-8	9-10	11-12	13-14	15-16	17-18	19-20	21-22	23-24	25-26	27-28	29-30	31-32	33-34	35-36	37-38	39-40	41-42	43-44	45-46	47-48	49-50	51-52	53-54	55-56	57-58	59-60	61-62	63-64	65-66	67-68	69-70	71-72	73-74	75-76	77-78	79-80	81-82	83-84	85-86	87-88	89-90	91-92	93-94	95-96	97-98	99-100
Level	1-2	3-4	5-6	7-8	9-10	11-12	13-14	15-16	17-18	19-20	21-22	23-24	25-26	27-28	29-30	31-32	33-34	35-36	37-38	39-40	41-42	43-44	45-46	47-48	49-50	51-52	53-54	55-56	57-58	59-60	61-62	63-64	65-66	67-68	69-70	71-72	73-74	75-76	77-78	79-80	81-82	83-84	85-86	87-88	89-90	91-92	93-94	95-96	97-98	99-100
Area	1-2	3-4	5-6	7-8	9-10	11-12	13-14	15-16	17-18	19-20	21-22	23-24	25-26	27-28	29-30	31-32	33-34	35-36	37-38	39-40	41-42	43-44	45-46	47-48	49-50	51-52	53-54	55-56	57-58	59-60	61-62	63-64	65-66	67-68	69-70	71-72	73-74	75-76	77-78	79-80	81-82	83-84	85-86	87-88	89-90	91-92	93-94	95-96	97-98	99-100
Level	1-2	3-4	5-6	7-8	9-10	11-12	13-14	15-16	17-18	19-20	21-22	23-24	25-26	27-28	29-30	31-32	33-34	35-36	37-38	39-40	41-42	43-44	45-46	47-48	49-50	51-52	53-54	55-56	57-58	59-60	61-62	63-64	65-66	67-68	69-70	71-72	73-74	75-76	77-78	79-80	81-82	83-84	85-86	87-88	89-90	91-92	93-94	95-96	97-98	99-100
Area	1-2	3-4	5-6	7-8	9-10	11-12	13-14	15-16	17-18	19-20	21-22	23-24	25-26	27-28	29-30	31-32	33-34	35-36	37-38	39-40	41-42	43-44	45-46	47-48	49-50	51-52	53-54	55-56	57-58	59-60	61-62	63-64	65-66	67-68	69-70	71-72	73-74	75-76	77-78	79-80	81-82	83-84	85-86	87-88	89-90	91-92	93-94	95-96	97-98	99-100
Level	1-2	3-4	5-6	7-8	9-10	11-12	13-14	15-16	17-18	19-20	21-22	23-24	25-26	27-28	29-30	31-32	33-34	35-36	37-38	39-40	41-42	43-44	45-46	47-48	49-50	51-52	53-54	55-56	57-58	59-60	61-62	63-64	65-66	67-68	69-70	71-72	73-74	75-76	77-78	79-80	81-82	83-84	85-86	87-88	89-90	91-92	93-94	95-96	97-98	99-100

Note

During hair removal operation, do not repeat or miss any area in a single treatment. When the operation is completed, the operation is completed, it is suggested to use again for a total of 2-3 times.

4.9 After hair removal is completed, cool the skin with ice water or frozen water.

4.10 Moisturize the skin, use lotion to moisturize the skin.

4.11 When high energy level is used, the time of operations should be appropriately reduced.

4.12 The skin appeared red swollen and feels a sense of burn heat, please stop operation in the area.

4.13 Do not sunbath within 5 weeks before and after the hair removal operation.

4.3 Clean the flash window and its glass with a cotton swab or plug-free cloth dipped in medical alcohol.

- 4.4 Plug the adapter to power supply and then plug its other end to the rear socket of the device. LCD display shows "OFF" and LCD display background light turns on.
- 4.5 Press down "Power button level button" and hold for 3 seconds till the product is turned on. The display shows current remaining amount of lighting, energy level and other information and operation mode.
- 4.6 Short press "Power" button to select an appropriate energy level. "1-5" indicates the five levels. Level 1 is minimum, level 5 is maximum. "1-5" indicates the five levels. Level 1 is minimum and level 5 is maximum.
- 4.7 Wear the protective glasses and let lighting lens area apply to the skin vertically. When the display background light flashes, short press the flash launch button to send the light to the skin and remove the hair.
- 4.8 Long press and hold launch button for 3 seconds when the light flashes, please stop operation in the area.
- 4.9 The device is set to Auto mode.

4.10 After hair removal is completed, cool the skin with ice water or frozen water.

4.11 Moisturize the skin, use lotion to moisturize the skin.

4.12 When high energy level is used, the time of operations should be appropriately reduced.

4.13 The skin appeared red swollen and feels a sense of burn heat, please stop operation in the area.

4.14 Do not sunbath within 5 weeks before and after the hair removal operation.

Operation Level and Circle Table

Part	1-2	3-4	5-6	7-8	9-10	11-12	13-14	15-16	17-18	19-20	21-22	23-24	25-26	27-28	29-30	31-32	33-34	35-36	37-38	39-40	41-42	43-44	45-46	47-48	49-50	51-52	53-54	55-56	57-58	59-60	61-62	63-64	65-66	67-68	69-70	71-72	73-74	75-76	77-78	79-80	81-82	83-84	85-86	87-88	89-90	91-92	93-94	95-96	97-98	99-100
Area	1-2	3-4	5-6	7-8	9-10	11-12	13-14	15-16	17-18	19-20	21-22	23-24	25-26	27-28	29-30	31-32	33-34	35-36	37-38	39-40	41-42	43-44	45-46	47-48	49-50	51-52	53-54	55-56	57-58	59-60	61-62	6																		