

troubleshooting



Make sure the device is fully charged.



Please choose a lower output of heat at the start of use.



Do not participate in any strenuous exercise after wearing the product.



Make sure the massager is properly fasten to the knee.

technical specifications

Power Supply	USB Type C DC5V
Battery	Lithium battery 3.7V 3000mAh (full charge about 2.5-3H)
Working time	1-1,5 hours in automatic mode
Charging current	5V/2.0A, the blue indicator light breathes when charging
Rated power	10W
Power-on Method	Long press the touch switch K1 button, short press to switch infrared, laser mode, or mode radio
Boot default:	Infrared photon skin rejuvenation, and laser skin rejuvenation
Short press once:	Photorejuvenation
Short press twice:	Infrared laser skin rejuvenation
Short press thrice:	Turn off infrared, photon, laser

When the laser indicator flashes, it means the laser is not turned on;


Vibration Mode The default automatic mode when starting up.

- Short press once: Low-grade massage
- Short press twice: High-end massage
- Short press thrice: Turn off the massage mode;

Heating mode the default low-grade hot compress when starting up

- Short press once: Mid-range hot compress
- Short press twice: High-grade hot compress

Temperature range in each gear heating mode:



Three vibration modes and three heating modes can be selected at will, and nine different experiences can be combined.

Intelligent Mode

Every time you turn on the machine, if you don't put the massager on your lap within 60 seconds, the machine will Make a beeping sound to remind you to use it. If not on the lap, the machine will automatically shut down.

The system automatically counts the time every time it is turned on and runs, and it will automatically shut down after 10 minutes of each run.

Indicator light display description


- Sufficient power: Normal use
- Insufficient battery: The blue light flashes three times within 5 seconds (when <3.4V), indicating that the battery is low, pay attention to charging.
- Charging: The blue indicator light shows breathing
- Fully Charged: The blue indicator light is always on
- Working temperature: 5°C ~ 40°C



nooro
WIRELESS KNEE MASSAGER

USER GUIDE

button panel



charging & maintenance:

- For the first use, please charge the product fully. Charging takes around 2-3 hours.
- The flashing blue light indicates low charging. When fully charged, the blue blinking light will stop.
- Input current: 5V2A Power supply: 110V-220V
- Turn the device off every time you stop using it by pressing the off button until the light goes off.
- Keep the instrument out of the reach of children, away from humidity, high temperature, direct sunlight and water splashes.

say hello

to your new nooro
knee massager.



4 simple steps on how to use the knee massager



1.
Firstly, fasten the nooro
massager to your knee.



2.
Turn it on or off, simply press
the center button with the
power icon. Long press the
power button to enter the
automated mode. When
the device is turned on,
the blue light illuminates.



3.
Use the button on the left w/
the {{{}} jagged lines to select
between the three vibration
modes. Having selected
what feels best for you.



4.
It's time to select the heat
level of your massage!, you
can do this by pressing
the heat button | | | |
on the right.

Note: Use it 1-2 times daily for 10-15 minutes per session.