

Paiaite Dog Recovery Suit User Manual

1. Product Introduction

This dog recovery suit is specially designed for dogs after surgery (such as neutering, spaying, or other skin surgeries). It is made of soft, skin-friendly, breathable and **elastic fabric**, which can effectively protect the dog’ s surgical incision, prevent the dog from licking or scratching the wound, and promote wound healing. The suit is easy to put on and take off, with a row of buttons on the back for adjustment, and fits most dog body types with its elastic design.

2. Product Specifications

Brand Size	Pet Weight (lb)	Chest Girth (in)	Neck Girth (in)	Back Length (in)
XS-(1~3lb) Chest(~11.1")	1 - 3	11.1	9.1	8.8
S-(3~8lb) Chest(~13.7")	3 - 8	13.8	10.2	10.2
M-(8~12lb) Chest(~15.5")	8 - 12	15.8	11.8	12.6
L-(11~16lb) Chest(~17.5")	11 - 16	17.7	13.4	14.5
XL-(17~23lb) Chest(~19.5")	17 - 23	19.7	15.8	16.3
XXL-(22~32lb) Chest(~21.5")	22 - 32	21.7	16.9	18.3

Note: Please measure your dog’ s chest circumference, neck circumference and back length carefully before purchasing to choose the most suitable size. The suit is made of elastic fabric, which has a certain stretching range. If your dog is between two sizes, we recommend choosing the larger size for more comfortable wearing.

3. How to Use

3.1 Preparation

1. Check the product for any damage (such as loose threads, broken buttons, etc.) before use. If there is any defect, please do not use it and contact the seller in time.
2. For male dogs: Please cut off the reserved circled part on the crotch area of the suit before wearing (see Section 4.1 for details). For female dogs, no cutting is required.
3. Ensure the dog's surgical incision is clean and dry before putting on the suit.

3.2 Wearing Steps

1. Unbutton all the buttons on the back of the suit first.
2. Let the dog stand or lie down comfortably. Put the dog's front paws into the two armholes of the suit respectively.
3. Pull the suit gently to cover the dog's back and abdomen, making sure the surgical incision is fully covered by the suit.
4. Fasten the buttons on the back one by one to adjust the fit. Thanks to the elastic fabric, the suit can fit snugly but not tightly. You should be able to insert 1-2 fingers between the suit and the dog's body, which means the tightness is appropriate.
5. Check again to make sure the suit does not rub the dog's wound, and the dog can move freely (such as walking, sitting, lying down) without discomfort.

3.3 Taking Off Steps

1. Unbutton all the buttons on the back of the suit.
2. Gently pull the suit off from the dog's front paws, and then take it off completely.

4. Important Notes

4.1 Special Reminder for Male Dogs

This suit is designed with a reserved circled part on the crotch for male dogs. **Please cut off this circled part with scissors before wearing.** If you wear it directly without cutting, the crotch area will be tight and uncomfortable for the male dog, which may cause discomfort and affect the dog's movement.

Cutting steps: 1. Find the pre-marked circle on the crotch of the suit. 2. Use clean scissors to cut along the edge of the circle. 3. Check the cut edge to ensure there are no sharp threads that may scratch the dog's skin.

4.2 Daily Care

1. Change the suit regularly: It is recommended to change the suit every 2-3 days, or when it is dirty, wet, or stained. Keep the suit clean to avoid bacterial growth and infection of the wound.
2. Cleaning method: Hand wash with mild detergent in cold or warm water (below 30°C). Do not bleach, do not dry clean, do not iron at high temperature. Hang it to dry naturally in a well-ventilated place. Note: The elastic fabric should be washed gently to avoid damaging the elasticity.
3. Do not let the dog chew or bite the suit, which may damage the product and pose a choking hazard.

4.3 Wearing Time

The suit should be worn until the dog's surgical incision is completely healed (usually 7-14 days, depending on the dog's recovery situation). Please consult your veterinarian for the specific wearing time.

If the dog shows obvious discomfort (such as restlessness, difficulty breathing, skin redness and swelling) when wearing the suit, please take off the suit immediately and contact your veterinarian.

4.4 Storage

After use, clean and dry the suit completely, then fold it and store it in a clean, dry and ventilated place, away from direct sunlight and high temperature. Avoid extrusion to prevent the elastic fabric from losing elasticity.

5. Troubleshooting

Problem	Possible Reason	Solution
The dog feels uncomfortable at the crotch (male dog)	Did not cut the reserved circled part on the crotch	Take off the suit, cut the circled part as required, and then put it on again
The suit is too tight	Incorrect size selection or over-fastened buttons; insufficient use of fabric elasticity	Replace with a larger size, or loosen the buttons appropriately to make full use of the elastic fabric
The suit is too loose, cannot protect the wound	Incorrect size selection or loose buttons	Replace with a smaller size, or fasten the buttons tightly to fit the dog's body

The suit is easy to fall off	Loose buttons or incorrect size	Fasten the buttons tightly to make it fit snugly, or replace with the correct size
The fabric loses elasticity	Excessive cleaning temperature or violent washing	Wash according to the correct method; replace with a new suit if necessary