RENPHO

User Manual

Pillow Massager

Model: RF-NM067

Contents

English..... 01-05

Deutsch..... 06-11

Français..... 12-17

Italiano..... 18-23

Español..... 24-28

IMPORTANT SAFETY INSTRUCTIONS

READ ALL INSTRUCTIONS CAREFULLY BEFORE USING.

Please consult your doctor if you are any of the following user types:

- Pregnant and postpartum women
- Persons with a pacemaker or implanted medical devices
- Persons with diabetes or other sensory impairments
- Persons with open wounds, discolored areas, or any part of the body that is swollen, burned, inflamed, or has skin eruptions or sores
- Persons with increased or reduced sensitivity to heat
- Persons with varicose vein, phlebitis, thrombosis, joint dysfunction, inflammation, or injuries
- Persons with osteoporosis, weak bones, abnormal or curved spine, or other acute pains
- Persons who have recently undergone cardiac or visceral surgery

If you feel unwell or experience any unusual discomfort, stop using it immediately.

CAUTION

1. The device is intended for household use only, not for medical or commercial purposes.
2. This device can be used by children aged 8 years and above and persons with reduced physical, sensory or cognitive capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the device in a safe way and understand the hazards involved.
3. Children should not play with the device. Cleaning and user maintenance should not be done by children without supervision.

DANGER - To reduce the risk of electric shock:

1. Always unplug this device from the electrical outlet immediately after use and before cleaning.
2. Do not use while bathing or in a shower. Do not use while driving.
3. Do not place or store the device where it can fall or be pulled into a tub or sink, and do not place it in or drop it into water or any other liquid.

WARNING - To reduce the risk of fire, electric shock, or injury to persons:

1. The device should never be left unattended when plugged in.
2. Do not operate under a blanket or pillow. Excessive heating can occur and cause fire, electric shock, or injury to persons.
3. Close supervision is necessary when this device is used by, on, or near children, invalids, or disabled persons.
4. Use this device only for its intended use as described in this manual. Do not use attachments not recommended by the manufacturer.
5. Never operate this device if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water. Return the device to the service center for examination and repair.
6. Do not carry this device by supply cord, or use the cord as a handle.
7. Keep the cord away from any heated surfaces.
8. Do not operate where aerosol (spray) products are being used or where oxygen is being administered.
9. Keep long hair, loose clothing, and jewelry away from the device while in use.
10. The device has a heated surface. Persons insensitive to heat must be careful when using it.
11. Use heated surfaces carefully, as they may cause serious burns. Do not use over insensitive skin areas or in the presence of poor circulation. The unattended use of heat by children or incapacitated persons may be dangerous.
12. Temperatures sufficiently high to cause burns may occur regardless of the control setting. Do not use on an infant or invalid or on a sleeping or unconscious person. Check the skin in contact with the heated area of the device frequently to reduce the risk of blistering.

SAVE THESE INSTRUCTIONS


WHAT'S IN THE BOX

Pillow Massager x 1
Pillow Cover x 1
Power Adapter x 1
User Manual x 1

SPECIFICATIONS

Product Name: Massage Pillow
Model: RF-NM067
Input: 12V $\equiv$ 2A
Rated Power: 24W
Heat Range: Up to 118°F/48°C
Auto-off Timer: 15mins

OVERVIEW

OPERATION

Power/Heat Button

Press and hold for 2 seconds to turn on the massager. (Heat On, Low Speed)

Press again to enable/disable the heating function.

Press and hold for 2 seconds to turn off the massager.

Speed/Direction Button

Press to cycle through speed levels among Low/Medium/High.

Press and hold for 2 seconds to switch the rotation direction.

Note:

1. Direction of massage nodes change every minute.
2. The device has a built-in auto-off timer that will shut it off after 15mins of continuous use.

CLEANING & MAINTENANCE

1. Turn off the device, unplug the power adapter, and allow it to cool completely before cleaning.
2. Wipe down the outer surfaces with a soft, slightly damp cloth. Do not use abrasives, solvents, chemical cleaners, or harsh detergents.
3. Store the device in a cool and dry location. Do not expose it to direct sunlight or high temperatures. Keep away from humidity, heat, and/or corrosive gas.
4. Do not attempt to disassemble or repair the device by yourself. Otherwise, it may cause damage and void the manufacturer's warranty.

WARRANTY POLICY

Your RENPHO product purchase is covered by a one year limited manufacturer warranty from the date of delivery. For warranty terms and conditions, please visit:

<https://renpho.com/pages/warranty-terms-and-conditions>

Note: Product registration is not required for the warranty. If you choose not to register your product, it will not diminish the product warranty.

CUSTOMER SERVICE

Please feel free to contact us if you have any questions or concerns. RENPHO Customer Service Team guarantees a quick response and hassle-free solutions to any issue you may have within business hours.

 **Tel:** +44(0)7434-666088 (UK ONLY)

Monday-Friday 9:00AM-5:00PM(GMT)

 **Email:** support-eu@renpho.com(UK&EU)

*For defective products or the return of items, please contact us with your order number within the specified warranty period. DO NOT dispose of any product parts as they may be required for inspection/repair.

Notice:

This marking indicates that this product should not be disposed with other household wastes throughout the EU. To prevent possible harm to the environment or human health from uncontrolled waste disposal, recycle it responsibly to promote the sustainable reuse of material resources. To return your used device, please use the return and collection systems or contact the retailer where the product was purchased. They can take this product for environmental safe recycling.

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EU Importer: Renpho GmbH
Manufacturer: Shenzhen Ruiyi Business Technology Co., Ltd.
C&E Connection E-Commerce (DE) GmbH
Richers & Vision Ltd.