GB

AIRMAX NASAL DILATOR

The Aimmax generates a scientifically proven increase in airflow of 76.1% on average by widening the nasal passage at the narrowest point of the nose. In this way, it helps reduce snoring and improve sleep quality and increases the intake of air during exercise. The product has been developed by airflow experts in collaboration with ENT specialists. The Aimmax is made of soft, odourless and medically certified plastic (S.E.B.S.), making it extremely comfortable and responsible to wear.

INSTRUCTIONS OF USE

Rinse the Aimmax with water before and after use. Do not use harsh cleaning agents or disinfectants on it. Use a mirror to correctly position the Aimmax and check the illustration for reference. Use the Aimmax for at least two consecutive weeks to get used to it. The Aimmax may fall out at night; this will decrease over time as use increases. Note: The Sport version is fitted with small 'barbs' to prevent it from falling out of your nose during exercise.

RECOMMENDED DURATION OF USE OF 2 MONTHS

Preferably store the product in a dry and dark place to preserve its durability of use. We recommend replacing the Aimmax every 2 months for optimal hygiene, comfort and effectiveness.

CAUTIONS

The Aimmax is only suitable for positioning in the nose and not suitable for children under the age of 12 years. Keep the Aimmax out of the reach of small children. The Aimmax should be used by one person only to avoid contamination. If you experience any difficulties using the Aimmax, we recommend that you consult an ENT specialist. No allergic reactions or other contraindications are known so far. Any serious incident in connection with this device should be reported to the manufacturer and the competent authorities of the Member State where the user is located. When disposing of this device, the general legal requirements for waste.

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Airmax B.V.
UKCA Experts Ltd.