READY SET,PRO

- USER MANUAL
- MANUALE UTENTE
- BENUTZERHANDBUCH
- MANUEL D'UTILISATION
- MANUAL DE USUARIO

R2 & R3

EASTMOUNT®

R2

R3

Hello, I'm Lukas

A sports lover changed by spinal injury. Rebuilding strength from a wheelchair taught me this: true recovery meets you where you are—whether you're a beginner or a seasoned athlete.

I tested 18 devices. Most asked too much, too soon. None offered exercise that was both simple and truly protective.

So I created tools for everyone in recovery. Our design delivers stable, gradual resistance. Start safe and build confidence. From first-timers to athletes, everyone gets an effective, low-impact workout.

Every product is science-backed, durable, and comes with 24/7 support and a 1-year warranty.

This isn't just equipment.
It's smarter, safer training for real life.
Join us at EASTMOUNT
Progress starts here — steady, safe, sure.

Lukas Wn

FIRST STEPS

WHAT'S IN THE BOX

- EAST MOUNT Twister Arm Exerciser
- User Manual

UNPACKING YOUR TWISTER ARM EXERCISER

1. Remove the components from the shipping box and inspect all items for shipping damage.
2. Carefully read all safety warnings and operating instructions before using the Twister Arm Exerciser .


SAFETY TIPS

- Never allow the handles to open in an uncontrolled manner. Always hold both grips firmly with one hand as you remove the lock bar with your other hand to avoid injury or property damage that may be caused by allowing the handles to swing open in an uncontrolled manner. Always unlock the Twister Arm Exerciser away from your body, other people, and surrounding objects to avoid injury or property damage.
- Before starting an exercise, always ensure the resistance is set to your desired weight level. Always return the resistance dial to the lowest setting after each workout to avoid future injury. Should the product get stuck, do not try to force it closed, as this may damage it. Return the resistance to the lowest level and then press the handles closed.
- If you are prone to dizziness or leg weakness, use the product while seated.
- To avoid injury, do not wear accessories such as rings, bracelets, or watches during use.
- Keep long hair and loose clothing away from the product to avoid getting it caught in the hinge.
- Adjust resistance to the lowest level and lock the product before storage or travel.

⚠ CAUTION:

- Stop immediately at the first sign of pain or discomfort, or if you begin to feel unwell.
- Keep the product away from children and pets

PARTS OF THE TWISTER ARM EXERCISER

EXERCISE TIPS

- Spend at least five minutes warming up, including stretching, before exercise and at least five minutes cooling down after exercise to avoid pulled muscles and other sports injuries.
- Always maintain control and move slowly to avoid injury.
- Always listen to your body and do not overwork your muscles to avoid injury. Take sufficient breaks.
- We recommend beginning with the lowest resistance, level 2 (22 lb). Once you get used to the exercise, increase the resistance level and number of repetitions by a little bit each week.
- Use the Twister Arm Exerciser every other day to allow your muscles sufficient time to recover from your workout.
- Tighten your abdominal muscles during exercise to prevent injury and provide a core workout.
- Breathe naturally during exercise. Never hold your breath.
- You may wish to use a pair of weightlifting gloves to prevent hand chafing, blistering, and callouses
- Drink plenty of water before and after exercise to stay hydrated

INSTRUCTIONS FOR USE

Please read these instructions before using the Twister Arm Exerciser .
Keep this manual for reference.

UNPACKING YOUR TWISTER ARM EXERCISER

1. Turn the resistance dial to the lowest setting, level 2(22 lb).
2. Position the Twister Arm Exerciser so that the lock bar faces away from you.
3. Grasp both grips with one hand. Hold the Twister Arm Exerciser away from your body, and hold the grips tightly.
4. Rotate the lock bar to unlock.

⚠ CAUTION: Ensure the product is positioned so it will not hit you, others,or surrounding objects when you remove the lock bar.

USING THE TWISTER ARM EXERCISER

1. Turn the resistance dial to the desired setting, level 2 (22 lb) through level 20 (440 lb). Use the chart below for reference.

	Easy		Moderate		Difficult		Advanced		Elite	
Level	2	4	6	8	10	12	14	16	18	20
Weight(kg)	20	40	60	80	100	120	140	160	180	200
Weight(lb)	22	88	132	176	220	264	308	352	396	440

- 2.Grasp each grip with your hands.
- 3.Push the grips towards each other slowly while keeping your abdominals tight.
- 4.Slowly pull the grips away from each other.
- 5.Repeat steps 3 and 4, opening and closing the Twister Arm Exerciser like a pair of scissors until you have completed the desired number of repetitions.
- 6.Attach the lock bar immediately after use and before cleaning or storing.

SAFETY INSTRUCTIONS

Keep this manual in a safe place for future reference.
Before using this product, read all safety warnings and operating instructions.Failure to do so may result in serious injury or death, or damage to the product or property.

CONTRAINDICATIONS

Users should consult with a medical professional before beginning a new exercise program. This is especially important for those over the age of 35and those who have pre-existing health conditions.

Your physician should assist you in determining the heart rate zone that is appropriate for your age and physical condition, They can help determine ifyou have any physical limitations that could create a health risk or prevent you from properly using this product.

Your doctor can help you determine the appropriate length of strength training sessions, how often you should exercise, and how intensely you should exercise based on your age, fitness level, and any existing health conditions you may have. Always exercise within the parameters discussed with your doctor to avoid injury or damage to your health.

Consult your physician before using this product if you:

- Are over the age of 35 or have pre-existing health conditions.
- Have osteoporosis, spinal injuries, fractures, herniated or slipped discs,chronic neck or back pain,sciatica, or other orthopedic conditions or disorders affecting the back,neck, or spinal column.
- Have joint pain, disorders, or conditions, including arthritis or rheumatoid arthritis, or if you have muscle or soft tissue damage or inflammation.
- Have severe anemia, a heart or vascular condition, hypertension, blood clots, or if you are taking medication that affects heart rate.
- Have a pacemaker, artificial heart,or another medical implant.
- Are pregnant, if you have recently had surgery, if you have cancer or malignant tumors,or if you have lung disease.
- Experience dizziness, have balance impairments, or if you are takingdrugs that affect balance, In

this case, you may need to use the product while seated.

DISCLAIMER

- The health benefits suggested or implied in this manual, other product literature, and company website are not certified or endorsed by any regulatory authority or medical institute.
- The use of this product is entirely at the user's discretion. Please read all instructions and safety information carefully before using this product. We assume no responsibility and expressly disclaim any and all liability for personal injury, death, or property damage or loss sustained through improper use, failure to follow the instructions and warning in this manual, neglect, or inadequate or improper maintenance of this product.
- It is the responsibility of the owner to ensure that all users of this product are adequately informed of all warnings and cautions.
- Use this product only as instructed in this manual.

WARNING

To reduce the risk of serious injury or death, or damage to the product or property:

- Excessive exercise or overexertion can be damaging to your health. Always exercise within the intensity and time limits recommended by your doctor.
- Stop exercising immediately and consult your doctor if you begin to feel faint, dizzy, nauseous, or experience pain, irregular heartbeat, or shortness of breath.
- Do not use the product in or near swimming pools or other bodies of water.
- Only use the provided accessories with the product. Do not modify, tamper with, damage, or attempt to disassemble the product. This product has no user-serviceable parts.
- Keep the product away from flammable or explosive items, chemicals, fire, and heat sources.
- Do not use in high temperatures or direct sunlight. Do not operate under a blanket or pillow.
- Keep away from liquids. Do not use outdoors, in high-moisture environments or any place where the product may become wet or damp. Do not immerse in water. If the product becomes wet, do not use.
- Do not run or walk around or use other exercise equipment while using the product.

- Always inspect the product for damage, wear and tear, or loose or missing screws before each use. Tighten any loose screws before use. Do not use if damaged, malfunctioning, or missing screws. Do not use the product again until the problem has been fixed.
- For adult use only. Keep away from children and pets. Children should be supervised to ensure that they do not play with the product.
- This product is not intended for use by people with diminished mental sensory, or physical capacity unless they are closely supervised by a responsible adult.

CAUTION

To reduce the risk of minor or moderate injury, or damage to the product:

- Should the product get stuck, do not try to force it closed, as this may damage it. Return the resistance to the lowest level and then press the handles closed.
- To avoid injury, do not wear accessories such as rings, bracelets, or watches while using this product.
- Keep long hair and loose clothing away from the product to avoid getting it caught in the hinge.
- Never use the product while under the influence of drugs or alcohol. Do not use during or immediately after eating or when you feel tired. Do not use if injured, ill, or if you have a fever.
- Do not allow more than one person to use the product at a time.
- Ensure adequate clear space around the product during use (at least seven feet on all sides). Remove any obstacles in the area, in case of falls or accidents.
- Always wear appropriate athletic clothing and shoes. Appropriate clothing should be comfortable but snug and allow you to move freely.

NOTICE

- Excessive force or weight applied to the product can crush or damage it. Keep away from sharp objects that could damage the product or accessories.
- This product is intended for in-home use only. Do not use in any commercial, rental, institutional, or therapeutic setting.
- Dispose of the product and all accessories according to state and federal regulations.

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

WAYS TO USE THE TWISTER ARM EXERCISER

You can strengthen your entire upper body with the Twister Arm Exerciser . Here are just a few examples.

OR


**STRENGTHEN YOUR
SHOULDER AND BACK MUSCLES**

**STRENGTHEN YOUR
BACK AND ARM MUSCLES**

TO REPLACE THE BATTERY (R3 ONLY)

1. Ensure the counter is powered off.
2. Use a screwdriver to remove the four small screws on the back of the repetition counter.
3. Remove the back portion of the counter.
4. Turn the ShredForce over and carefully remove the front portion of the repetition counter.
5. Slip the spent battery out (you may need to use a small wooden dowel or the end of the screwdriver) and properly dispose of it.
6. Insert a new LR44 button cell battery.
7. Replace the front and back portions of the repetition counter. Screw the four screws back into place.

PRODUCT SPECIFICATIONS

Model R2 & R3

Folded Size: 43 cm × 16 cm × 5 cm (17" x 6.30" x 1.97")

Net Weight: 2.3 kg (5 lb)

Resistance Levels(10): 10 - 200 kg (22 - 440 lb)

Operating Temperature: -10°C +50°C (14-122°F)

Materials: Carbon Steel

Frame: Carbon Steel

Handles: PU Rubber, ABS

STRENGTH TRAINING

TWISTER ARM EXERCISER

CHEST & ARM TWISTER

ONE-YEAR WARRANTY

EASTMOUNT® creates quality, durable exercise products and we stand by that quality with one-year warranty on all of our products. If your Twister Arm Exerciser ever breaks (and we doubt it will), we'll send your replacement parts and show you how to repair it. If your Twister Arm Exerciser cannot be repaired, we'll replace it-free of charge.

More Products Available

R2

R3

R2 PRO

R8 PLUS

YOU'VE GOT THIS

chinacheers@hotmail.com