User Manual

Model No.:TK-HP2412F
Product Type:Heating Pad
Serial No: 6789202405051524

Model No.: TK-HP2420F
Product Type: Heating Pad
Serial No: 6789202405053589

EN

This is an electronic product. Please be sure to read the instructions below before using the heating pad. Failure to do so may result in electric shock or fire, even serious physical injury.

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
DE-5

This is an electronic product. Please be sure to read the instructions below before using the heating pad. Failure to do so may result in electric shock or fire, even serious physical injury.
This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety. Children must not play with the appliance.
Cleaning and user maintenance must not be carried out by children, unless they are aged 8 or over and are supervised.
This appliance can be used by children aged from 8 years and above and by persons with reduced physical, sensory or mental capabilities or lack of experience and/or knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved.
The appliance is not to be used for warming animals.
This appliance has a polarized plug (one blade is wider than the other). To reduce the risk of electric shock, this plug is intended to fit in a polarized outlet only one way. If the plug does not fit fully in the outlet, reverse the plug. If it still does not fit, contact a qualified electrician. Do not modify the plug in any way.

WARNING! Electric shock Risk
Important instructions! Retain for future use.

- Read all instructions before using.
- For household and indoor use only.
- Do not use it if there is any obvious damage. (cracks, torn or exposed wires and heating elements).
- If the item is overheating or there is burnt smell during use, use, please unplug it immediately to turn off the heating pad.
- Patients with severe heart disease, those who rely on pacemakers and pregnant women should consult a doctor before using.
- This appliance must not be used by persons insensitive to heat and other very vulnerable persons who are unable to react to overheating.
- Do not make any modifications to the switch buttons, wires and heating elements.
- Do not reconnect the heated product or use it after cleaning until it is completely dry.
- Do not swing the wires. Do not use cord as handle.
- Do not sit or lie on the heating pad, put it on your body instead of underneath, the heat may increase when the pad stays close to the body.
- Do not pierce it or use any sharp objects to secure the heating pad.
- Do not use folded or rucked the heating pad during use.
- Do not heap or cover anything (incl. mattress or quilt) on the heating pad during use.
- If the supply cord is damaged, it must be replaced by the manufacturer, its service agent or similarly qualified persons in order to avoid a hazard.
- Do not insert pins.
- Do not use any other electronic heating devices while using the heating pad.
- Unplug it when not in use.
- Do not use if wet, ensure that the power cord never gets wet or damp when the pad is in use. Route the power cord in such a way that it cannot be pinched or damaged. Risk of electric shock!
- Do not dry clean the heating pad. Dry cleaning solvents may damage the wires.
- Do not use bleach or cleaning solution
- Do not use a wringer and do not dry in the dryer.
- Do not iron the heating pad.

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
- Do not use any mothballs or sprays.
- Very sensitive wash at 30°C.
- Do not use while sleeping or better set the timer before sleeping. When the appliance is likely to be used for an extended period, for instance by the user falling asleep, controls should be adjusted to a setting recommended for continuous use.
- Children under the age of three are not to use this appliance due to their inability to react to overheating.
- The appliance is not to be used by young children over the age of three unless the controls have been pre-set by a parent or guardian, and unless the child has been adequately instructed on how to operate the controls safely.
- Do not crease the appliance by placing items on top of it during storage.
- Examine the appliance frequently for signs of wear or damage. If there are such signs, if the appliance has been misused or does not work, return it to the supplier before switching it on again.
- Always store the appliance indoors. Do not place or store the appliance where it can fall or be pulled into a tub or sink. To prevent accidents, keep the appliance in a dry location when not in use.
- When storing the appliance, allow it to cool down before folding.
- This product should not be disposed with other household wastes. To prevent possible harm to the environment or human health from uncontrolled waste disposal, recycle it responsibly to promote the sustainable reuse of material resource.
- This appliance is not intended for medical use in hospitals.

PRODUCT OVERVIEW

UNIT FEATURES

【Extra-Large Pad】Offers maximum coverage for the lower back, shoulders, legs and other large muscle groups.
【Timer】4 settings range from 30 to 120 minutes, help to save energy and prevent burns caused by overheating.
【7-Foot Power Cord】Makes it easier to use on your comfort bed or favorite chair.
【Machine Washable】For easy care, keeping the pad as beautiful as new.
【6 Temperature Settings】

Number	Level	Temperature
1	Low	~104°F/40°C
2	Medium Low	~112°F/45°C
3	Default	~119°F/48°C
4	Warm	~126°F/52°C
5	Medium High	~133°F/56°C
6	High	~140°F/60°C

Caution: The heating pad is not intended for hospital or commercial use. It can't be used for infants, toddlers, heat sensitive person.

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OPERATING INSTRUCTIONS

- Plug into a safety tested power outlet.
- Power on: Push the switch key to ON, the heating pad goes into the default mode, which is 119°F/48°C (1/2/3 heating indicator and 60 minutes timing indicator is on).
- Temperature Adjustment: 1-6 grades corresponds to temperature range from 104°F-140°F/40°C-60°C. You can press the heating button to adjust the temperature to meet your requirements, the heating indicator will light to the corresponding temperature.
- Timing: The default time is 60min. Press the timer button to set the heating time. The timer indicator will be on to the corresponding timer--4 levels (30/60/90/120min). When the timer is off, all LED lights will be off and the heating pad will turn off automatically, in this case, press the heating button or the timing button to return to the default mode.
- For continuous use, readjust the timer and temperature settings as described above again. (The highest level is for 2 hours, after 2 hours it will turn off automatically).
- Power Off: Push the switch button to OFF.

PARAMETER

Material	Super Plush Flannel Fabric
Temperature	6 Temperature Levels (104°F-140°F/40°C-60°C)
Timer	0.5/1/1.5/2hour Auto-off
Safety Protection	Over-heating Protection

CLEANING & STORAGE TIPS

- It is recommended to clean the heating pad only in necessary conditions. If it were washed frequently, the product may wear out quickly.
- Cleaning: We suggest use a cloth or a damp sponge to remove small stains and add a little liquid detergent to clean the fabric.
- Washing: Turn off the heating pad, unplug the pad from the wall socket then remove the controller and power supply. We recommend to put the heating pad in a laundry bag, don't spin-dry it in washing machine or expose it in the strong sunlight. Do not use chemical cleaners when cleaning. Do not rotate, wring out or press hard the heating pad after cleaning. The heating pad must be completely dried before use to prevent the risk of electric shock.
- Storage: Turn off the heating pad, and unplug it. Please store it in a dry, safe, clean and cool place.

SPECIFICATIONS

Model name:	Heating pad	
Heating pad Model No.:	TK-HP2412	TK-HP2420
Controller Model No.:	TK-HP2412F	TK-HP2420F
Size:	60x30cm	60x50cm
Power:	75 Watts	120 Watts
Input:	EU/UK version:AC 220~240V 50/60Hz CA/US version:AC 100~120V 50/60Hz	
The appliance is only to be used with the types that are marked on the appliance.		

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FR-13

This marking indicates that this product should not be disposed with other household wastes throughout the EU. To prevent possible harm to the environment or human health from uncontrolled waste disposal, recycle it responsibly to promote the sustainable reuse of material resources. To return your used device, please use the return and collection systems or contact the retailer where the product was purchased. They can take this product for environmental safe recycling.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

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