
**Bob
AND
Brad®**

User Manual

iNeck Pro
Neck And Shoulder Massager

Model:206DM

Meet Bob and Brad

Bob and Brad are dedicated to creating a healthier world. We bring to you some amazing products that we have personally vetted and truly believe in! We want only the best for you and the best is right here! Enjoy!

NEED HELP?

SUPPORT@BOBANDBRAD.COM

www.bobandbrad.com

Tel: +1(612) 567-3035

M-F 9:00 AM- 5:30 PM

Your guide to being healthy, fit, and pain-free

B&B Innovation LLC
3310 S Fairview St, Santa Ana, CA 92704

Warnings and Safety Cautions

To reduce the risk of malfunction, electric shock, injury, fire, death, and equipment damage, please carefully read all instructions and warning notices in this manual before using the product.

Precautions

- Do not pull the charging cable to avoid wire breakage or poor contact.
- When unplugging the power adapter, hold the adapter and pull it out gently.
- If the product stops working or exhibits any abnormalities, stop using it immediately.
- Do not use the product on an empty stomach, after overeating, while intoxicated, or after intense exercise.
- Infants or children are strictly prohibited from using or playing with the product.
- Do not insert fingers or toes into gaps during use to prevent injury.
- Do not operate the control panel with wet hands or feet.
- If the charging cable is damaged, stop using the product.
- If you feel discomfort during use, discontinue use immediately.
- Do not immerse the entire product in water for cleaning. If the product accidentally falls into water, it must be sent for repair before further use.
- Do not attempt to disassemble the product. Modifying the product in any way may cause personal injury or property damage.
- This product is prohibited from being used in the following locations: bathrooms, kitchens, or other high-temperature, high-pressure, humid, or near-water and near-fire areas.
- Keep the product away from sharp objects. If the exterior is damaged, discontinue use and have it replaced or repaired at an authorized service center before further use.
- Do not clean the product while it is charging.
- Read the user manual carefully. Do not use the product continuously for extended periods to prevent muscle soreness or skin damage caused by excessive massage.
- Do not use the product if you have undergone surgery on the massage area or have skin injuries, redness, inflammation, bruising, or other discomfort.
- This product is a household electronic healthcare device and is not a medical instrument. It is not intended for diagnosis, treatment, or other commercial purposes.

FCC Statement

This device complies with Part-15 of the FCC Rules.

Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference.
- (2) This device must accept any interference received, including interference that may cause undesired operation.

FCC Statement NOTE

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in residential installation. This equipment generates, uses, and can radiate radio frequency energy, and if not installed and used in accordance with the instructions, may cause harmful interference to radio communications.

However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause any harmful interference to radio or television receptions, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.
- Do not hang the unit by the power cord.

Disclaimer

Please contact your physician before using the Bob and Brad Massager if you have any physical health concerns or have any of the following conditions:

- Individuals currently under medical treatment or experiencing physical discomfort.
- Patients with diabetes, hypertension, heart disease, cerebral arteriosclerosis, or cerebral thrombosis.
- Individuals with malignant tumors.
- Pregnant women and women during menstruation.
- Individuals with osteoporosis.
- Individuals using pacemakers or implanted medical devices.
- Individuals suffering from acute diseases.
- Minors and individuals with physical or mental disabilities.


THIS PRODUCT IS NOT A MEDICAL DEVICE.

IC Warning

This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s). Operation is subject to the following two conditions:

- (1) This device may not cause interference.
- (2) This device must accept any interference, including interference that may cause undesired operation of the device.


Main Structure



Usage Instructions

Note: The default working time of the product is 10 minutes. To use it again, the device must be restarted.


1. Place the product around your neck.

2. Tighten the pull strap or secure it with the shoulder strap.

3. Turn on the function to begin massage.

Button Operations

Power Button (Turn The Massager On/Off):

Turn on: Long-press the power button for 2 seconds. The buzzer will sound for 1 second, the heating function will activate automatically, and the indicator light will turn white.

Turn off: Long-press the power button for 2 seconds. The buzzer will sound for 1 second, the product will shut down, and the indicator light will turn off.

Power Button (Switch The Massage Mode):

Short-press the power button to switch between massage modes. The buzzer will beep for 0.3 seconds, and the indicator light will briefly turn off and then back on, indicating a successful mode switch (default mode at startup: Mode 1).

Mode 1: Constant Speed Mode

Mode 2: Rhythmic Mode

Heating Button (Turn The Heating Function On/Off):

Short-press the heating button to turn the heating function on/off. The buzzer will beep for 0.3 seconds each time (Heating is enabled by default upon startup).

Heating Lamp Cover Temperature: 109.4°F (43°C)

Heating Massage Head Temperature: 100.4°F (38°C)

Intensity Button (Switch The Massage Intensity):

Short-press the intensity button to adjust massage strength. The buzzer will beep for 0.3 seconds (Low Intensity: 1 beep; High Intensity: 2 beeps), and the indicator light will briefly turn off and then back on, indicating a successful intensity switch (default at startup: Low Intensity).

Low Intensity


High Intensity

Indicator Light Instructions

Device Power	Low Battery Warning	Red Flashing Light
	Auto Shutdown	The Light Flashes Red Three Times + The Buzzer Beeps Three Times
Charging	Charging in Progress	Pulses Red Light
	Fully Charged	Solid White Light
Note: During Charging, Massage And Heating Functions Are Automatically Disabled, And Button Operations Are Inactive.		

Charging Instructions

1. Use the included charging cable to connect the product to a USB port for charging.
2. The product requires approximately 4 hours to fully charge.
3. When the battery level is low, the power indicator will flash red as a low battery warning.
4. The default working time of the product is 10 minutes. To use it again, the device must be restarted.
5. The product cannot be turned on during the charging process.
6. During charging, the indicator light on the power button will display a breathing red light until it turns solid white when fully charged. Adapter is not included.
7. Supports multiple charging methods, including computer USB ports, power banks, and 5V/2A power adapters.

Use with 5V  2A charger output

5V  2A

Do not use the product when it is being charged

Technical Details

- Name: iNeck Pro Neck And Shoulder Massager
- Model No.: 206DM
- Input: 5V  2A
- Rated Power: 10W
- Dimensions: 9.1" x 6.3" x 5.5" / 230 x 160 x 140 mm
- Weight: 1.1kg
- Battery Rated Voltage: 7.4V 

- Battery Capacity: 2500 mAh
- Charging Time: 4 Hours

Troubleshooting

PROBLEM	CAUSES	SOLUTION
Fails to power on	Low battery	Charge the product for over 4 hours
	Product malfunction	Contact support
Weak massage intensity	Intensity setting at low level	Press "Intensity" button to adjust levels
	Excessive body pressure on massage head	Reduce body pressure on massage head
Sudden shutdown during use	Low battery/power interruption	Recharge the product
	Automatic timer shutdown	Restart the device
	Overheating protection from prolonged continuous use	Allow to cool before reuse
No heating function	Heating function not activated	Press "Heating" button to activate function
Friction sound from massage head/internal components or motor operation noise	These are normal operational sounds from mechanical components, not product defects	Normal operational condition
If you encounter any other problems, contact: support@bobandbrad.com		

Guarantee

This product comes with a two-year warranty from the date of purchase. Warranty will be voided if:

- Damaged due to improper usage or storage.
- Damaged by attempting to repair or disassemble the device.
- Period of warranty exceeded.

WEEE Instructions

Correct Disposal of this product	
	This marking indicates that this product should not be disposed with other household wastes throughout the EU. To prevent possible harm to the environment or human health from uncontrolled waste disposal. Recycle it responsibly to promote the sustainable reuse of material resources. To return your used device, please use the return and collection systems or contact the retailer where the product was purchased. They can take this product for environmental safe recycling.

Battery Instructions

- Do not submerge the product into water or get it wet.
- Do not charge the product near an open flame or in any extremely hot conditions.
- Do not use or store the product near heat sources like furnaces or heaters.
- If the battery leaks or emits an odor, remove it immediately from close proximity to an open flame.
- Do not put the product into a fire or heat it up in any way.
- Do not disassemble the battery in any way.
- Do not use the product if, while the battery is in use or charging, it emits an odor; gives off excessive heat; shows signs of deformation, discoloration, or any other abnormalities. If the battery is in use or being charged, disconnect the battery from the power source.

RoHs Declaration

Our Product is compliant with RoHS 2.0 (2011/65/EU Directive) regarding restriction of hazardous substances, the raw materials and components of our product meet the following toxic and hazardous substance limits.

Cadmium (Cd): < 100 ppm	Lead (Pb): < 1000 ppm
Mercury (Hg): < 1000 ppm	Hexavalent Chromium: (Cr VI) < 1000 ppm
Polybrominated Biphenyls (PBB): < 1000 ppm	Polybrominated Diphenyl Ethers (PBDE): < 1000 ppm
Bis(2-Ethylhexyl) phthalate (DEHP): < 1000 ppm	Benzyl butyl phthalate (BBP): < 1000 ppm
Dibutyl phthalate (DBP): < 1000 ppm	Diisobutyl phthalate (DIBP): < 1000 ppm


Storage:

- When not in use, please place the product in its packaging box and store it in a safe, dry, ventilated, and cool place.
- Do not hang the product using the charging cable.
- Avoid scratching the product surface with sharp objects.
- Do not place the product in direct sunlight or high-temperature environments.
- Keep the product clean and away from high-humidity areas.

Cleaning:

- Before cleaning, ensure the product is powered off and disconnected from the charging cable.
- Note: Do not use corrosive cleaners such as benzene or thinner to clean the product.

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IT

EC Declaration of Conformity

We(UK Importer is responsible for this declaration

KOVA ASSOCIATES LTD

(Company name)

72a Kingston Road, Portsmouth, Hampshire, PO2 7PA, United Kingdom

(Company address)

KOVA.Compliance@outlook.com

(Email Address)

declare under our sole responsibility that the product

Tradename: BOB AND BRAD

Equipment: Neck And Shoulder Massager

Model No.: 206DM

ASIN: B0DTYGK6YX

to which this declaration relates is in conformity with the essential requirement set out in the Council Directive on the Approximation of the laws of the Member States relating to LVD Directive(2014/35/EU) & EMC Directive(2014/30/EU) & RoHS(2002/95/EC) product is responsible to affix CE marking, the following standards were applied:

EN 60335-1:2012+A11:2014+A13:2017

EN 60335-2-32:2003+A1:2008+A2:2015

EN 62233:2008

EN 55014-1:2017

EN 55014-2:2015

EN 61000-3-2:2014

EN 61000-3-3:2013

Full Name: Taozirui

Position: Manager

Signature: Taozirui

Date: 2021.5.6

EC Declaration of Conformity

We (EU Importer is responsible for this declaration)

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xdh.tech@outlook.com

(Email Address)

declare under our sole responsibility that the product

Tradename: BOB AND BRAD

Equipment: Neck And Shoulder Massager

Model No.: 206DM

ASIN: B0DTYGK6YX

to which this declaration relates is in conformity with the essential requirement set out in the Council Directive on the Approximation of the laws of the Member States relating to LVD Directive(2014/35/EU) & EMC Directive(2014/30/EU) & RoHS(2002/95/EC) product is responsible to affix CE marking, the following standards were applied:

EN60335-1:2012+A11:2014+A13:2017

EN60335-2-32:2003+A1:2008+A2:2015

EN 62233:2008

EN 55014-1:2017

EN 55014-2:2015

EN 61000-3-2:2014

EN61000-33:2013

Full Name: Taozirui Position: Manager

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Date: 2021.5.6

We Stand Behind Our Products!

2-Year Limited Warranty & Troubleshooting **Support@bobandbrad.com**

At Bob and Brad, your happiness comes first. Our products come with a **2-year** warranty, and our Customer Happiness Team is always here to help--anytime you need us. Happiness, Guaranteed.

Troubleshoot

Help you solve any issues with your device.

Product Knowledge

Answer questions about your device and best practices.

Initiate Returns

Provide instructions for replacements, returns and refunds.

Place New Orders

Help with placing orders online or via authorized resellers.