THERABAND

THERABAND™

ELEVATE PERFORMANCE. ACCELERATE RECOVERY.

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THERABAND
Trusted
Progression™
System

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Band mittig in der Tür befestigen und die Bandschleife um den Bauch legen.
Fasten the band in the center of the door and place the band loop around the belly.

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Pull the ends of the band forcefully downwards and forwards and very slowly give way again without completely losing the band tension. Variation: Pull the ends of the straps next to the body. 10 times per leg side.

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Rotate your upper body and pull the ends of the bands alternately to the left and right. The feet/legs remain in the middle position as a fixed point. 15-20 times.

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Raise both arms simultaneously over the side against the resistance of the band. Slowly lower to the starting position without completely losing the band tension. 15-20 times.

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Forcefully kick/ extend the band leg backwards against the resistance and then slowly return to the starting position. 10 times per side.

PRODUKTINFORMATIONEN

PRODUCT INFORMATION

This document is also available in other languages

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What began with a few colorful latex bands in 1978 in Akron, Ohio, USA, has developed into the epitome of resistance products and exercise programs to accompany various physiotherapies, rehabilitation measures and also generally to increase strength, mobility and flexibility.

Over the years, the TheraBand® range has grown in line with the wishes and needs of the growing health and fitness market. Today, therapists, professionals and recreational users can find a complete range of products and exercise programs under the TheraBand® brand to get and stay healthy and fit.

A key factor in the success of the TheraBand® brand is the ongoing, close dialog with experts. The TheraBand® Academy was founded in 1999 as a platform for research, education, training and exchange of experience. The aim is to transfer the latest research findings into products and scientifically based exercise programs.

The elastic exercise band in particular is used in a wide range of exercise programs and has been proven to lead to positive training and therapy results.

Study results show the following effects, among others:

- Increase in strength
- Pain reduction
- Balance improvement
- Improving endurance
- Fall prevention
- Lowering blood pressure
- Improving posture
- Reduction of movement restrictions
- Improved functionality
- Improved mobility and flexibility

You can always find exact references and current studies at:
www.thera-bandacademy.com/research

The exercise band is a natural product made of pure latex*. The material is coated with talcum powder after manufacture and should be regularly treated with talcum powder or powder during use. People with latex allergies can use latex-free bands from TheraBand®. The durability and safe training with the exercise bands essentially depend on proper handling:

- Check the band for small holes, lateral tears and porosity before every training session. Never train with a damaged or porous band. Avoid mechanical impacts such as rings, long fingernails or rough sneaker soles.
- Use the TheraBand® AssistTM to create a sling and do not use knots or plastic clips.
- Clamps in particular can become dangerous objects if they come loose during training.
- Store the exercise band rolled up in a bag or box and protect it from direct sunlight or heat, as these dry out the natural material.
- Regular powdering prevents the band from sticking due to moisture or sweat. Dirty exercise bands can be easily cleaned with soap and water and are ready for use again after drying and powdering.
- Disinfect smooth surfaces and medical devices after each use with ready-to-use

disposable wipes. Improper or careless use of the band can lead to serious injuries. Therefore, only perform safe and pain-free exercises. Avoid exercises that could cause the band to kick back towards the head.

- The band is not a plaything for children. Children should only practice under supervision.
- Avoid using straps that are too short. The ideal band length is 2.50 - 3 m. This prevents high increases in force, uncontrolled movements and, above all, tearing of the band (when stretched by more than 300%).
- The tape can be used in water. After use, however, the tape should be rinsed with chlorine-free water and dried thoroughly.
- The TheraBand® exercise band is a natural product made of latex whose elasticity decreases over time. With regular and frequent use, the TheraBand® should generally be replaced after approx. one year of use. If there is visible damage such as tears, holes, etc., you should replace your band immediately.

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The choice of band color and therefore the correct band thickness depends on the individual training requirements. It should be possible to repeat an exercise about 15 times with the elastic band.

Practical experience has shown that beige and yellow exercise bands are mainly used for rehabilitation after injuries and for senior citizens.

For preventive training, red and green exercise bands are ideal for untrained people.

Trained people generally get on well with green and blue exercise bands and should only use the black exercise band when they are well trained.

The colors silver and gold are mainly used in competitive sports.

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Winding and fixing the exercise band

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The exercise band is always wrapped twice around the hand or foot to prevent it from slipping.

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Lay the band flat in your hand with the end towards the small finger. Now wrap the band around your hand and secure it with your thumb.

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
Place the band flat, more than shoulderwidth apart and with the ends towards the little fingers, in both hands. Turn your palms inwards and wrap the band twice around your hands.

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
Stand with one foot on the center of the band. Loop the band around your foot and secure it with the other foot.

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Stand with one foot on the center of the band. Wrap one end of the band around the foot.

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TheraBand® training with the door anchor

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

Please observe the following instructions when using the door anchor:

- Always use your door anchor against the opening direction of the door.
- Make sure that the door is locked before you start training.
- Do not use the door anchor on glass doors.
- Only use the original TheraBand® in conjunction with the door anchor.
- If these instructions are followed, nothing stands in the way of successful TheraBand® training with the door anchor.

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

Instructions folding mechanism door anchor pocket

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Turn the bag into a door anchor with one hand movement. Fold the bag lengthwise so that the black drawstring folds outwards. A loop is now formed through which you can pass the Thera-Band®.

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Now hold the band behind the door leaf and slowly close the door. The loop is now clamped in the door. The pocket on the other side serves as an anchor and prevents it from slipping. Make absolutely sure that you always attach the door anchor against the opening direction of the door.

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Training with the exercise band

Each exercise should be repeated 15 - 20 times. The complete exercise sequence can be repeated 1 - 3 times.

The load should never be perceived as very heavy during the entire training session. If this is the case, you should switch to a lighter band. Once you have created a certain basis and find the training easy, you can increase the resistance, number of exercise sequences and training frequency.

Always keep the following in mind when training:

- Correct posture and slow, controlled execution of the Movement in the starting position, during the execution and at the end of the movement.
- Each exercise begins with a slightly pre-tensioned exercise band. This band tension is primarily used to secure the joints.
- The direction in which the exercise band is pulled has a significant influence on the muscle groups and joint system being exercised. The angle between the band and the exercising body segment (arm, leg) must therefore be observed. Resistance is at its maximum at 90°. There is no effective training resistance below 30°.
- The training should be carried out without pain. If pain occurs, take a break. If pain occurs repeatedly, seek medical advice.
- Physiotherapy and any rehabilitation treatment must be carried out under expert supervision.
- Always train the right and left side of the body, even if only one side is shown in the illustration.

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Ludwig Artzt GmbH