Troubleshooting

Issue	Possible Cause	Solution
Water leaking from nozzle	Nozzle not installed correctly	Ensure the nozzle is properly installed until it clicks into place
Device won't power on	Battery is completely drained	Recharge the device
Unable to charge/Charging issues	No output from adapter or loose charging cable	Replace the adapter or charging cable
	Dirty or oxidized charging contacts	Clean the charging base regularly and keep it dry
No water flow	Valve seal dry due to first use or long-term disuse	Add a small amount of water to the nozzle installation port before use
	Nozzle blocked by debris	Clean the nozzle or replace it
	Cracked silicone water tube	Replace the silicone tube or contact customer support for repair
	Blocked filter mesh at the silicone tube connector	Clean the filter mesh
Device turns on automatically	Water ingress into the mainboard	Contact customer support for repair
Device turns off automatically	Auto shut-off feature activated	Turn the device back on
	Battery is low	Recharge the device

If the problems cannot be solved according to the above methods, please consult the selling store or go to the designated after-sales repair point for repair.

FlossJet

Portable Oral Irrigator

User Manual

FC5240


Table of Contents

Accessories and Specifications	02
EN–Product usage steps–English instructions.....	03–06
DE–Schritte zur Produktverwendung–deutsche Anleitung	07–10
FR–Étapes d'utilisation du produit–instructions en français.....	11–14
IT–Passaggi di utilizzo del prodotto–Istruzioni in italiano	15–18
ES–Pasos de uso del producto–instrucciones en español	19–22
JP–商品使用手順–日本説明	23–26

FCC STATEMENT

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Warning: Changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures

—Reorient or relocate the receiving antenna.

—Increase the separation between the equipment and receiver.

—Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.

—Consult the dealer or an experienced radio/TV technician for help.

DANGERS—To reduce the risk of electrocution:


1. Always unplug this product immediately after using (Exception: Battery chargers while recharging batteries and permanently connected products such as toilet seat assemblies.)
2. Do not use while bathing.
- Exception Hydromassage units and similar products.
3. Do not place or store product where it can fall or be pulled into a tub or sink.
4. Do not place in or drop into water or other liquid.
5. Do not reach for a product that has fallen into water. Unplug immediately.

WARNING—To reduce the risk of burns, electrocution, fire, or injury to persons.

1. A product should never be left unattended when plugged in. (Exception: Contact lens disinfectors, toilet seat assemblies, and similar products.)
2. Close supervision is necessary when this product is used by, on, or near children or invalids.
3. Use this product only for its intended use as described in this manual. Do not use attachments not recommended by the manufacturer.
4. Never operate this product if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water. Return the product to a service center for examination and repair.
5. Keep the cord away from heated surfaces.

Accessories and Specifications

Product Accessories:

Accessories can be purchased from the vendor

Product Specifications:

Product Name		FlossJet Oral Irrigator
Oral Irrigator	Product No.	FC5240
	Rated Voltage	3.7V ===
	Power Rating	8W
	Waterproof Rating	IPX7
	Charge Time	≤4h
	Input	5V === 1A
	Adaptor Rating Value	Input: 100-240V~50/60Hz 0.3A Output: 5V === 1A 5.0W
Place of origin		China


I. Steps

Before Starting: Very Important! Ensure the device is fully charged before using it for the first time by connecting the USB cable to the charging port and turning on the power. This device cannot be used while charging.

Step 1: Fill the Water Tank


Rotate the water tank counterclockwise and pull it down to remove the water tank. Then add water into the water tank. After that, align the water tank with the buckle point and rotate the water tank clockwise to install it.

Note: Ensure that the device is turned off before filling the water tank.

Step 2: Install the Nozzle

Insert the nozzle into the hole on the top of the controller. You should hear a "click" sound on successful installation. Rotate the nozzle to adjust the direction of water outflow. If you need to replace the irrigator, press and hold the uninstall button, then remove the irrigator.


Installation

Replacement

Step 3: Select Mode

This product has 4 modes. You can select the required mode by gently pressing the mode button when the equipment is off. The functions of each mode are as follows:

1. Press the on/off button, and the product will enter the working mode. Press the mode button to switch circularly among four modes, namely normal, medium, soft and pulse modes. When the product is working in any mode, press the on/off button and the equipment will shut down.

2. During charging, the function indicator will display the charging status, the battery will display the current battery level, and the charging indicators will flash in turn. When the battery has been fully charged, the charging indicator will be on. Press and hold the mode button to display the current battery level, which is indicated by the number of lights that are on. If the charging indicator flashes during use, it reminds the user to charge the battery in time. When the battery is about to be used up, the charging indicator will flash for 4 times and then the equipment will shut down automatically.

3. Memory mode: When the device is on, it will memorize the last mode you use for more than 2 minute before you turn off the device.

4. Hold the device upright during use to make sure that it is in the optimal working condition. If the device operates with any mode for more than 120 seconds, it will be off automatically.

Step 4: Before Use

Before using the device, lean slightly over the sink, extend the nozzle into your mouth and aim it at the teeth, gums, or between the teeth. Rotate the nozzle to adjust the direction of the water outflow so that it is perpendicular to the part of your mouth you wish to clean. Ensure that your mouth is open enough so that any water flows right out.

Note:


When using the oral irrigator for the first time, start from the lowest gear. Water pressure can be gradually increased as the oral cavity adapts. You may experience itching or minor gum bleeding during use. This is normal and should typically stop after 1–2 weeks.

Step 5: Use

Press and hold the on/off button for 0.6 second to start the equipment. If you need to pause during use, press the on/off button. If you want to continue to use it, press and hold the on/off button again. This product has 120-second timed shutdown function. Before the equipment shuts down automatically, the water in the water tank will be used up. If you need to continue to use it, press the on/off button to shut down, add water and then start the equipment again. Turn off the equipment after using.

Techniques:

For best results, start with the molars (posterior teeth) and start cleaning towards the front teeth. While cleaning, move the nozzle 90° along the gum line. Stay between each tooth for a short period. Close your lips slightly to avoid the water from spraying out, but leave enough space so that water can flow freely from your mouth to the sink.

Step 6: After Use

After use, open the water tank and pour out any remaining liquid. Turn on the device and drain any liquid left until the irrigator no longer sprays out any liquid, then turn off the device.

Step 7: Charging

The indicator is blue. During use, the indicator is normally on. When the indicator flashes for 5 times, the product will be shut down, indicating low voltage. In the low-voltage state, press the on/off button, the product will not be started, the indicator will flash for 4 times, and the equipment needs to be charged in time. During charging, the indicator will flash. When the battery has been fully charged, the indicator will be normally on.

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