Vibrating and Rolling Foot Massager

Model: R-D200R

CONTENTS

- ENGLISH 01-11
-

- DEUTSCH 12-22
-

- FRANÇAIS 23-33
-

- ITALIANO 34-44
-

- ESPAÑOL 45-55
-

Thank you for your purchase of the RENPHO Vibrating and Rolling Foot Massager. This device provides a powerful massage. To ensure safety and avoid any property damage, please read this manual carefully and use as indicated. Non-compliance may cause serious injuries or severe property damage.

Important Safety Instructions

When using an electric appliance, basic precautions should always be followed, including the following:

Read all instructions before using (this appliance).

- For household use only, not for medical or commercial purposes.
- This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or cognitive capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved.
- Children shall not play with the appliance.
- Cleaning and user maintenance shall not be made by children without supervision.

DANGER – To reduce the risk of electric shock:

- 1) Always unplug this appliance from the electrical outlet immediately after using and before cleaning.
- 2) Never use pins or other metallic fasteners with this appliance.
- 3) Carefully examine the covering before each use. Discard the appliance if the covering shows any sign of deterioration,

such as checking, blistering or cracking.

4) Do not use or store the appliance in wet environments (e.g. spa, bathroom, sink) or moist conditions, high temperatures or dusty places.

5) Do not reach for the appliance if it is submerged in water. Unplug immediately.

6) Do not place or store the appliance where it can fall or be pulled into a tub or sink.

7) Do not pour water into the appliance. Do not immerse the appliance in water or other liquids.

8) Do not operate the appliance or handle the plug with wet hands or wet feet.

9) Please attach the plug fully for use to avoid unnecessary electric shock or injury.

WARNING – To reduce the risk of burns, fire, electric shock, or injury to persons:

1. If the appliance is to be used by the following user groups, please consult a doctor, otherwise, accident, injury or discomfort may occur:

1) Persons suffering from thrombus, severe aneurysm, malignant tumor, venous tumor, skin infections, etc.

2) Persons suffering from osteoporosis, weak bones, abnormal or curved spine, or other acute pains

3) Persons suffering from diabetes or other sensory impairment

4) Persons have a pacemaker or medical devices implanted

5) Pregnant or postpartum women

6) Persons currently suffering from an injury to the lower extremities

7) Persons being treated for a medical condition or physical

injury

8) Persons with a heart condition, or abnormal blood circulation on feet

9) Persons with open wounds, discolored areas, or any part of the body that is swollen, burned, inflamed or where skin eruptions or sores are present

10) Persons had received cardiac or visceral surgery

11) Persons insensitive to heat or with heat-feeling loss

12) Persons suffering from a fever

13) Persons feeling discomfort

14) Persons with carsick, seasick feeling, or reacting vigorously to other means of transport

15) Persons that are easily susceptible to getting foot sprained or suffering injury

If you have any discomfort (such as inflamed skin, redness, itching, etc.), please stop using the appliance immediately and consult your doctor.

2. The appliance should never be left unattended when plugged in. Unplug from outlet when not in use, and before cleaning, putting on or taking off parts.

3. Do not use the appliance while covered in quilt or blanket or near heat resources (e.g. radiators or stoves). Excessive heating can occur and cause fire, electric shock, or injury to persons.

4. Close supervision is necessary when this appliance is used by, on, or near children, invalids, or disabled persons.

5. Never operate this appliance if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water. Return the appliance to a service center for examination and repair.

6. Do not carry this appliance by supply cord or use cord as a handle. Keep the cord away from heated surfaces.
7. Never operate the appliance with the air openings blocked. Keep the air openings free of lint, hair, etc.
8. Never drop or insert any objects into any opening.
9. Do not use outdoors. Place and use the appliance on a dry, hard and flat surface indoors.
10. Do not operate where aerosol (spray) products are being used or where oxygen is being administered.
11. Do not crush – Avoid sharp folds. Keep away from sharp objects to avoid scratching.
12. Do not stand on or in appliance, use while seated only.
13. Do not place excessively heavy objects on the appliance to avoid damage.
14. Do not allow pets and animals to lie or chew on this appliance or on the electric cord.
15. Do not use while sleeping/feeling asleep or after intake of alcohol.
16. Any modification or disassembly of the appliance without authorization is prohibited.
17. If a high temperature or foul odor occurs, please stop using the appliance immediately and cut off power.
18. To disconnect, turn all controls to the off position, then remove plug from outlet.
19. Only use accessories officially provided by RENPHO.

SAVE THESE INSTRUCTIONS

About the Product


1. What's in the Box

- 1 x Vibrating and Rolling Foot Massager
- 1 x Remote Control
- 1 x Power Adapter
- 2 x TPR Foot Pedal
- 1 x User Manual

2. Specifications

Product name: Foot Massager
Model: R-D200R
Input: 24V $\equiv$ 1.5A
Rated power: 30W
Auto-off timer: 10/15/20min


3. Product Overview

Using the Product

Connect the power adapter into the power port and plug into socket.

Control Panel

⏻ : Press to turn on the machine. P1 mode and 15min auto-off timer will be activated by default. TIME/MODE window will alternate between displaying P1 and countdown auto-off time. The SPEED window will display 05. Press again to turn off the machine.

P : Press to cycle through P1/P2 preset programs. TIME/MODE window will display corresponding mode and countdown auto-off time, and SPEED window will display 05.

P1: Vibration only

P2: Vibration and Rotation on

+ : Press to increase the vibration and rotation speed from 1 to 20.

- : Press to decrease the vibration and rotation speed from 20 to 1.

A : Press to cycle through A1/A2/A3 modes.

A1: Vibration only. Speed increases from 1 to 20 and decreases 20 to 1 (in cycle).

A2: Vibration and Rotation. Speed increases from 1 to 20 and decreases from 20 to 1 (in cycle).


A3: Vibration always on, Rotation on at 1-minute intervals. Speed will run in the cycle of 1→15→10→15→10→15→10→15→1.

T : Press to cycle through 10 / 15 (Default)/ 20-minute auto-off time. The timer will restart each time you press the button.

Note:

1. Each press is accompanied by a beep sound and description of the setting will show on the screen.
2. The speed can be adjusted only when the Program mode is on. If you press the Speed "+" or Speed "-" button when the Auto mode is on, the machine will switch to Program (P1) and Speed 5 by default.
3. During the massage, the timer will not restart unless the Time button is pressed.
4. The machine is programmed with a memory function. It will recall the last setting used once the Power button is turned on again. The memory function will not work if the power supply is unexpectedly cut off.
5. It is recommended to enjoy the massage for no more than 30 minutes at a time.

Remote Control

 : Press to turn on/off the machine.

 : Press to cycle through 10/15/20-minute auto-off time.



 : Press to increase the vibration and rotation speed from 1 to 20.

 : Press to decrease the vibration and rotation speed from 20 to 1.

 : Press to cycle through P1/P2 preset programs.

 : Press to cycle through A1/A2/A3 modes.

Tip: You can massage your feet or calf with the machine as shown.

Note: Do not stand on the machine, use while seated only.

Maintenance and Care

1. Unplug the appliance and allow it to cool before cleaning. Use a soft cotton cloth dampened with mild detergent to wipe it clean. Remove the foot pedals for easy cleaning.
2. It is prohibited for non-professional personnel to repair the appliance or replace parts.
3. If the power adapter is damaged, it must be replaced by the maintenance department or qualified personnel approved by the manufacturer so as to avoid potential hazards.
4. Do not use or store the appliance in wet environments or high temperatures.
5. To avoid damage, do not use sharp objects to handle the appliance.

Troubleshooting Guide

Common problems	Causes	Solutions
The massager does not work.	1. Disconnected power or loose plug causing poor contact. 2. No power or power cord is damaged. 3. The power button is off or the timer has been reset. 4. The motor auto shuts-off due to the device being on too long activating overheating protection. 5. Didn't use original adapter.	1. Remove plug from outlet and re-insert into power source. 2. Try another outlet. Check the power cord, replace if damaged. 3. Turn on the power button. 4. Use the massager after the motor cools off naturally. 5. Use original adapter provided in the package only.
The massage intensity is too high or too low.	The incorrect intensity level for your preference was selected.	The massager has 20 speed levels, choose the speed you prefer. We recommend low speed to begin with.
The speed can not be adjusted.	The speed can be adjusted only when the Program mode is on. If you press any speed button when the Auto mode is on, the machine will switch to Program (P1) and Speed 5 by default.	Activate the program function and adjust the speed.
Rotation doesn't work.	P1 mode or A1 mode is activated. Both don't include the rotation function.	Choose other modes other than P1 or A1 mode.

Warranty Policy

Your RENPHO product purchase is covered by a one-year limited manufacturer warranty from the date of delivery.

For warranty terms and conditions, please visit:

<https://renpho.com/pages/warranty-terms-and-conditions>

Note: Product registration is not required for the warranty.

If you choose not to register your product, it will not diminish the product warranty.

Customer Service

Please feel free to contact us if you have any questions or concerns. RENPHO Customer Service Team guarantees a quick response and hassle-free solutions to any issue you may have within business hours.

 **Tel:** +44(0)7434-666088 (UK ONLY)
Monday-Friday 9:00AM-5:00PM(GMT)

 **Email:** support-eu@renpho.com(UK&EU)

*For defective products or the return of items, please contact us with your order number within the specified warranty period. DO NOT dispose of any product parts as they may be required for inspection/repair.

Notice:

This marking indicates that this product should not be disposed with other household wastes throughout the EU. To prevent possible harm to the environment or human health from uncontrolled waste disposal, recycle it responsibly to promote the sustainable reuse of material resources. To return your used device, please use the return and collection systems or contact the retailer where the product was purchased. They can take this product for environmental safe recycling.

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Manufacturer: Shenzhen Ruiyi Business Technology Co., Ltd.
EU REP C&E Connection E-Commerce (DE) GmbH
UK REP Richers & Vision Ltd.