89mm


129mm

hetaida

Non Contact

Infrared Body Thermometer

HTD8813C




QUICK USER GUIDE


The Non Contact Infrared Body Thermometers are designed to be used for intermittent measurement and monitoring of human body temperature by consumers in the home setting and primary care setting as screening tool.

Measurement steps:

1. Aim the thermometer probe at the center of the forehead and keep the distance within 1.97in. Do not touch the forehead directly.



2. Simply press the measurement key to start the measurement.




3. You will hear a beeping sound. Read the number on the display. The whole process takes about 2 second.



*After the measurement is over, the thermometer will automatically turn off after 18 seconds.

The thermometer is a precision instrument. If there are small foreign objects such as floating dust and droplets falling into the probe, it should be cleaned up in time. It is easier to obtain accurate values by keeping the probe clean.




* Please wipe it with a damp cotton swab first when cleaning, then dry it with a dry cotton swab, then let it sit for a few minutes for normal use.

How to switch measurement mode?

In the power-on state, short press the mode key (MODE) to switch between the " Body " mode and the " Surface " mode. The measurement range of the surface mode is: 32~212℉. The Body mode is used to measure the body temperature in the human body. (The factory default is " Body " mode).

How to install battery?

Push the battery door down and back, insert 2 AAA alkaline batteries, and close the battery cover. When installing, please carefully check the positive and negative signs inside the battery cover, and note that the positive and negative electrodes of the battery cannot be installed in reverse.



* The "+" sign on the convex part of the battery is the positive electrode, and the "-" sign is the negative electrode.

Special Reminder

1. The measurement position of the thermometer is the center of the forehead, please do not measure other parts.
2. When the thermometer is brought into a new environment (such as bringing it home from the outdoors), please let it stand for 30 minutes to reach an equilibrium state before measuring.
3. When measuring for a baby, crying and bathing will cause body temperature fluctuations, and it is not recommended to measure immediately at this time.
4. When the forehead is covered with bangs, antipyretic stickers, and sweat, the local skin temperature will be abnormal. At this time, the measurement results are not for reference. Please keep the forehead bare and clean for at least 30 minutes before measuring.
5. When measuring with a thermometer, please keep away from electromagnetic radiation products, such as mobile phones that are talking.

*If there are any questions, please contact us before leave a feedback.
ttt

267mm