PISOOKA X

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projectmars.shop

WARRANTY

The warranty period of this product is 12 months, and the right is NOT transferable. From the purchase date, you enjoy free warranty for non-artificial damage within the warranty period of 12 months with valid warranty card. Please contact Project Mars Customer Service Provider first if needed maintenance. The provider will indicate to return the device, the original copy of warranty card, the proof of purchase and the other maintenance information back to the maintenance center, and the customer will be responsible for the back freight. The maintenance should take 5-7 working days after the provider receives the device, and there is no same-day service.

The warranty is limited to the device, excluding consumables and accessories, such as massage heads and charging cables, etc. The warranty also does NOT cover damage under the following circumstances:

1. Artificial damage by accident, improper usage or transportation.
2. Operation without following the manual instructions or maintain improperly.
3. Disassembled by engineers who were not approved by the Project Mars.
4. Unable to provide valid proof of purchase.

Project Mars CO LTD

ADR:

Rm. 3, 17F., No. 193, Jingxing Rd., Wenshan Dist.,

Taipei City 116051, Taiwan (R.O.C.)

1F, Newmarket House, Newmarket Square, The Coombe, Dublin 8. D08VAP9

EML: mars-service@projectmars.com.tw

projectmars.info

Trouble Shooting	
Product does not work	Check if the device is sufficiently charged. If the device hasn't been used for a long time, it must be charged for 8 hours.
Abnormal noise during use	Slight noise is a normal phenomenon when the longer the time of use, and you can shake the device for preliminary inspection. However, if there is a knock, there may be loose parts, please send the device back for examination.

MAINTENANCE

1. If it is not used for a long time, it's highly suggested to store the device in the original package in a safe, cool, and dry place for save-keeping with routine recharging around once every 2 to 3 months to keep 75% of the power.
2. Please choose neutral detergent and soft cloth for routine maintenance. Do not use engine oil, organic solvents or other chemicals for cleaning.
3. Please turn off the power after using the product or before cleaning the product.
4. The product should be secured and safeguarded during transportation to avoid damage from collision or falling.

ENVIRONMENTAL PROTECTION

Equipment name: Mini Massage Gun				Type designation(Type):PISOOKA X		
Unit	Restricted substances and its chemical symbols					
	Lead (Pb)	Mercury (Hg)	Cadmium (Cd)	Hexavalent chromium (Cr^{+6})	Polybrominated biphenyls (PBB)	Polybrominated diphenyl ethers (PBDE)
Shell	○	○	○	○	○	○
Circuit board	-	○	○	○	○	○
Metallic parts	-	○	○	○	○	○
Plastic parts	○	○	○	○	○	○

Note 1: "Exceeding 0.1 wt %"and "exceeding 0.01 wt %"indicate that the percentage content of the restricted substance exceeds the reference percentage value of presence condition.

Note 2: "○"indicates that the percentage content of the restricted substance does not exceed the percentage of reference value of presence.

Note 3: The" - "indicates that the restricted substance corresponds to the exemption.


The device has a built-in charged battery, and the charged battery contains substances that may pollute the environment. Please dispose of the replaced battery in the recycling bin designated by the government.

CHARGING INTRODUCTION

1. As discharge test has been performed when leaving the factory, the product can be used upon purchase.
2. For recharging, turn off the device and connect the cable to the Type-C charging port. The status indicator will flash in sequence to indicate the charging process has started. DO NOT move the device when charging. The status indicator will show a steady full bar to indicate that the device has been fully charged.

REPLACING THE MASSAGE HEAD

To replace the massage head, turn off the device and slowly rotate the head to insert it into the barrel. As shown in the figure, the red line position of the massage head should be aligned with the red line position of the machine.

SPECIFICATIONS

Net weight	350 g	**Size**	145x 78x 40 mm
Rated speed	1600-3200 rpm	**Stroke**	9 mm
Battery capacity	800 mAh	**Rated input**	Type-C, 5V ⎓ 2A
Charging time	2.5 hours	**Endurance time**	4 hours*

*This is the result of no-load test in the laboratory. Actual time of usage may vary in accordance of how much pressure is applied to the device by different user.

FAQ

1. Low temperature:

 Temperature below 0°C (32°F) may intervene the electrochemical reaction of the battery, resulting in failure when trying to turn on the power or start vibration.

 Solutions: Keep the device in temperature above 0°C (32°F) for two hours, or using hair drier for quick result.

2. Low battery power:

 When the battery level is low or haven't been charged for a long time, it will enter low power protection mode. The product cannot operate normally in said mode.

 Solutions: Turn off the product and charge it. It will work properly after fully charged.

PRODUCT INTRODUCTION

MASSAGE HEAD

Spinal Head
For spine massage and relaxation.

Round Head
For whole body muscle massage and relaxation.


Flat Head
For large muscle group massage and relaxation.

Bullet Head
For pain point massage and relaxation.

INSTRUCTION FOR USE

OPERATING INTRODUCTION

1. Press and hold the multi-function button for 3 seconds to turn on the power. The status indicator will show the current battery level.
2. Short press the multi-function button to start the vibration. The status indicator will show the current vibration strength. Each press will increase the vibration strength by one gear.
3. Press and hold the muti-function button for 3 second to shut down the power.
4. The device is equipped with automatic power-off protection mechanism, which will stop the vibration after vibrating for approximately 10 minutes. Putting the device on standby approximately 10 minutes, as in turn on the power but without vibrating, will also automatically turn the power off.

PRECAUTIONS

To avoid harm, fault, electric shock, fire, death and other equipment damage, please read the instructions and warnings carefully.

PROHIBTION

1. When using the device, Do Not press the massage gun. If the motor is squeezed, the use period will be reduced.
2. To prevent injury from being caught between, please ensure that NO finger, hair, or other body part is between the attachment and the device.
3. The device should NOT be near fire or immersed in water or any liquid. Please return the device for examination immediately if dropped into water by accident.
4. DO NOT drop the device intentionally or attempt to disassemble the device.
5. When charging the device, DO NOT use anything other than the original charger. The device should be operated on an electric power source matched to the safety extra low voltage indicated on the product. Failure to do so may cause injury, electric shock or fire.
6. The device should NEVER be used on bones or around genitals. Use this product for soft tissues and muscles only. Please operate with caution when used on the head, the neck, or around the heart.
7. The device should NEVER be an alternative to any medical treatment. For those suffering muscle pain, under other medical treatment, or have other health conditions, please consult a doctor before using.
8. If the power is on without a massage head attached to the device, DO NOT touch the massage head socket so as to prevent injury.

9. This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety. Children should be supervised to ensure that they do not play with the appliance.
10. It is recommended that the time per usage DOES NOT exceed 20 minutes in total. Usage on the same body part SHOULD NOT exceed 5 minutes.

DISCLAIMER

Please consult your doctor before using the device if you have any contraindication or are under any circumstance listed below:

1. Pregnancy, diabetes, chronic diseases, neurological diseases, retinal detachment, and other complications that may raise concern.
2. Epilepsy, migraine, prolapse of lumbar intervertebral disc, anterior displacement of bone cone and other symptoms.
3. Recent surgeons, using pacemakers, recent joint replacement surgery, Intrauterine contraceptive device placing surgery, metal parts implantation, or other concerns about your own body.

Researches (and ongoing researches) on the effect of vibration therapy toward specific diseases showed that vibration therapy may ease many symptoms. However, employing massage gun or vibration therapy as physiotherapy must be under advice, agreement, or guidance of a doctor or physiotherapist.

Thank you for choosing our PISOOKA massage gun. Before using it, please read this manual carefully to ensure proper operation, and keep this manual for future reference.

Main functions:
Relieving muscle soreness, stimulating the blood circulation, relieving stiffness of cervical vertebra, and improving muscle recovery.

CONTENTS

Precautions / 01

Product Introduction / 03

Instruction for Use / 04

Specifications / 06

FAQ / 06

Maintenance / 07

Environment Protection / 08

Warranty / 09