FAQ

After applying Cushion Grip, how long do I have to wait to eat, drink, or take dentures out?

Answer

At least 2 hours. The longer the better. Cushion Grip takes an average of 2-6 hours to fully cure.

Tip: Our brand ambassadors like to apply it at night before bed so it can cure the whole night uninterrupted.

FAQ

Can I soak my dentures with Cushion Grip on them?

Answer

Yes, but do not soak for the first 24 hours after application. After that, you can clean and soak like normal because Cushion Grip is 100% waterproof and will not wash off or break down in water.

How much Cushion Grip should I use for my denture?

Answer

The quantity needed to secure your dentures properly is based on the looseness of the denture.

For instance, it can depend on if you have a single or full set of dentures and how often you'd like to change your application. Most people get 3-4 applications per tube.

Application use is dependent on the user's number of plates (an upper, bottom) or partial. as well. The instructions for application state to squeeze out a ball of Cushion Grip 3/4 inch in diameter per plate after the tube has been warmed in hot tap water.

FAQ

Is Cushion Grip a glue adhesive?

Answer

Cushion Grip is not a glue adhesive; it works similarly to a soft reline by filling in the gaps and spaces between your gums and dentures. It adheres to your dentures, not to your gums; it only secures the fit with your gums.

When used correctly, Cushion Grip will give you more comfort and a better fitting denture along with saving you money.

What do I do when the Cushion Grip has hardened in the tube?

Answer

A thermoplastic dental adhesive requires heat to produce a softer pliable adhesive that provides a soft cushion for gums when applied to a denture.

It hardens as it cools creating a stronger hold when compared to ordinary cream/paste dental adhesives. Cushion Grip Thermoplastic Denture Adhesive holds up to 4 days versus other dental adhesives that last 8-12 hours.

Ensure tube is placed in hot water upside down with cap tightly closed to facilitate easy removal from the tube. At room temperature, Cushion Grip is a highly viscous semi-solid material that becomes pliable with heat.

FAQ

I am having a difficult time with the application and/or removal process, what should I do?

Answer

Please call 1-800-277-6151 and ask for a Brand Ambassador. They will be able to help you with the application and/or removal process and answer questions regarding Cushion Grip.

Why does it tingle? What to expect after applying Cushion Grip to your dentures?

Answer

A temporary slight stinging or burning sensation like mouthwash may be experienced at first application even with healthy gums but should last no more than 5 minutes.

Cushion Grip contains antibacterial ingredients that inhibit the formation of biofilm. The tingling sensation is normal and only lasts a few minutes.

This product contains a small amount of alcohol. Do not use it if you have any open sores or cuts in your mouth or gums. Please wait to use until after your extraction sites are healed.

How do you store and what is the expiration date of Cushion Grip?

Answer

Please store Cushion Grip at room temperature. Keep out of cold climates to avoid hardening. Keep the cap on tightly to prevent product hardening at the neck.

When stored appropriately Cushion Grip will not expire.

Can you use Cushion Grip on a partial?

Answer

Yes, Cushion Grip Denture Adhesive is great for both dentures and partials.

FAQ

Can you use other adhesive/reliner with Cushion grip?

Answer

The product is recommended to be used without additional adhesives or reliners; however, some customers prefer to use Cushion Grip with other products for extra security.