Kopfkissen
Sleeping pillow

DE EN

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EN

Intended Use:

- Ergonomic support for the neck, head and shoulders while you sleep
- Promotes anatomically correct alignment of the Spine
- May help relieve neck and back pain

Intended user group: Persons of adult age

Indications:

- For people with neck or back pain and unfavourable sleeping posture
- Supports incorrect posture of the spine and muscular tension
- For postoperative support after consultation with the treating doctor

Contraindications:

- Not suitable for acute injuries to the neck or back that require special medical treatment
- Do not use if you have known allergies to the materials of the pillow (e.g. memory foam)
- Do not use if a specific sleeping pillow or sleeping position has been prescribed by your doctor
- Be careful with respiratory conditions that require special head or neck positioning; consult your doctor

Instructions for use:

- Determine your optimum pillow height for restful sleep
- Open the zip and remove any foam layers for your personal feel-good sleeping pillow
- The head should be in a straight line to the rest of the body and the neck must not be bent or overbe stretched
- This maintains the natural double S-shape of the spine and ensures healthy sleep for your back
- Cover your sleeping pillow again and place it in the optimum position for you

Safety Instructions: ⚠

- Use the pillow only for its intended purpose
- If symptoms persist or worsen, stop use and seek medical advice
- Check the cushion regularly for damage and replace if necessary
- Not suitable for children under 3 years due to choking hazard
- Care for and clean the pillow according to the manufacturer's instructions to ensure hygiene and longevity
- Store the pillow away from open flames orsources of strong heat

Specification of technical data:

- Dimensions (width x height x depth): 400 x 600 x 120/100 mm
- Weight: 1250 g (soft), 910 g (medium), 1250 g (hard)
- Density: 45 kg/m^3 (soft), 30 kg/m^3 (medium), 45 kg/m^3 (hard)
- Cores: 100% Polyurethane
- Outer cover: 70% Polyester + 30% Viskose

Information on cleaning:

- Open the zipper and remove the foam cores
- Wash the cover at 40 °C in the washing machine (core not washable)

Instructions for safe disposal:

- Not recyclable, disposal as plastic waste

Storage and transport conditions:

- Constant conditions, ideally at temperatures of 18 to 24 °C and relative humidity of 40 to 60%.

Intended service life:

- 5 years

Information to authorities in case of "events":

- All serious incidents that have occurred in connection with the wedge pad must be reported to the manufacturer and the competent authority:

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Further information on the application:

- In case of problems with the application, please ask your specialist dealer or the responsible feela. customer advisor for advice. Phone: +49 991 9992267

Scope of delivery:

- 1 x feela.-sleeping pillow, consisting of 3 elements
- 1 x cushion cover
- Also: instructions for use, protective packaging

DE / EN / FR / IT / ES
Waschkennzeichnungen für den Bezug
Washing indicators for the pillowcase
Indications quant au lavage de l`enveloppe extérieure
Etichette di lavaggio per la copertina
Etiquetas de lavado para la cubierta

	Maschinenwäsche 40°C / Machine wash 40°C / Lavage en machine à 40°C / Lavare in lavatrice a 40°C / Lavado a máquina a 40°C
	Nicht bleichen / Do not bleach / Ne pas javelliser / Non candeggiare / No blanquear
	Bei niedriger Temperatur im Wäschetrockner trocknen, bis es vollständig trocken ist / Tumble dry at a low setting until completely dry / Sécher à feu doux jusqu'à séchage complet / Asciugare in asciugatrice a fuoco basso fino a completa asciugatura / Secar en secadora a fuego lento hasta que esté completamente seco sécher en machine / Non asciugare in asciugatrice / No secar en secadora
	Nicht bügeln / Do not iron / Ne pas repasser / Non stirare / No planchar
	Nicht trockenreinigen / Do not dry clean / Ne pas rettoyer à sec / Non lavare a secco / No limpiar en seco

DE / EN / FR / IT / ES
Erklärung der verwendeten Piktogramme
Explanation of the pictograms used
Explication des pictogrammes utilisés
Spiegazione dei pittogrammi utilizzati
Explicación de los pictogramas utilizados

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Typesetting and printing errors expected.
Changes in colour, material and construction possible.