Disinfectant

Retainer Brite®

EN

INSTRUCTIONS FOR USE

INTENDED USE/PURPOSE:

A disinfectant and cleaner used on a variety of removable, custom-made dental devices.

INDICATIONS FOR USE:

Product is indicated for patients who are using removable dental appliances such as retainers, aligners, dentures, sports mouthguards, nightguards, bruxing and TMJ appliances, which need to be cleaned and disinfected at home.

PERFORMANCE CHARACTERISTICS AND CLINICAL BENEFITS:

Retainer Brite® is specially formulated to help protect removable appliances by disinfecting*, removing plaque and preserving clarity and cleanliness of the appliance.

Clinical benefits of orthodontic treatment may include overall improvement of tooth health and correction of malocclusion.

CONTRAINDICATIONS:

This product is contraindicated for patients allergic to Persulfates.

WARNINGS:

Contains persulfates which are a known allergen.

Keep out of reach of children.

Use only as directed. If under the care of a dental professional, follow the instructions provided by your dentist and notify them if any skin irritation occurs.

If product comes into contact with eyes, wash eyes thoroughly with water for several minutes.

Do not put tablets or cleaning solution in mouth. If ingested (swallowed) or inhaled, seek medical assistance immediately. In the USA, please call Poison Control Center at 1-800-222-1222. In the UK, please call NHS 111.

Report serious incidents that have occurred in relation to the device to device manufacturer and the Competent Authority in Member State of residence, or other applicable regulatory agency.

PRECAUTIONS:

**Do not allow wire appliances to soak for more than 15 minutes. Remove promptly and rinse thoroughly with water.

Wash hands thoroughly before and after handling tablets.

Do not mix with other cleaners.

DO NOT REUSE, once dissolved tablet cannot be reused. Solution will not clean as intended if reused.

ADVERSE REACTIONS:

Allergic reactions.

STEP BY STEP INSTRUCTIONS:

Use daily for best results.

1

Drop one tablet into 100mL (0.4 cups or 3.3 oz) water at or above 50°C (122° F).

2

Place appliance into bubbling blue solution.

3

After 15 minutes** appliance is clean and odor free!

4

Rinse appliance and hands thoroughly with running water.

STORAGE:

Keep Dry

SAFE DISPOSAL OF THE DEVICE:

Dispose of down the drain, as you would your mouthwash or toothpaste. Flush drain with plenty of running water.

**Bactericidal according to ASTM E2315-031[1], EN 13727, EN 14561; Yeasticidal according to EN 13624 and EN 14562; when dirty dental devices are immersed in 100 mL water/tablet solution at 50 °C for 15 minutes (data on file).*

[1] for US market

INGREDIENTS: Sodium Bicarbonate, Corn Syrup Solids, Sodium Carbonate, Citric Acid, Potassium Persulfate Compound, Sodium Percarbonate, Sodium Sulfate (may contain silica), Sorbitol, Tetraacetylethylenediamine (TAED), PEG-180, Sodium Lauryl Sulfoacetate, PEG-8, Aroma (Flavor), Magnesium Stearate, FD&C Blue #1 Aluminum Lake (CI 42090), FD&C Blue #2 (CI 73015).

Manufacturer:
Dentsply Sirona Orthodontics Inc.
7290 26th Court East
Sarasota, FL 34243 USA

www.dentsplysirona.com/ifu

MKT-DFU0048-ML Rev. 10 Issued 12-2024

© 2024 Dentsply Sirona

CE 0123

MD

EC REP **Authorized EU Representative**
Emergo Europe
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands

CH REP **Authorized Switzerland Representative**
Maillefer Instruments
Holding Sàrl
Chemin du verger 3
CH-1338 Ballaigues

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