PILL MILL

Thank you for purchasing the Pill Mill Pill Crusher.
Please read the manual below thoroughly before use.

INSTRUCTIONS:


It is recommended to wash the Crusher before first use (Instructions are on page 2).

1- Take off the LID (B) and place the pills you wish to crush inside the GRINDER BODY (G), and then place the LID (B) back on and mount the HANDLE (A) on top of the METAL ROD (C).

2- Hold the crusher firmly with one hand and turn the handle clockwise to start crushing. Crush until you don't feel any resistance on the handle anymore. To prevent pills from getting stuck, turn the handle counterclockwise for a turn every now and then while crushing to release the pills.

3- To adjust the final size of the crushed particles, find the ADJUSTMENT KNOB (I) that is located under the BURR/GRINDER (H). Turn the knob clockwise for finer powder, or counterclockwise for a coarser powder. Before adjusting the grade of crushing, make sure there are no bigger particles between the burrs that are left over from the previous crushing.

Important: Do not tighten the ADJUSTMENT KNOB (I) too much or the ceramic blades will come in contact with each other, resulting in the handle becoming difficult to turn. The blades may also become damaged if they make contact. Stop tightening just before this point for the finest grade of crushing.

HOW TO CLEAN THE CRUSHER?

It is recommended that you clean the burrs every other week with a soft brush (toothbrush or similar) and water to prevent residue buildup on the burrs. Overtime, buildup can cause a poorer grind and result in crusher components becoming stuck. **This routine cleaning is very important for your crusher's ability to stay in working order.**

Please make sure to dry all parts very well in order to prevent the risk of rusting. For even more detailed instructions on cleaning and reassembling your crusher, go to www.pill-mill.com or contact us at **customers@pill-mill.com** for help. You are also welcome to watch our "How to wash your grinder video" on YouTube by following this link **https://www.youtube.com/watch?v=mWyDg_NvcJE** or by scanning this QR Code.

HOW TO DISASSEMBLE AND REASSEMBLE YOUR CRUSHER FOR CLEANING:

Disassembling:

Turn the ADJUSTMENT KNOB (I) counterclockwise while holding the handle until the knob comes off, followed by the SPRING (E). There is no need to take out the METAL ROD (C). If you wish to take the rod out, you will need to slide WASHER 2 (F) off the METAL ROD (C), by aligning the grooves of the washer with the grooves of the rod. Then, pull WASHER 2 (F) off and the rod will come out. Then, remove WASHER 1 (D) from the rod if you wish.

Reassembling:

(Follow STEP 1 and 2 if you took the METAL ROD (C) out, if not skip to STEP 3)

STEP 1: Slide plastic WASHER 1 (D) on METAL ROD (C) with the square side of the washer pointing toward the bottom of the rod. Make sure the washer grooves line up so that it slides over the METAL ROD(C). Push Washer to the top of METAL ROD (C).

STEP 2: Put METAL ROD (C) into the top half, slide down and slide WASHER 2 (F) on METAL ROD (C), with the square side of the washer pointing towards the top of the rod. Slide it up from below, again aligning the grooves, then push the washer up as far up as it goes.

STEP 3: Put SPRING (E) onto the bottom of the METAL ROD (C) in the top section.

STEP 4: Slide BURR/GRINDER (H) onto metal rod under SPRING (E).

STEP 5: Hold top of METAL ROD (C) (square top) either with your hand or with the HANDLE/CRANK (A) and press the BURR/GRINDER (H) onto the METAL ROD (C).

STEP 6: Screw on ADJUSTMENT KNOB (I) with one hand while holding the HANDLE/CRANK (A) or top of the rod with your other hand. Turn until tight. Turn crank to check if it's too tight.

STEP 7: Place top half of grinder onto BOTTOM (J).

STEP 8: Put on LID (B), then attach HANDLE/CRANK (A).

FAQ:

What to do if a pill gets stuck in the crusher?

In case a pill becomes stuck in the crusher, simply turn the handle in a counterclockwise direction until the pill becomes loose. You can remove the LID (B) to more easily notice when the pill becomes dislodged. If this does not work, turn the ADJUSTMENT KNOB (I) towards the coarser size setting and crush the pill. You can also use a sharp object to help release the stuck.

Crushing a bigger pill:

There are two ways to crush a larger pill. One option is splitting the pill into smaller pieces first and then crushing the pieces. The second option is letting your crusher do the breaking for you by setting the size of the crushing to the coarsest setting and crushing the pill into smaller pieces. Once the pill is broken into smaller sizes, adjust the size to the finest setting and crush the coarsely crushed pill again.

How does one consume crushed pills?

If your medication allows for it (consult your doctor first or check the medication specification), the easiest way to consume crushed pills is to mix the powder with food. It is best to mix the powder with any kind of fruit puree (like applesauce), yogurt, or pudding. This ensures an easy intake and covers up the bitter taste some pills might have, making your daily routine more pleasurable.

Is the crusher able to crush all type of pills?

It is not recommended to crush gelatin medication, capsules, gel caps and soft medication or gummies. It is able to crush pills that have a protective layer, and may leave the protective layer in bigger pieces. That is normal and desired by some.

If you desire a finer powder than what is produced, simply pour the contents of the jar back into the crusher and repeat the process.

Disclaimer: Consult your physician or a specialist before crushing your medication. We do not take any responsibility for misuse of your medication.

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