Seat Cushion for pressure relief

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DE EN ES FR IT

Instructions of use

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Healthfix®

Healthfix GmbH
Bentheimer Str. 66
48529 Nordhorn
Germany
service@gohealthfix.com

REF

Art. 100 110

HFXSITZ01

MD

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SEAT CUSHION
WITH COCCYX AND SIT BONES CUTOUT

EN

Purpose of use:

The **Healthfix® seat cushion** is used to correct sitting posture and improve sitting comfort. It supports the natural curvature of the back, promotes an upright posture and can help reduce back pain and discomfort. The U-shaped neckline relieves pressure on the tailbone. The two middle cutouts help reduce pressure on the sit bones.

Ideal for office workers, people who spend a lot of time sitting and anyone who wants to improve their sitting posture ergonomically. Suitable for various seating furniture such as office chairs, dining chairs or in the car.

Operating instructions:

-Place the seat cushion on a flat seat surface with the matt anti-slip nubs facing down and the recess for the tailbone facing backwards.

-Make sure the U-shaped cutout is level with your spine.

-Placement: Place the cushion on a chair or seat of your choice. It should fit tightly and securely without slipping.

After just 30 seconds of using the **Healthfix® seat cushion**, the MEMORY-FOAM will adapt to your body shape with the help of your body heat.

We still recommend that you get up regularly.

Indications:

Poor posture, back pain, coccyx and ischial tuberosity problems, piriformis syndrome, sciatica problems, prevention of postural problems, rehabilitation after operations or injuries, pressure sores

Contraindications:

Lymphatic drainage disorders, sensory and critical circulatory disorders in the supplied body regions, PAD

unstable fractures in the spine and pelvic area.

Medical conditions for which the treating doctor or therapist advises against.

Storage and transport conditions:
Constant conditions, ideally temperatures between 18 to 25° degrees and a humidity of 35 to 65%.

Cleaning information:
Open the zipper on the bottom and one inside of the middle recess and take out the foam core. The foam core should not be washed as this could cause it to lose its effectiveness.
The cover can be washed in the washing machine at 40°.

Scope of delivery:
1x Healthfix® seat cushion
1x seat cushion cover
Instruction manual and brochure

Warranty and customer service
The Healthfix® seat cushion comes with a 12-month guarantee.
If you have any questions or problems, please contact our customer service at the email address provided.

Intended useful life
- 5 years

Safety instructions: ⚠
Keep the seat cushion away from fire-inducing heat, sparks or open flames.
Only use the product as described in these instructions.
Make sure that the pillow sits firmly and does not slip.
The pillow is not a toy. Keep the pillow out of the reach of children.

Instructions for safe disposal:
Dispose of the seat cushion as plastic waste

Further information about the application:
If you have any difficulties in use, please contact your dealer or the responsible Healthfix® customer advisor.
Email: service@gohealthfix.com

Technical data:
Dimensions: W45 x H7.8 x L39 cm
Weight: 0.8 kg
User weight: 50-135 kg

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All rights reserved.
Healthfix GmbH | Bentheimer Str. 66 | 48529 Nordhorn | Germany