

What are the features of this product ?

GOQii Smart Vital a cutting-edge smart device with a host of features such as Integrated Pulse Oximeter (SpO2), 24\*7 Body Temperature & Heart Rate Monitor, Blood Pressure, Step Count, Auto Sleep Tracking, 18 exercise modes, music control, and phone notifications. Additionally, you get 3 months personal health coaching with one doctor consultation.

How accurate is GOQii ?

GOQii Smart Vital is tested by various labs, it is International Organization for Standardization (ISO) certified with multiple grades and is approved as a registered medical device by the Central Drugs and Standard Control Organization (CDSCO). The readings on Heart Rate, Body Temperature, Blood Oxygen (SpO2) are approximately 97% accurate & Blood Pressure close to 94% accurate. We use PPG (photoplethysmography) which is a proven standard. If you compare it to devices using Oscillometric method or to the traditional Sphygmomanometer, you will notice a variance. Do read & research more on each method on the internet and why the medical fraternity recommends all. You need to always use devices like ours to help monitor trends.

IS GOQii certified ?

Yes. GOQii Smart Vital is ISO certified under the following grades:

**ISO 9001** standard is based on a number of quality management principles including a strong customer focus, the motivation and implication of top management, the process approach and continual improvement. These principles are explained in more detail in ISO's quality management principles. Using ISO 9001 helps ensure that customers get consistent, good-quality products and services, which in turn brings many business benefits.

**ISO 13485:2016** specifies requirements for a quality management system where an organisation needs to demonstrate its ability to provide medical devices and related services that consistently meet customer and applicable regulatory requirements.

**ISO 27001** is widely known, providing requirements for an information security management system (ISMS), though there are more than a dozen standards in the ISO/IEC 27000 family. Using them enables organizations of any kind to manage the security of assets such as financial information, intellectual property, employee details or information entrusted by third parties.

Is GOQii medical registered?

Yes. GOQii Smart Vital is approved as a registered medical device by the Central Drugs and Standard Control Organization (CDSCO) (GOQii-MUMBA-MH/M/MD/001104).

Battery life of GOQii ?

GOQii Smart Vital has a 230 mAh battery which provides a battery life of upto 7 days, depending on the usage. If all functions are enabled, it will last for upto 3-4 days. You can extend the same to 7 days by lowering brightness, disabling notifications, continuous Heart Rate and Temperature Monitoring.

How to connect GOQii and app

Ensure your mobile phone's Bluetooth, Internet & GPS is ON to link the Tracker. Open Apps & Devices >> Tap on GOQii Smart Vital >> Tap "SETUP NOW" Button. Keep the Tracker near to Phone >> Tap on "START SEARCH" button. Once the device is searched, it will say "GOQii Device found." >> Tap on the "SEND A LINK REQUEST" button. Link request code will be sent to your GOQii Smart Vital Tap on the screen of GOQii Smart Vital to accept the Link request. After tapping on the Link request, it will verify the device availability and it will link the GOQii Smart Vital with the phone successfully.

Charger pin type ?

Yes. Magnetic two pin charger is provided with this device.